



Photometer-System MultiDirect



GB Instruction manual

Safety precautions



Reagents are formulated exclusively for chemical analysis and must not be used for any other purpose. Reagents must not get into the hands of children. Some of the reagents contain substances which are not entirely harmless environmentally. Be aware of the ingredients and take proper care when disposing of the test solution.



Please read this instruction manual before unpacking, setting up or using the photometer. Please read the method description completely before performing the test. Be aware of the risks of using the required reagents by reading the MSDS (Material Safety Data Sheets). Failure could result in serious injury to the operator or damage to the instrument.

MSDS: www.tintometer.de



Use the charger unit only with rechargeable batteries. Failure can result in serious injury to the operator or damage to the instrument.

Do not use charger with non rechargeables batteries.



The accuracy of the instrument is only valid if the instrument is used in an environment with controlled electromagnetic disturbances according to DIN 61326. Wireless devices, e.g. wireless phones, must not be used near the instrument.

Revision_7 08/2007

Table of contents

Part 1 Methods	7
1.1 Table of Methods	8
Acid demand to pH 4.3.....	12
Alkalinity-total (Alkalinity-m, m-Value).....	14
Alkalinity-p (p-value)	16
Aluminium with tablets.....	18
Aluminium (powder pack).....	20
Ammonium with tablets	22
Ammonium (powder pack)	24
Ammonium, low range (LR)	26
Ammonium, high range (HR)	28
Boron	30
Bromine.....	32
Chloride.....	34
Chlorine.....	36
Chlorine with tablet	
differentiated determination (free, combined, total)	38
free Chlorine.....	40
total Chlorine.....	41
Chlorine with liquid reagent	
differentiated determination (free, combined, total)	42
free Chlorine.....	44
total Chlorine.....	45
Chlorine (powder pack)	
differentiated determination (free, combined, total)	46
free Chlorine.....	48
total Chlorine.....	49
Chlorine dioxide.....	50
in presence of Chlorine	52
in absence of Chlorine	55
Chlorine HR (KI)	56
COD, low range (LR)	58
COD, middle range (MR).....	60
COD, high range (HR)	62
Copper with tablet.....	64
differentiated determination (free, combined, total)	61
free Copper	66
total Copper	67
Copper PP.....	68
Cyanide	70

Cyanuric acid	72
DEHA T.....	74
DEHA PP	76
Fluoride	78
Hardness, Calcium	80
Hardness, total.....	82
Hardness, total HR	84
Hydrazine	86
Hydrazine with liquid reagent	88
Hydrazine with Vacu-vials.....	90
Hydrogen peroxide	92
Iodine	94
Iron.....	96
Iron with tablet	98
Iron (powder packs)	100
Iron (TPTZ) (powder packs)	102
Mangenes with tablet	104
Manganese LR (powder packs).....	106
Manganese HR (powder packs).....	108
Molybdate with tablet.....	110
Molybdate HR (powder pack).....	112
Nitrate	114
Nitrite with tablet.....	116
Nitrite LR (powder pack).....	118
Nitrogen, total LR.....	120
Nitrogen, total HR.....	122
Oxygen, active	124
Oxygen, dissolved	126
Ozone.....	128
in presence of Chlorine	130
in absence of Chlorine	132
PHMB (Biguanide).....	134
Phosphate.....	136
Phosphate, ortho LR with tablet	138
Phosphate, ortho HR with tablet	140
Phosphate, ortho (powder packs).....	142
Phosphate, ortho (tube test).....	144
Phosphate 1, ortho C	146
Phosphate 2, ortho C	148
Phosphate, hydrolysable (tube test)	150
Phosphate, total (tube test)	152

pH-Value LR with tablet	154
pH-Value with tablet	156
pH-Value with liquid reagent.....	158
pH-Value HR with tablet.....	160
Potassium	162
Silica	164
Silica LR PP.....	166
Silica HR PP.....	168
Sodium hypochlorite	170
Sulfate T	172
Sulfate (powder pack).....	174
Sulfide	176
Sulfite.....	178
Urea	180
Zinc	182
1.2 Important notes.....	184
1.2.1 Correct use of reagents	184
1.2.2 Cleaning vials and accessories for analysis	185
1.2.3 Guidelines for photometric measurements	185
1.2.4 Sample dilution techniques	187
1.2.5 Correcting for volume additions	187
Part 2 Operating manual.....	189
2.1 Operation.....	190
2.1.1 Commissioning	190
2.1.2 Saving data – Important Notes	190
2.1.3 Replacement of rechargeable batteries resp. Lithium-battery	190
2.1.4 Charging the rechargeable batteries.....	191
2.1.5 Fuse	191
2.1.6 Protective caps	192
2.2 Overview of function keys.....	193
2.2.1 Overview.....	193
2.2.2 Displaying time and date	193
2.2.3 User-countdown.....	194
2.3 Operation mode	195
2.3.1 Automatic switch off.....	195
2.3.2 Selecting a method	195
2.3.2.1 Method-information.....	195

2.3.2.2	Chemical Species Information.....	196
2.3.3	Differentiation.....	196
2.3.4	Performing Zero.....	196
2.3.5	Performing Test.....	197
2.3.6	Ensuring reaction periods (countdown).....	197
2.3.7	Changing chemical species.....	198
2.3.8	Storing results.....	198
2.3.9	Printing results.....	199
2.3.10	Perform additional measurements.....	199
2.3.11	Selecting a new method.....	200
2.3.12	Measure absorbance.....	200
2.4	Photometer settings <Mode-Menu>.....	201
2.4.1	blank because of technical requirements	
2.4.2	Instrument basic settings 1.....	202
2.4.3	Printing of stored results.....	206
2.4.4	Recall / delete stored results.....	211
2.4.5	Calibration.....	215
2.4.6	Lab function.....	221
2.4.7	User operations.....	222
2.4.8	Special functions.....	232
2.4.9	Instrument basic settings 2.....	234
2.4.10	Instrument special functions / Service.....	234
2.5	Data transfer.....	235
2.5.1	Connection to a printer.....	235
2.5.2	Data transfer to a personal computer (PC).....	235
2.5.3	Internet-Updates.....	235
2.6	blank because of technical requirements	
Part 3	Enclosure.....	237
3.1	Unpacking.....	238
3.2	Delivery content.....	238
3.3	blank because of technical requirements	
3.4	Technical data.....	239
3.5	Abbreviations.....	240
3.6	Troubleshooting.....	241
3.6.1	Operating messages in the display / error display.....	241
3.6.2	General problems.....	243
3.7	Declaration of CE-Conformity.....	244

Part 1

Methods

1.1 Methods

1.1 Table of Methods

No.	Analysis	Reagent	Range	Displayed as	Method	λ [nm]	Page
20	Acid demand to pH 4.3 T	tablet	0.1-4	mmol/l	Acid/Indicator ^{1,2,5}	610	12
30	Alkalinity, total T	tablet	5-200	mg/l CaCO ₃	Acid/Indicator ^{1,2,5}	610	14
35	Alkalinity-p T	tablet	5-500	mg/l CaCO ₃	Acid/Indicator ^{1,2,5}	560	16
40	Aluminium T	tablet	0.01-0.3	mg/l Al	Eriochrome Cyanine R ²	530	18
50	Aluminium PP	PP + liquid	0.01-0.25	mg/l Al	Eriochrome Cyanine R ²	530	20
60	Ammonium T	tablet	0.02-1	mg/l N	Indophenol blue ^{2,3}	610	22
62	Ammonium PP	PP	0.01-0.8	mg/l N	Salicylate ²	660	24
65	Ammonium LR TT	tube test	0.02-2.5	mg/l N	Salicylate ²	660	26
66	Ammonium HR TT	tube test	1-50	mg/l N	Salicylate ²	660	28
85	Boron T	tablet	0.1-2	mg/l B	Azomethine ³	430	30
80	Bromine T	tablet	0.05-13	mg/l Br ₂	DPD ⁵	530	32
90	Chloride T	tablet	0.5 -25	mg/l Cl	Silver nitrate/ turbidity	530	34
100	Chlorine T *	tablet	0.01-6	mg/l Cl ₂	DPD ^{1,2,3}	530	36, 38
101	Chlorine L *	liquid	0.02-4	mg/l Cl ₂	DPD ^{1,2,3}	530	36, 42
110	Chlorine PP *	PP	0.02-2	mg/l Cl ₂	DPD ^{1,2}	530	36, 46
120	Chlorine dioxide T	tablet	0.05-11	mg/l ClO ₂	DPD, Glycine ^{1,2}	530	50
105	Chlorine HR (KI) T	tablet	5-200	mg/l Cl ₂	KI/Acid ⁵	530	56
130	COD LR TT	tube test	0 -150	mg/l O ₂	Dichromate/H ₂ SO ₄ ^{1,2}	430	58
131	COD MR TT	tube test	0 -1500	mg/l O ₂	Dichromate/H ₂ SO ₄ ^{1,2}	610	60
132	COD HR TT	tube test	0 -15	g/l O ₂	Dichromate/H ₂ SO ₄ ^{1,2}	610	62
150	Copper T *	tablet	0.05-5	mg/l Cu	Biquinoline ⁴	560	64
153	Copper PP	PP	0.05-5	mg/l Cu	Bicinchoninate	560	68
157	Cyanide	PP + liquid	0.01-0.5	mg/l CN	Pyridine- barbituric acid ¹	580	70
160	Cyanuric acid T	tablet	2-160	mg/l Cys	Melamine	530	72
165	DEHA T	tablet + liquid	20-500	μ g/l DEHA	PPST ³	560	74
167	DEHA PP	PP + liquid	20-500	μ g/l DEHA	PPST ³	560	76
170	Fluoride L	liquid	0.05-2	mg/l F	SPADNS ²	580	78
190	Hardness, Calcium T	tablet	50-900	mg/l CaCO ₃	Murexide ⁴	560	80

* = free, combined, total; PP = powder pack; T = tablet;
L = liquid; TT = tube test; LR = low range; MR = middle range; HR = high range;
Vacu-vial® is a registered trade mark of CHEMetrics Inc.

1.1 Methods

1.1 Table of Methods

No.	Analysis	Reagent	Range	Displayed as	Method	λ [nm]	Page
200	Hardness, total T	tablet	2-50	mg/l CaCO ₃	Metallphthalein ³	560	82
201	Hardness, total HR T	tablet	20-500	mg/l CaCO ₃	Metallphthalein ³	560	84
205	Hydrazine P	powder	0.05-0.5	mg/l N ₂ H ₄	4-(Dimethyl-amino)-benzaldehyde ³	430	86
206	Hydrazine L	liquid	0.01-0.6	mg/l N ₂ H ₄	4-(Dimethyl-amino)-benzaldehyde ³	430	88
207	Hydrazine C	Vacu-vial	0.01-0.7	mg/l N ₂ H ₄	PDMAB	430	90
210	Hydrogen peroxide	tablet	0.03-3	mg/l H ₂ O ₂	DPD/catalyst ⁵	530	92
215	Iodine T	tablet	0.05-3.6	mg/l I	DPD ⁵	530	94
220	Iron T	tablet	0.02-1	mg/l Fe	PPST ³	560	96, 98
222	Iron PP	PP	0.02-3	mg/l Fe	1,10-Phenanthroline ³	530	96, 100
223	Iron (TPTZ) PP	PP	0.02-1.8	mg/l Fe	TPTZ	580	96, 102
240	Manganese T	tablet	0.2-4	mg/l Mn	Formaloxime	530	104
242	Manganese LR PP	PP + liquid	0.01-0.7	mg/l Mn	PAN	560	106
243	Manganese HR PP	PP + liquid	0,1-18	mg/l Mn	Periodate oxidation ²	530	108
250	Molybdate T	tablet	1-50	mg/l MoO ₄	Thioglycolate ⁴	430	110
252	Molybdate HR PP	PP	0.5-66	mg/l MoO ₄	Mercaptoacetic acid	430	112
265	Nitrate TT	tube test	1-30	mg/l N	Chromotropic acid	430	114
270	Nitrite T	tablet	0.01-0.5	mg/l N	N-(1-Naphthyl)-ethylenediamine ^{2,3}	560	116
272	Nitrite LR PP	PP	0.01-0.3	mg/l N	Diazotization	530	118
280	Nitrogen, total LR TT	tube test	0.5-25	mg/l N	Persulfate digestion method	430	120
281	Nitrogen, total HR TT	tube test	5-150	mg/l N	Persulfate digestion method	430	122
290	Oxygen, active T	tablet	0.1-10	mg/l O ₂	DPD	530	124
292	Oxygen, dissolved	Vacu-vial	10-800	μ g/l O ₂	Rhodazine D TM	530	126
300	Ozone (DPD) T	tablet	0.02-1	mg/l O ₃	DPD/Glycine ⁵	530	128
70	PHMB T	tablet	2-60	mg/l PHMB	Buffer/Indicator	560	134

* = free, combined, total; PP = powder pack; T = tablet;
 L = liquid; TT = tube test; LR = low range; MR = middle range; HR = high range;
 Vacu-vial[®] is a registered trade mark of CHEMetrics Inc.

1.1 Methods

1.1 Table of Methods

No.	Analysis	Reagent	Range	Displayed as	Method	λ [nm]	Page
320	Phosphate, T ortho LR	tablet	0.05-4	mg/l PO ₄	Ammonium-molybdate ^{2,3}	660	136, 138
321	Phosphate, ortho HR T	tablet	1-80	mg/l PO ₄	Vanado-molybdate ²	430	136, 140
323	Phosphate, PP ortho	PP	0.06-2.5	mg/l PO ₄	Ascorbic acid ²	660	136, 142
324	Phosphate, ortho TT	tube test	0.06-5	mg/l PO ₄	Ascorbic acid ²	660	136, 144
327	Phosphate 1 C, ortho	Vacu-vial	5-40	mg/l PO ₄	Vanado-molybdate ²	430	136, 146
328	Phosphate 2 C, ortho	Vacu-vial	0.05-5	mg/l PO ₄	Stannous chloride ²	660	136, 148
325	Phosphate, hydr. TT	tube test	0.02-1.6	mg/l P	Acid digestion, Ascorbic acid ²	660	136, 150
326	Phosphate, total TT	tube test	0.02-1.1	mg/l P	Acid persulf digestion, Ascorbic acid ²	660	136, 152
329	pH-Value LR T	tablet	5.2-6.8	—	Bromocresolpurple ⁵	560	154
330	pH-Value T	tablet	6.5-8.4	—	Phenolred ⁵	560	156
331	pH-Value L	liquid	6.5-8.4	—	Phenolred ⁵	560	158
332	pH-Value HR T	tablet	8.0-9.6	—	Thymolblue ⁵	560	160
340	Potassium T	tablet	0.7-12	mg/l K	Tetraphenylborate-Turbidity ⁴	430	162
350	Silica T	tablet	0.05-4	mg/l SiO ₂	Silicomolybdate ^{2,3}	660	164
351	Silica LR PP	PP	0.1-1.6	mg/l SiO ₂	Heteropolyblue ²	660	166
352	Silica HR PP	PP	1-90	mg/l SiO ₂	Silicomolybdate ²	430	168
212	Sodium hypochlorite T	tablet	0.2-16	% NaOCl	Potassium iodide ⁵	530	170
355	Sulfate T	tablet	5-100	mg/l SO ₄	Bariumsulfate-Turbidity	610	172
360	Sulfate PP	PP	5-100	mg/l SO ₄	Bariumsulfate-Turbidity ²	530	174
365	Sulfide	tablet	0.04-0.5	mg/l S	DPD/Catalyst ^{3,4}	660	176
370	Sulfite T	tablet	0.1-5	mg/l SO ₃	DTNB	430	178
390	Urea T	tablet + liquid	0.1-3	mg/l Urea	Indophenol/Urease	610	180
400	Zinc T	tablet	0.02 -1	mg/l Zn	Zincon ³	610	182

* = free, combined, total; PP = powder pack; T = tablet;
 L = liquid; TT = tube test; LR = low range; MR = middle range; HR = high range;
 Vacu-vial® is a registered trade mark of CHEMetrics Inc.

1.1 Methods

Literature

The reagent formulations are based on internationally recognised test methods. Some are described in national and/or international guidelines.

- 1) Deutsche Einheitsverfahren zur Wasser-, Abwasser- und Schlammuntersuchung
- 2) Standard Methods for the Examination of Water and Wastewater; 18th Edition, 1992
- 3) Photometrische Analysenverfahren, Schwedt, Wissenschaftliche Verlagsgesellschaft mbH, Stuttgart 1989
- 4) Photometrische Analyse, Lange / Vejdelek, Verlag Chemie 1980
- 5) Colorimetric Chemical Analytical Methods, 9th Edition, London

Notes for searching:

Active Oxygen	->	Oxygen, activ
Alkalinity-m	->	Alkalinity, total
Alkalinity, total	->	Alkalinity, total
Biguanide	->	PHMB
Calcium Hardness	->	Hardness, Calcium
Total Hardness	->	Hardness, total
m-Value	->	Alkalinity, total
p-Value	->	Alkalinity-p
Silicon dioxide	->	Silica
total Alkalinity	->	Alkalinity, total
total Hardness	->	Hardness, total

Langelier Saturation Index (Water Balance)	->	Mode function 70
---	----	-------------------------

1.1 Methods

2

0

Acid demand to pH 4.3 with Tablet

0.1 – 4 mmol/l



1. Fill a clean vial (24 mm Ø) with **10 ml of water sample**, close the vial with the cap tightly.

2. Place the vial in the sample chamber making sure that the \times marks are aligned.

prepare Zero
press ZERO

3. Press **ZERO** key.

4. Remove the vial from the sample chamber.

5. Add **one ALKA-M-PHOTOMETER tablet** straight from the foil to the water sample and crush the tablet using a clean stirring rod.

6. Close the vial with the cap tightly and swirl the vial several times until the tablet is dissolved.

7. Place the vial in the sample chamber making sure that the \times marks are aligned.

Zero accepted
prepare Test
press TEST

8. Press **TEST** key.

The result is shown in the display as Acid demand to pH 4.3 in mmol/l.

1.1 Methods

Notes:

1. The terms total Alkalinity, Alkalinity-m, m-Value and Acid demand to pH 4.3 are identical.
2. For accurate results exactly 10 ml of water sample must be taken for the test.

1.1 Methods

3

0

Alkalinity, total = Alkalinity-m = m-Value with Tablet

5 – 200 mg/l CaCO₃



prepare Zero
press ZERO

1. Fill a clean vial (24 mm Ø) with **10 ml of water sample**, close the vial with the cap tightly.
2. Place the vial in the sample chamber making sure that the Σ marks are aligned.
3. Press **ZERO** key.
4. Remove the vial from the sample chamber.
5. Add **one ALKA-M-PHOTOMETER tablet** straight from the foil to the water sample and crush the tablet using a clean stirring rod.
6. Close the vial with the cap tightly and swirl the vial several times until the tablet is dissolved.
7. Place the vial in the sample chamber making sure that the Σ marks are aligned.
8. Press **TEST** key.

Zero accepted
prepare Test
press TEST

The result is shown in the display as total Alkalinity.

1.1 Methods

Notes:

1. The terms total Alkalinity, Alkalinity-m, m-Value and Alkalinity to pH 4.3 are identical.
2. For accurate results exactly 10 ml of water sample must be taken for the test.
3. Conversion table:

	Acid demand to pH 4.3 DIN 38 409 (K _s 4.3)	German °dH*	English °eH*	French °fH*
1 mg/l CaCO ₃	0.02	0.056	0.07	0.1

*Carbonate hardness (reference = Hydrogencarbonate-anions)

Example:

$$10 \text{ mg/l CaCO}_3 = 10 \text{ mg/l} \times 0.056 = 0.56 \text{ mg/l } ^\circ\text{dH}$$

$$10 \text{ mg/l CaCO}_3 = 10 \text{ mg/l} \times 0.02 = 0.2 \text{ mmol/l}$$

4. ▲ CaCO₃
°dH
°eH
°fH
▼ °aH

1.1 Methods

3

5

Alkalinity-p = p-value with Tablet

5 – 500 mg/l CaCO₃



prepare Zero
press ZERO

1. Fill a clean vial (24 mm Ø) with **10 ml of water sample**, close the vial with the cap tightly.
2. Place the vial in the sample chamber making sure that the \times marks are aligned.
3. Press **ZERO** key.
4. Remove the vial from the sample chamber.
5. Add **one ALKA-P-PHOTOMETER tablet** straight from the foil to the water sample and crush the tablet using a clean stirring rod.
6. Close the vial with the cap tightly and swirl the vial several times until the tablet is dissolved.
7. Place the vial in the sample chamber making sure that the \times marks are aligned.

Zero accepted
prepare Test
press TEST

8. Press **TEST** key.

The result is shown in the display as Alkalinity-p.

1.1 Methods

Notes

1. The terms Alkalinity-p, p-Value and Alkalinity to pH 8.2 are identical.
2. For accurate test results exactly 10 ml of water sample must be taken for the test.
3. This method was developed from a volumetric procedure for the determination of Alkalinity-p. Due to undefined conditions, the deviations from the standardised method may be greater.
4. Conversion table:

	mg/l CaCO ₃	°dH	°fH	°eH
1 mg/l CaCO ₃	----	0,056	0,10	0,07
1 °dH	17,8	----	1,78	1,25
1 °fH	10,0	0,56	----	0,70
1 °eH	14,3	0,80	1,43	----

▲ CaCO₃
°dH
°eH
°fH
▼ °aH

5. By determining Alkalinity-p and Alkalinity-m it is possible to classify the alkalinity as Hydroxide, Carbonate and Hydrogencarbonate.

The following differentiation is only valid if:

- a) no other alkalis are present and
- b) Hydroxide und Hydrogen are not present in the same water sample.

If condition b) is not fulfilled please get additional information from "Deutsche Einheitsverfahren zur Wasser-, Abwasser- und Schlammuntersuchung, D 8".

Case 1: Alkalinity-p = 0

Hydrogen carbonate = m

Carbonate = 0

Hydroxide = 0

Case 2: Alkalinity-p > 0 and Alkalinity-m > 2p

Hydrogen carbonate = m – 2p

Carbonate = 2p

Hydroxide = 0

Case 3: Alkalinity-p > 0 and Alkalinity-m < 2p

Hydrogen carbonate = 0

Carbonate = 2m – 2p

Hydroxide = 2p – m

1.1 Methods

4

0

Aluminium with Tablet

0.01 – 0.3 mg/l Al



1. Fill a clean vial (24 mm Ø) with **10 ml of water sample**, close the vial with the cap tightly.

2. Place the vial in the sample chamber making sure that the Σ marks are aligned.

prepare Zero
press ZERO

3. Press **ZERO** key.

4. Remove the vial from the sample chamber.

5. Add **one ALUMINIUM No. 1 tablet** straight from the foil to the water sample and crush the tablet using a clean stirring rod (dissolve the tablet).

6. Add **one ALUMINIUM No. 2 tablet** straight from the foil to the same water sample and crush the tablet using a clean stirring rod.

7. Close the vial with the cap tightly and swirl the vial gently several times until the tablets are dissolved.

8. Place the vial in the sample chamber making sure that the Σ marks are aligned.

Zero accepted
prepare Test
press TEST

9. Press **TEST** key.
Wait for a **reaction period of 5 minutes**.

Count-Down
5:00

After the reaction period is finished the reading starts automatically.

The result is shown in the display in mg/l Aluminium.

1.1 Methods

Notes:

1. Before using clean the vials and the measuring beaker with Hydrochloric acid (approx. 20%). Rinse then thoroughly with deionized water.
2. To get accurate results the sample temperature must be between 20°C and 25°C.
3. A low test result may be given in the presence of Fluorides and Polyphosphates. The effect of this is generally insignificant unless the water has fluoride added artificially. In this case, the following table should be used:

Fluoride [mg/l F]	Displayed value: Aluminium [mg/l Al]					
	0,05	0,10	0,15	0,20	0,25	0,30
0,2	0,05	0,11	0,16	0,21	0,27	0,32
0,4	0,06	0,11	0,17	0,23	0,28	0,34
0,6	0,06	0,12	0,18	0,24	0,30	0,37
0,8	0,06	0,13	0,20	0,26	0,32	0,40
1,0	0,07	0,13	0,21	0,28	0,36	0,45
1,5	0,09	0,20	0,29	0,37	0,48	---

Example: If the result of Aluminium determination is 0.15 mg/l Al and the Fluoride concentration is known to be 0.4 mg/l F, the true concentration of Aluminium is 0.17 mg/l Al.

4. ▲ Al
▼ Al₂O₃

1.1 Methods

5

0

Aluminium with Vario Powder Pack

0.01 – 0.25 mg/l Al




Ø 24 mm

Use two clean vials (24 mm Ø) and mark one as blank for zeroing.

1. Fill **20 ml of water sample** in a 100 ml beaker.
2. Add **one Vario Aluminum ECR F20 powder pack** straight from the foil to the water sample.
3. Dissolve the powder using a clean stirring rod.
4. Press [**↵**] key.
Wait for a **reaction period of 30 seconds**.

Countdown 1
0:30
start: ↵

After reaction period is finished proceed as follows:



5. Add **one Vario Hexamine F20 powder pack** straight from the foil to the same water sample.
6. Dissolve the powder using a clean stirring rod.
7. Add **1 drop of Vario Aluminum ECR Masking Reagent** in the vial marked as blank.
8. Add 10 ml of the prepared water sample to the vial (this is the blank).
9. Add the remaining 10 ml of the prepared water sample in the second clean vial (this is the sample).
10. Close the vials with the caps tightly and swirl the vials several times to mix the contents.

Countdown 2
5:00
start: ↵

11. Press [**↵**] key.
Wait for a **reaction period of 5 minutes**.

1.1 Methods

After reaction period is finished proceed as follows:

12. Place the vial (the blank) in the sample chamber making sure that the Σ marks are aligned.

**prepare Zero
press ZERO**

13. Press **ZERO** key.

14. Remove the vial from the sample chamber.

15. Place the vial (the sample) in the sample chamber making sure that the Σ marks are aligned.

**Zero accepted
prepare Test
press TEST**

16. Press **TEST** key.

The result is shown in the display in mg/l Aluminium.

Notes:

1. Before using clean the vials and the measuring beaker with Hydrochloric acid (approx. 20%). Rinse then thoroughly with deionized water.
2. To get accurate results the sample temperature must be between 20°C and 25°C.
3. A low test result may be given in the presence of Fluorides and Polyphosphates. The effect of this is generally insignificant unless the water has fluoride added artificially. In this case, the following table should be used:

Fluoride [mg/l F]	Displayed value: Aluminium [mg/l Al]					
	0,05	0,10	0,15	0,20	0,25	0,30
0,2	0,05	0,11	0,16	0,21	0,27	0,32
0,4	0,06	0,11	0,17	0,23	0,28	0,34
0,6	0,06	0,12	0,18	0,24	0,30	0,37
0,8	0,06	0,13	0,20	0,26	0,32	0,40
1,0	0,07	0,13	0,21	0,28	0,36	0,45
1,5	0,09	0,20	0,29	0,37	0,48	---

Example: If the result of Aluminium determination is 0.15 mg/l Al and the Fluoride concentration is known to be 0.4 mg/l F, the true concentration of Aluminium is 0.17 mg/l Al.

4. ▲ Al
▼ Al₂O₃


1.1 Methods

6

0

Ammonium with Tablet

0.02 – 1 mg/l N



prepare Zero
press ZERO

1. Fill a clean vial (24 mm Ø) with **10 ml of water sample**, close the vial with the cap tightly.
2. Place the vial in the sample chamber making sure that the Σ marks are aligned.
3. Press **ZERO** key.
4. Remove the vial from the sample chamber.
5. Add **one AMMONIA No. 1 tablet** straight from the foil to the water sample and crush the tablet using a clean stirring rod.
6. Add **one AMMONIA No. 2 tablet** straight from the foil to the same water sample and crush the tablet using a clean stirring rod.
7. Close the vial with the cap tightly and swirl the vial several times until the tablets are dissolved.
8. Place the vial in the sample chamber making sure that the Σ marks are aligned.
9. Press **TEST** key.
Wait for a **reaction period of 10 minutes**.

Zero accepted
prepare Test
press TEST

Countdown
10:00

After the reaction period is finished the reading starts automatically.

The result is shown in the display in mg/l Ammonium.

1.1 Methods

Notes:

1. The tablets must be added in the correct sequence.
2. The AMMONIA No. 1 tablet will only dissolve completely after the AMMONIA No. 2 tablet has been added.
3. The temperature of the sample is important for full colour development.
At a temperature below 20°C the reaction period is 15 minutes.
4. Sea water samples:
Ammonia conditioning reagent is required when testing sea water or brackish water samples to prevent precipitations of salts.
Fill the test tube with the sample to the 10 ml mark and add one level spoonful of Conditioning Powder. Mix to dissolve, then continue as described in the test instructions.
5. Conversion:
 $\text{mg/l NH}_4 = \text{mg/l N} \times 1.29$
 $\text{mg/l NH}_3 = \text{mg/l N} \times 1.22$
6. ▲ N
 NH₄
 ▼ NH₃


1.1 Methods

6

2

Ammonium mit Vario Powder Pack

0.01 – 0.8 mg/l N



Countdown 1
3:00
start: ↵

Countdown 2
15:00
start: ↵

prepare Zero
press ZERO

Zero accepted
prepare Test
press TEST

Use two clean vials (24 mm Ø) and mark one as blank for zeroing.

1. Fill a clean vial (24 mm Ø) with **10 ml of deionised water** (this is the blank).
2. Fill the other clean vial (24 mm Ø) with **10 ml of water sample** (this is the sample).
3. Add **one Vario Ammonium Salicylate F10 powder pack** straight from the foil to each vial.

4. Close the vials with the caps and shake to mix the contents.

5. Press [↵] key.
Wait for a **reaction period of 3 minutes**.

After reaction period is finished proceed as follows:

6. Add **one Vario Ammonium Cyanurate F10 powder pack** straight from the foil to each sample.

7. Close the vials with the caps tightly and shake to mix the contents.

8. Press [↵] key.
Wait for a **reaction period of 15 minutes**.

After reaction period is finished proceed as follows:

9. Place the vial (the blank) in the sample chamber making sure that the Σ marks are aligned.

10. Press **ZERO** key.

11. Remove the vial from the sample chamber.

12. Place the vial (the sample) in the sample chamber making sure that the Σ marks are aligned.

13. Press **TEST** key.

The result is shown in the display in mg/l Ammonium.

1.1 Methods

Notes:

1. Extremely basic or acidic water samples should be adjusted with 0.5 mol/l (1 N) Sulfuric acid solution or 1 mol/l (1 N) Sodium hydroxide solution to pH 7.
2. Interferences:

Interfering substance	Interference levels and treatments
Calcium	greater than 1000 mg/l CaCO ₃
Iron	Interferes at all levels. Correct as follows: a) determine the concentration of iron present in the sample by performing a total Iron test b) add the same iron concentration as determined to the deionised water (step 1). The interference will be blanked out successfully.
Magnesium	greater than 6000 mg/l CaCO ₃
Nitrate	greater than 100 mg/l NO ₃ -N
Nitrite	greater than 12 mg/l NO ₂ -N
Phosphate	greater than 100 mg/l PO ₄ -P
Sulfate	greater than 300 mg/l SO ₄
Sulfide	intensifies the colour
Glycine, Hydrazine, Colour, Turbidity	Less common interferences such as Hydrazine and Glycine will cause intensified colours in the prepared sample. Turbidity and colour will give erroneous high values. Samples with severe interferences require distillation.

3. ▲ N
NH₄
▼ NH₃


1.1 Methods

6


5

Ammonium LR with Vario Tube Test

0.02 – 2.5 mg/l N



Ø 16 mm



1. Open the white cap of one reaction vial and add **2 ml deionised water** (this is the blank).
2. Open the white cap of another reaction vial and add **2 ml water sample** (this is the sample).
3. Add **one Vario AMMONIA Salicylate F5 powder pack** straight from the foil into each vial.
4. Add **one Vario AMMONIA Cyanurate F5 powder pack** straight from the foil into each vial.
5. Close the vials with the caps tightly and swirl the vials several times to dissolve the powder.
6. Press **[↓]** key.
Wait for a **reaction period of 20 minutes**.

Countdown 1
20:00
start: ↵

After reaction period is finished proceed as follows:

7. Place the vial (the blank) in the sample chamber making sure that the marks are **Δ** aligned.
Place the cover on the adapter.
8. Press **ZERO** key.
9. Remove the vial from the sample chamber.
10. Place the vial (the sample) in the sample chamber making sure that the marks are **Δ** aligned.
Place the cover on the adapter.
11. Press **TEST** key.

prepare Zero
press ZERO

Zero accepted
prepare Test
press TEST

The result is shown in the display in mg/l Ammonium.

1.1 Methods

Notes:

1. Strong alkaline or acidic water samples must be adjusted to approx. pH 7 before analysis (use 1 mol/l Hydrochloric acid resp. 1 mol/l Sodium hydroxide).
2. If chlorine is known to be present, add one drop of 0.1 mol/l Sodium thiosulfate for each 0.3 mg/l Cl_2 in a one litre water sample.
3. Iron interferes with the test. The interferences will be eliminated as follows:
Determine the amount of total iron present in the water sample. To produce the blank add an iron standard solution with the same iron concentration to the vial (point 1) instead of deionised water
4. Conversion:
 $\text{mg/l NH}_4 = \text{mg/l N} \times 1.29$
 $\text{mg/l NH}_3 = \text{mg/l N} \times 1.22$
5. ▲ N
 NH₄
 ▼ NH₃



1.1 Methods

6

6

Ammonium HR with Vario Tube Test

1 – 50 mg/l N



1. Open the white cap of one reaction vial and add **0.1 ml deionised water** (this is the blank).
2. Open the white cap of another reaction vial and add **0.1 ml water sample** (this is the sample).
3. Add **one Vario AMMONIA Salicylate F5 powder pack** straight from the foil into each vial.
4. Add **one Vario AMMONIA Cyanurate F5 powder pack** straight from the foil into each vial.
5. Close the vials with the caps tightly and swirl the vials several times to dissolve the powder.

Countdown 1
20:00
start: ↵

6. Press **[↵]** key.
Wait for a **reaction period of 20 minutes**.

After reaction period is finished proceed as follows:

7. Place the vial (the blank) in the sample chamber making sure that the marks are **Δ** aligned.
Place the cover on the adapter.

prepare Zero
press ZERO

8. Press **ZERO** key.

9. Remove the vial from the sample chamber.

10. Place the vial (the sample) in the sample chamber making sure that the marks are **Δ** aligned.
Place the cover on the adapter.

Zero accepted
prepare Test
press TEST

11. Press **TEST** key.

The result is shown in the display in mg/l Ammonium.

1.1 Methods

Notes:

1. Strong alkaline or acidic water samples must be adjusted to approx. pH 7 before analysis (use 1 mol/l Hydrochloric acid resp. 1 mol/l Sodium hydroxide).
2. If chlorine is known to be present, add one drop of 0.1 mol/l Sodium thiosulfate for each 0.3 mg/l Cl₂ in a one litre water sample.
3. Iron interferes with the test. The interferences will be eliminated as follows:
Determine the amount of total iron present in the water sample. Add an iron standard solution with the same concentration to the vial (point 1) instead of deionised water to produce the blank.
4. Conversion:
 $\text{mg/l NH}_4 = \text{mg/l N} \times 1.29$
 $\text{mg/l NH}_3 = \text{mg/l N} \times 1.22$
5. ▲ N
 NH₄
 ▼ NH₃


1.1 Methods

8

5

Boron with Tablet

0.1 – 2 mg/l B



prepare Zero
press ZERO

1. Fill a clean vial (24 mm Ø) with **10 ml of water sample**, close the vial with the cap tightly.
2. Place the vial in the sample chamber making sure that the Σ marks are aligned.
3. Press **ZERO** key.
4. Remove the vial from the sample chamber.
5. Add **one BORON No. 1 tablet** straight from the foil to the water sample and crush the tablet using a clean stirring rod and dissolve the tablet.
6. Add **one BORON No. 2 tablet** straight from the foil to the same water sample and crush the tablet using a clean stirring rod.
7. Close the vial with the cap tightly and swirl the vial several times until the tablets are dissolved.
8. Place the vial in the sample chamber making sure that the Σ marks are aligned.

Zero accepted
prepare Test
press TEST

9. Press **TEST** key.

Countdown
20:00

Wait for a **reaction period of 20 minutes**.

After the reaction period is finished the reading starts automatically.

The result is shown in the display in mg/l Boron.

1.1 Methods

Notes:

1. The tablets must be added in the correct sequence.
2. The sample solution should have a pH value between 6 and 7.
3. Interferences are prevented by the presence of EDTA in the tablets.
4. The rate of colour development depends on the temperature. The temperature of the sample must be $20^{\circ}\text{C} \pm 1^{\circ}\text{C}$.
5. ▲ B
▼ H_3BO_3


1.1 Methods

8

0

Bromine with Tablet

0.05 – 13 mg/l Br₂



prepare Zero
press ZERO

1. Fill a clean vial (24 mm Ø) with **10 ml of water sample**, close the vial with the cap tightly.
2. Place the vial in the sample chamber making sure that the X marks are aligned.
3. Press **ZERO** key.
4. Remove the vial from the sample chamber, **empty the vial leaving a few drops in**.
5. Add **one DPD No. 1 tablet** straight from the foil and crush the tablet using a clean stirring rod.
6. Add water sample to the 10 ml mark.
7. Close the vial with the cap tightly and swirl the vial several times until the tablet is dissolved.
8. Place the vial in the sample chamber making sure that the X marks are aligned.

Zero accepted
prepare Test
press TEST

9. Press **TEST** key.

The result is shown in the display in mg/l Bromine.

1.1 Methods

Notes:

1. Vial cleaning:

As many household cleaners (e.g. dishwasher detergent) contain reducing substances, the subsequent determination of Bromine may show lower results. To avoid any measurement errors, only use glassware free of Chlorine consumption.

Preparation: Put all applicable glassware into Sodium hypochlorite solution (0.1 g/l) for one hour, then rinse all glassware thoroughly with deionized water.

2. Preparing the sample:

When preparing the sample, the escape of Bromine gases, e.g. by pipetting or shaking, must be avoided. The analysis must take place immediately after taking the sample.

3. The DPD colour development is carried out at a pH value of 6.3 to 6.5. The reagent tablet therefore contains a buffer for the pH adjustment. Strong alkaline or acidic water samples must be adjusted between pH 6 and pH 7 before the reagent is added (use 0.5 mol/l Sulfuric acid resp. 1 mol/l Sodium hydroxide).

4. Exceeding of the measuring range:

Concentrations above 22 mg/l Bromine can produce results within the measuring range up to 0 mg/l. In this event, the water sample must be diluted with water free of Bromine. 10 ml of the diluted sample will be mixed with the reagent and the measurement repeated.

Oxidizing agents such as Chlorine, Ozone etc. interfere as they react like Bromine.


1.1 Methods

9

0

Chloride with Tablet

0.5 – 25 mg/l Cl



prepare Zero
press ZERO

1. Fill a clean vial (24 mm Ø) with **10 ml of water sample**, close the vial with the cap tightly.
2. Place the vial in the sample chamber making sure that the X marks are aligned.
3. Press **ZERO** key.
4. Remove the vial from the sample chamber.
5. Add **one CHLORIDE T1 tablet** straight from the foil to the water sample, crush the tablet using a clean stirring rod and dissolve the tablet.
6. Add **one CHLORIDE T2 tablet** straight from the foil to the same water sample and crush the tablet using a clean stirring rod.
7. Close the vial with the cap tightly and swirl the vial gently several times until the tablet is dissolved (Note 1).
8. Place the vial in the sample chamber making sure that the X marks are aligned.

Zero accepted
prepare Test
press TEST

Count-Down
2:00

9. Press **TEST** key.
Wait for a **reaction period of 2 minutes**.
After the reaction period is finished the reading starts automatically.
The result is shown in the display in mg/l Chloride.

1.1 Methods

Notes:

1. Ensure that all particles of the tablet are dissolved – Chloride causes an extremely fine distributed turbidity with a milky appearance.
Heavy shaking leads to bigger sized particles which can cause false readings.
2. High concentrations of electrolytes and organic compounds have different effects on the precipitation reaction.
3. Ions which also form deposits with Silver nitrate in acidic media, such as Bromides, Iodides and Thiocyanates, interfere with the analysis.
4. Highly alkaline water should - if necessary - be neutralised using Nitric acid before analysis.

1.1 Methods

1 **0** **0**

Chlorine with Tablet

0.01 – 6 mg/l Cl₂

1 **0** **1**

Chlorine with Liquid Reagent

0.02 - 4 mg/l Cl₂

1 **1** **0**

Chlorine with Vario Powder Pack

0.02 - 2 mg/l Cl₂

Chlorine
>> **diff**
free
total

The following selection is shown in the display:

>> **diff**

for the differentiated determination of free, combined and total Chlorine.

>> **free**

for the determination of free Chlorine.

>> **total**

for the determination of total Chlorine.

Select the desired determination with the arrow keys [▲] and [▼]. Confirm with [↵] key.

1.1 Methods

Notes:

1. Vial cleaning:

As many household cleaners (e.g. dishwasher detergent) contain reducing substances, the subsequent determination of Chlorine may show lower results. To avoid any measurement errors, only use glassware free of Chlorine consumption.

Preparation: Put all applicable glassware into Sodium hypochlorite solution (0.1 g/l) for one hour, then rinse all glassware thoroughly with deionized water.
2. For individual testing of free and total Chlorine, the use of different sets of glassware is recommended (EN ISO 7393-2, 5.3)
3. Preparing the sample:

When preparing the sample, the escape of Chlorine gases, e.g. by pipetting or shaking, must be avoided. The analysis must take place immediately after taking the sample.
4. The DPD colour development is carried out at a pH value of 6.3 to 6.5. The reagents therefore contain a buffer for the pH adjustment.

Strong alkaline or acidic water samples must be adjusted between pH 6 and pH 7 before the reagent is added (use 0.5 mol/l Sulfuric acid resp. 1 mol/l Sodium hydroxide).
5. Exceeding of the measuring range:

Concentrations above

 - 10 mg/l Chlorine using tablets
 - 4 mg/l Chlorine using liquid reagents
 - 2 mg/l using powder packs

can produce results within the measuring range up to 0 mg/l. In this event, the water sample must be diluted with water free of Chlorine. 10 ml of the diluted sample will be mixed with the reagent and the measurement repeated.
6. Turbidity (lead to errors):

The use of the DPD No. 1 tablet (method 100) in samples with high Calcium ion content* and/or high conductivity* can lead to turbidity of the sample and therefore incorrect measurements. In this event, the reagent tablet DPD No. 1 High Calcium should be used as an alternative. Even if the turbidity does occur after the DPD No. 3 tablet has been added, this can be prevented by using the DPD No. 1 HIGH CALCIUM tablet.

** it is not possible to give exactly values, because the development of turbidity depends on nature and ingredients of the sample.*
7. If ??? is displayed at a differentiated test result see page 242.

Oxidizing agents such as Bromine, Ozone etc. interfere as they react like Chlorine.

1.1 Methods

1 0 0

Chlorine, differentiated determination with Tablet

0.01 – 6 mg/l Cl₂



prepare Zero
press ZERO

1. Fill a clean vial (24 mm Ø) with **10 ml of water sample**, close the vial with the cap tightly.
2. Place the vial in the sample chamber making sure that the X marks are aligned.
3. Press **ZERO** key.
4. Remove the vial from the sample chamber, **empty the vial leaving a few drops in**.
5. Add **one DPD No. 1 tablet** straight from the foil and crush the tablet using a clean stirring rod.
6. Add water sample to the 10 ml mark.
7. Close the vial with the cap tightly and swirl the vial several times until the tablet is dissolved.
8. Place the vial in the sample chamber making sure that the X marks are aligned.

Zero accepted
prepare T1
press TEST

9. Press **TEST** key.
10. Remove the vial from the sample chamber.
11. Add **one DPD No. 3 tablet** straight from the foil to the same water sample and crush the tablet using a clean stirring rod.
12. Close the vial with the cap tightly and swirl the vial several times until the tablet is dissolved.

1.1 Methods

13. Place the vial in the sample chamber making sure that the Σ marks are aligned.

T1 accepted
prepare T2
press TEST

14. Press **TEST** key.
Wait for a **reaction period of 2 minutes**.

Countdown
2:00

After the reaction period is finished the reading starts automatically.

*,** mg/l free Cl
*,** mg/l comb Cl
*,** mg/l total Cl

The result is shown in the display in:

mg/l free Chlorine
mg/l combined Chlorine
mg/l total Chlorine

Notes:

See page 37.

1.1 Methods

1 0 0

Chlorine, free with Tablet

0.01 – 6 mg/l Cl₂



prepare Zero
press ZERO

1. Fill a clean vial (24 mm Ø) with **10 ml of water sample**, close the vial with the cap tightly.
2. Place the vial in the sample chamber making sure that the X marks are aligned.
3. Press **ZERO** key.
4. Remove the vial from the sample chamber, **empty the vial leaving a few drops in**.
5. Add **one DPD No. 1 tablet** straight from the foil and crush the tablet using a clean stirring rod.
6. Add water sample to the 10 ml mark.
7. Close the vial with the cap tightly and swirl the vial several times until the tablet is dissolved.
8. Place the vial in the sample chamber making sure that the X marks are aligned.

Zero accepted
prepare Test
press TEST

9. Press **TEST** key.

The result is shown in the display in mg/l free Chlorine.

Notes:


See page 37.

1.1 Methods

1 0 0

Chlorine, total with Tablet

0.01 – 6 mg/l Cl₂



prepare Zero
press ZERO

1. Fill a clean vial (24 mm Ø) with **10 ml of water sample**, close the vial with the cap tightly.
2. Place the vial in the sample chamber making sure that the X marks are aligned.
3. Press **ZERO** key.
4. Remove the vial from the sample chamber, **empty the vial leaving a few drops in**.
5. Add **one DPD No. 1 tablet** and **one DPD No. 3 tablet** straight from the foil and crush the tablets using a clean stirring rod.
6. Add water sample to the 10 ml mark.
7. Close the vial with the cap tightly and swirl the vial several times until the tablets are dissolved.
8. Place the vial in the sample chamber making sure that the X marks are aligned.

Zero accepted
prepare Test
press TEST

Countdown
2:00

9. Press **TEST** key.
Wait for a **reaction period of 2 minutes**.

After the reaction period is finished the reading starts automatically.

The result is shown in the display in mg/l total Chlorine.

Notes:

See page 37.

1.1 Methods

1 0 1

Chlorine, differentiated determination with Liquid Reagent

0.02 – 4 mg/l Cl₂



prepare Zero
press ZERO

1. Fill a clean vial (24 mm Ø) with **10 ml of water sample**, close the vial with the cap tightly.
2. Place the vial in the sample chamber making sure that the X marks are aligned.

3. Press **ZERO** key.

4. Remove the vial from the sample chamber and **empty the vial**.

5. Fill the vial with drops of the same size by holding the bottle vertically and squeeze slowly:

6 drops of DPD 1 buffer solution

2 drops of DPD 1 reagent solution

6. Add water sample to the 10 ml mark.

7. Close the vial with the cap tightly and swirl the vial several times to mix the contents.

8. Place the vial in the sample chamber making sure that the X marks are aligned.

Zero accepted
prepare T1
press TEST

9. Press **TEST** key.

10. Remove the vial from the sample chamber.

11. **Add 3 drops of DPD 3 solution** to the same water sample.

12. Close the vial with the cap tightly and swirl the vial several times to mix the contents.

1.1 Methods

T1 accepted
prepare T2
press TEST

Countdown
2:00

*,** mg/l free Cl
*,** mg/l comb. Cl
*,** mg/l total Cl

13. Place the vial in the sample chamber making sure that the \times marks are aligned.

14. Press **TEST** key.
Wait for a **reaction period of 2 minutes**.

After the reaction period is finished the reading starts automatically.

The result is shown in the display in:

mg/l free Chlorine
mg/l combined Chlorine
mg/l total Chlorine

Notes:


1. After use replace the bottle caps securely noting the colour coding.
2. **Store the reagent bottles in a cool, dry place ideally between 6°C and 10°C.**
3. Also see page 37.

1.1 Methods

1 0 1

Chlorine, free with Liquid Reagent

0.02 – 4 mg/l Cl₂



prepare Zero
press ZERO

1. Fill a clean vial (24 mm Ø) with **10 ml of water sample**, close the vial with the cap tightly.
2. Place the vial in the sample chamber making sure that the X marks are aligned.
3. Press **ZERO** key.
4. Remove the vial from the sample chamber and **empty the vial**.
5. Fill the vial with drops of the same size by holding the bottle vertically and squeeze slowly:

6 drops of DPD 1 buffer solution
2 drops of DPD 1 reagent solution

6. Add water sample to the 10 ml mark.
7. Close the vial with the cap tightly and swirl the vial several times to mix the contents.
8. Place the vial in the sample chamber making sure that the X marks are aligned.
9. Press **TEST** key.

Zero accepted
prepare Test
press TEST

The result is shown in the display in mg/l free Chlorine.

Notes (free and total Chlorine):


1. After use replace the bottle caps securely noting the colour coding.
2. **Store the reagent bottles in a cool, dry place ideally between 6°C and 10°C.**
3. Also see page 37.

1.1 Methods

1 0 1

Chlorine, total with Liquid Reagent

0.02 – 4 mg/l Cl₂



prepare Zero
press ZERO

1. Fill a clean vial (24 mm Ø) with **10 ml of water sample**, close the vial with the cap tightly.

2. Place the vial in the sample chamber making sure that the X marks are aligned.

3. Press **ZERO** key.

4. Remove the vial from the sample chamber and **empty the vial**.

5. Fill the vial with drops of the same size by holding the bottle vertically and squeeze slowly:

6 drops of DPD 1 buffer solution

2 drops of DPD 1 reagent solution

3 drops of DPD 3 solution

6. Add water sample to the 10 ml mark.

7. Close the vial with the cap tightly and swirl the vial several times to mix the contents.

8. Place the vial in the sample chamber making sure that the X marks are aligned.

Zero accepted
prepare Test
press TEST

9. Press **TEST** key.
Wait for a **reaction period of 2 minutes**.

Countdown
2:00

After the reaction period is finished the reading starts automatically.

The result is shown in the display in mg/l total Chlorine.

1.1 Methods


1 1 0

Chlorine, differentiated determination with Vario Powder Pack

0.02 – 2 mg/l Cl₂



prepare Zero
press ZERO



Zero accepted
prepare T1
press TEST

1. Fill a clean vial (24 mm Ø) with **10 ml of water sample**, close the vial with the cap tightly.
2. Place the vial in the sample chamber making sure that the Σ marks are aligned.
3. Press **ZERO** key.
4. Remove the vial from the sample chamber.
5. Add **one VARIO Chlorine FREE-DPD/F10 powder pack** straight from the foil to the water sample.
6. Close the vial with the cap tightly and swirl the vial several times to mix the contents (approx. 20 seconds).
7. Place the vial in the sample chamber making sure that the Σ marks are aligned.
8. Press **TEST** key.
9. Remove the vial from the sample chamber, empty the vial, rinse vial and cap several times and fill the vial with **10 ml of water sample**.
10. Add **one VARIO Chlorine TOTAL-DPD / F10 powder pack** straight from the foil to the water sample.
11. Close the vial with the cap tightly and swirl the vial several times to mix the contents (approx. 20 seconds).

1.1 Methods

T1 accepted
prepare T2
press TEST

Countdown
3:00

*,** mg/l free Cl
*,** mg/l comb. Cl
*,** mg/l total Cl

12. Place the vial in the sample chamber making sure that the \times marks are aligned.

13. Press **TEST** key.
Wait for a **reaction period of 3 minutes**.

After the reaction period is finished the reading starts automatically.

The result is shown in the display in:

mg/l free Chlorine
mg/l combined Chlorine
mg/l total Chlorine

Notes:

See page 37.

1.1 Methods

1 1 0

Chlorine, free with Vario Powder Pack

0.02 – 2 mg/l Cl₂




1. Fill a clean vial (24 mm Ø) with **10 ml of water sample**, close the vial with the cap tightly.

2. Place the vial in the sample chamber making sure that the X marks are aligned.

**prepare Zero
press ZERO**

3. Press **ZERO** key.

4. Remove the vial from the sample chamber.



5. Add **one VARIO Chlorine FREE-DPD / F10 powder pack** straight from the foil to the water sample.

6. Close the vial with the cap tightly and swirl the vial several times to mix the contents (approx. 20 seconds).

7. Place the vial in the sample chamber making sure that the X marks are aligned.

**Zero accepted
prepare Test
press TEST**

8. Press **TEST** key.

The result is shown in the display in mg/l free Chlorine.

Notes:

See page 37.

1.1 Methods

1 1 0

Chlorine, total with Vario Powder Pack

0.02 – 2 mg/l Cl₂




1. Fill a clean vial (24 mm Ø) with **10 ml of water sample**, close the vial with the cap tightly.
2. Place the vial in the sample chamber making sure that the X marks are aligned.

prepare Zero
press ZERO

3. Press **ZERO** key.

4. Remove the vial from the sample chamber.



5. Add **one VARIO Chlorine TOTAL-DPD / F10 powder pack** straight from the foil to the water sample.

6. Close the vial with the cap tightly and swirl the vial several times to mix the contents (approx. 20 seconds).

7. Place the vial in the sample chamber making sure that the X marks are aligned.

Zero accepted
prepare Test
press TEST

8. Press **TEST** key.
Wait for a **reaction period of 3 minutes**.

Countdown
3:00

After the reaction period is finished the reading starts automatically.

The result is shown in the display in mg/l total Chlorine.

Notes:

See page 37.

1.1 Methods

1 2 0

Chlorine dioxide with Tablet

0.05 – 11 mg/l ClO₂

Chlorine dioxide
>> **with Cl**
without Cl

The following selection is shown in the display:

>> **with Cl**

for the determination of Chlorine dioxide in the presence of Chlorine.

>> **without Cl**

for the determination of Chlorine dioxide in the absence of Chlorine.

Select the desired determination with the arrow keys [▲] and [▼]. Confirm with [↵] key.

1.1 Methods

Notes:

1. Vial cleaning:

As many household cleaners (e.g. dishwasher detergent) contain reducing substances, the subsequent determination of Chlorine dioxide may show lower results. To avoid any measurement errors, only use glassware free of Chlorine consumption.

Preparation: Put all applicable glassware into Sodium hypochlorite solution (0.1 g/l) for one hour, then rinse all glassware thoroughly with deionized water.

2. Preparing the sample:

When preparing the sample, the escape of Chlorine dioxide gases, e.g. by pipetting or shaking, must be avoided. The analysis must take place immediately after taking the sample.

3. The DPD colour development is carried out at a pH value of 6.3 to 6.5. The reagent tablet therefore contains a buffer for the pH adjustment.

Strong alkaline or acidic water samples must be adjusted between pH 6 and pH 7 before the tablet is added (use 0.5 mol/l Sulfuric acid resp. 1 mol/l Sodium hydroxide).

4. Exceeding of the measuring range:

Concentrations above 19 mg/l Chlorine dioxide can produce results within the measuring range up to 0 mg/l. In this event, the water sample must be diluted with water free of Chlorine dioxide. 10 ml of the diluted sample will be mixed with the reagent and the measurement repeated.

5. If ??? is displayed at a differentiated test result see page 242.

Oxidizing agents such as Chlorine, Ozone etc. interfere as they react like Chlorine dioxide.

1.1 Methods

1 2 0

Chlorine dioxide in the presence of Chlorine with Tablet

0.05 – 11 mg/l ClO₂



prepare Zero
press ZERO

1. Fill a clean vial (24 mm Ø) with **10 ml of water sample**, close the vial with the cap tightly.
2. Place the vial in the sample chamber making sure that the X marks are aligned.
3. Press **ZERO** key.
4. Remove the vial from the sample chamber, **empty the vial leaving a few drops in.**
5. Add **one DPD No. 1 tablet** straight from the foil and crush the tablet using a clean stirring rod.
6. **Fill a second clean vial with 10 ml of water sample.**
7. Add **one GLYCINE tablet** straight from the foil and crush the tablet using a clean stirring rod.
8. Close the vial with the cap tightly and swirl the vial several times until the tablet is dissolved.
9. **Transfer the content of the second vial into the prepared vial.**
10. Close the vial with the cap tightly and swirl the vial several times until the tablet is dissolved.
11. Place the vial in the sample chamber making sure that the X marks are aligned.

Zero accepted
prepare T1
press TEST

12. Press **TEST** key.

1.1 Methods

13. Remove the vial from the sample chamber, empty the vial, rinse vial and cap several times. Fill with **a few drops of water sample**.
14. Add **one DPD No. 1 tablet** straight from the foil and crush the tablet using a clean stirring rod.
15. Add water sample to the 10 ml mark.
16. Close the vial with the cap tightly and swirl the vial several times until the tablet is dissolved.
17. Place the vial in the sample chamber making sure that the \times marks are aligned.
18. Press **TEST** key.
19. Remove the vial from the sample chamber.
20. Add **one DPD No. 3 tablet** straight from the foil to the same water sample and crush the tablet using a clean stirring rod.
21. Close the vial with the cap tightly and swirl the vial several times until the tablet is dissolved.
22. Place the vial in the sample chamber making sure that the \times marks are aligned.
23. Press **TEST** key.
Wait for a **reaction period of 2 minutes**.

T1 accepted
prepare T2
press TEST

T2 accepted
prepare T3
press TEST

Countdown
2:00

***,** mg/l ClO₂ [Cl]**

***,** mg/l ClO₂**

***,** mg/l free Cl**
***,** mg/l comb. Cl**
***,** mg/l total Cl**

After the reaction period is finished the reading starts automatically.

The result is shown in the display in:

as Chlorine dioxide in mg/l Chlorine,

or

as Chlorine dioxide in mg/l ClO₂.

mg/l free Chlorine

mg/l combined Chlorine

mg/l total Chlorine

Notes:

See next page.

1.1 Methods

Notes: (Chlorine dioxide in the presence of Chlorine)

1. The conversion factor to convert Chlorine dioxide as Chlorine to Chlorine dioxide as ClO_2 is approximately 0.4 (more exactly 0.38).

$$\text{mg/l } \text{ClO}_2 = \text{mg/l } \text{ClO}_2 [\text{Cl}] \times 0.38$$

▲ $\text{ClO}_2[\text{Cl}]$

▼ ClO_2

(Chlorine dioxide displayed as Chlorine units $\text{ClO}_2[\text{Cl}]$ has its origin out of the swimming poolwater treatment according to DIN 19643.)


2. The total Chlorine result given includes the contribution by the Chlorine dioxide (as Chlorine) reading. For true total Chlorine value subtract the Chlorine dioxide (as Chlorine) reading from the quoted total Chlorine reading.
3. Also see page 51.

1.1 Methods

1 2 0

Chlorine dioxide in absence of Chlorine with Tablet

0.05 – 11 mg/l ClO₂



prepare Zero
press ZERO

1. Fill a clean vial (24 mm Ø) with **10 ml of water sample**, close the vial with the cap tightly.
2. Place the vial in the sample chamber making sure that the X marks are aligned.
3. Press **ZERO** key.
4. Remove the vial from the sample chamber, **empty the vial leaving a few drops in**.
5. Add **one DPD No. 1 tablet** straight from the foil and crush the tablet using a clean stirring rod.
6. Add water sample to the 10 ml mark.
7. Close the vial with the cap tightly and swirl the vial several times until the tablet is dissolved.
8. Place the vial in the sample chamber making sure that the X marks are aligned.

Zero accepted
prepare Test
press TEST

9. Press **TEST** key.

*,** mg/l ClO₂ [Cl]

*,** mg/l ClO₂

The result is shown in the display
as Chlorine dioxide in mg/l Chlorine,
or
as Chlorine dioxide in mg/l ClO₂.

Notes:


See page 51.

1.1 Methods


1 **0** **5**

Chlorine HR (KI) with Tablet

5 – 200 mg/l Cl₂



prepare Zero
press ZERO

1. Fill a clean vial (16 mm Ø) with **10 ml of water sample**, close the vial with the cap tightly.
2. Place the vial in the sample chamber making sure that the marks are  aligned. Place the cover on the adapter.


3. Press **ZERO** key.

4. Remove the vial from the sample chamber.

5. Add **one CHLORINE HR (KI) tablet** straight from the foil to the water sample and crush the tablet using a clean stirring rod.

6. Add **one ACIDIFYING GP tablet** straight from the foil to the same water sample and crush the tablet using a clean stirring rod.

7. Close the vial with the cap tightly and swirl the vial several times until the tablets are dissolved.

8. Place the vial in the sample chamber making sure that the marks are  aligned. Place the cover on the adapter.

Zero accepted
prepare Test
press TEST

9. Press **TEST** key.

The result is shown in the display in mg/l Chlorine.

1.1 Methods

Notes:


1. Oxidizing agents interfere as they react like Chlorine.



1.1 Methods

1 **3** **0**

COD LR with Vario Tube Test

0 – 150 mg/l O₂



1. Open the white cap of one reaction vial and add **2 ml deionised water** (this is the blank (Note 1)).
2. Open the white cap of another reaction vial and add **2 ml water sample** (this is the sample).
3. Close the vials with the cap tightly. Invert the vial gently several times to mix the contents.
(CAUTION: The vial will become hot during mixing!)
4. Heat the vials for **120 minutes** in the preheated reactor at a temperature of **150°C**.
5. **(CAUTION: The vials are hot!)**
Remove the tubes from the heating block and allow them to cool to 60°C or less. Mix the contents by carefully inverting each tube several times while still warm. Then allow the tubes to cool to ambient temperature before measuring. (Note 2).
6. Place the vial (the blank (Note 3, 4)) in the sample chamber making sure that the marks are  aligned. Place the cover on the adapter.
7. Press **ZERO** key.
8. Remove the vial from the sample chamber.
9. Place the vial (the sample (Note 3, 4)) in the sample chamber making sure that the marks are  aligned. Place the cover on the adapter.
10. Press **TEST** key.
The result is shown in the display in mg/l COD.

**prepare Zero
press ZERO**

**Zero accepted
prepare Test
press TEST**

1.1 Methods

Notes:


1. Run samples and blanks with the same batch of vials. The blank is stable when stored in the dark and can be used for further measurements with vials of the same batch.
2. Do not place the hot vials in the sample chamber. Cool the vials to room temperature for final measurements.
3. Suspended solids in the vial lead to incorrect measurements. For this reason it is important to place the vials carefully in the sample chamber. The precipitate at the bottom of the sample should be not suspended.
4. Clean the outside of the vials with a towel. Finger prints or other marks will be removed.
5. Samples can be measured when the Chloride content does not exceed 1000 mg/l.
6. In exceptional cases, compounds contained in the water cannot be oxidized adequate, what results in minimum findings, compared with the reference method.

1.1 Methods





COD MR with Vario Tube Test

0 – 1500 mg/l O₂



Ø 16 mm

1. Open the white cap of one reaction vial and add **2 ml deionised water** (this is the blank (Note 1)).
2. Open the white cap of another reaction vial and add **2 ml water sample** (this is the sample).
3. Close the vials with the cap tightly. Invert the vial gently several times to mix the contents.
(CAUTION: The vial will become hot during mixing!)
4. Heat the vials for **120 minutes** in the preheated reactor at a temperature of **150°C**.
5. **(CAUTION: The vials are hot!)**
Remove the tubes from the heating block and allow them to cool to 60°C or less. Mix the contents by carefully inverting each tube several times while still warm. Then allow the tubes to cool to ambient temperature before measuring. (Note 2).
6. Place the vial (the blank (Note 3, 4)) in the sample chamber making sure that the marks are  aligned. Place the cover on the adapter.
7. Press **ZERO** key.
8. Remove the vial from the sample chamber.
9. Place the vial (the sample (Note 3, 4)) in the sample chamber making sure that the marks are  aligned. Place the cover on the adapter.
10. Press **TEST** key.
The result is shown in the display in mg/l COD.

prepare Zero
press ZERO

Zero accepted
prepare Test
press TEST

1.1 Methods

Notes:


1. Run samples and blanks with the same batch of vials. The blank is stable when stored in the dark and can be used for further measurements with vials of the same batch.
2. Do not place the hot vials in the sample chamber. Cool the vials to room temperature for final measurements.
3. Suspended solids in the vial lead to incorrect measurements. For this reason it is important to place the vials carefully in the sample chamber. The precipitate at the bottom of the sample should be not suspended.
4. Clean the outside of the vials with a towel. Finger prints or other marks will be removed.
5. Samples can be measured when the Chloride content does not exceed 1000 mg/l.
6. In exceptional cases, compounds contained in the water cannot be oxidized adequate, what results in minimum findings, compared with the reference method.
7. For samples under 100 mg/l COD it is recommendable to repeat the test with the tube test for COD LR.

1.1 Methods





COD HR with Vario Tube Test

0 – 15 g/l O₂ (≙ 0 – 15 000 mg/l O₂)



Ø 16 mm

1. Open the white cap of one reaction vial and add **0.2 ml deionised water** (this is the blank (Note 1)).
2. Open the white cap of another reaction vial and add **0.2 ml water sample** (this is the sample).
3. Close the vials with the cap tightly. Invert the vial gently several times to mix the contents.
(CAUTION: The vial will become hot during mixing!)
4. Heat the vials for **120 minutes** in the preheated reactor at a temperature of **150°C**.
5. **(CAUTION: The vials are hot!)**
Remove the tubes from the heating block and allow them to cool to 60°C or less. Mix the contents by carefully inverting each tube several times while still warm. Then allow the tubes to cool to ambient temperature before measuring. (Note 2).
6. Place the vial (the blank (Note 3, 4)) in the sample chamber making sure that the marks are  aligned. Place the cover on the adapter.
7. Press **ZERO** key.
8. Remove the vial from the sample chamber.
9. Place the vial (the sample (Note 3, 4)) in the sample chamber making sure that the marks are  aligned. Place the cover on the adapter.
10. Press **TEST** key.
The result is shown in the display in **g/l** COD.

prepare Zero
press ZERO

Zero accepted
prepare Test
press TEST

1.1 Methods

Notes:

1. Run samples and blanks with the same batch of vials. The blank is stable when stored in the dark and can be used for further measurements with vials of the same batch.
2. Do not place the hot vials in the sample chamber. Cool the vials to room temperature for final measurements.
3. Suspended solids in the vial lead to incorrect measurements. For this reason it is important to place the vials carefully in the sample chamber. The precipitate at the bottom of the sample should be not suspended.
4. Clean the outside of the vials with a towel. Finger prints or other marks will be removed.
5. Samples can be measured when the Chloride content does not exceed 1000 mg/l.
6. In exceptional cases, compounds contained in the water cannot be oxidized adequate, what results in minimum findings, compared with the reference method.
7. For samples under 1 g/l COD it is recommendable to repeat the test with the test kit for COD MR or for samples under 0,1 g/l COD with the tube test COD LR.

1.1 Methods

1 5 0

Copper with Tablet

0.05 – 5 mg/l Cu

Copper
>> diff
free
total

The following selection is shown in the display:

>> diff

for the differentiated determination of free, combined and total Copper.

>> free

for the determination of free Copper.

>> total

for the determination of total Copper.

Select the desired determination with the arrow keys [▲] and [▼]. Confirm with [↵] key.

Note:

1. If ??? is displayed at the differentiated test result see page 242.

1.1 Methods

1 5 0



Ø 24 mm

prepare Zero
press ZERO

Zero accepted
prepare T1
press TEST

T1 accepted
prepare T2
press TEST

*,** mg/l free Cu
*,** mg/l comb Cu
*,** mg/l total Cu

Copper, differentiated determination with Tablet

0.05 – 5 mg/l Cu

1. Fill a clean vial (24 mm Ø) with **10 ml of water sample**, close the vial with the cap tightly.
2. Place the vial in the sample chamber making sure that the Σ marks are aligned.
3. Press **ZERO** key.
4. Remove the vial from the sample chamber.
5. Add **one COPPER No. 1 tablet** straight from the foil to the water sample and crush the tablet using a clean stirring rod.
6. Close the vial with the cap tightly and swirl the vial several times until the tablet is dissolved.
7. Place the vial in the sample chamber making sure that the Σ marks are aligned.
8. Press **TEST** key.
9. Remove the vial from the sample chamber.
10. Add **one COPPER No. 2 tablet** straight from the foil to the same water sample and crush the tablet using a clean stirring rod.
11. Close the vial with the cap tightly and swirl the vial several times until the tablet is dissolved.
12. Place the vial in the sample chamber making sure that the Σ marks are aligned.
13. Press **TEST** key.

The result is shown in the display in:

mg/l free Copper
mg/l combined Copper
mg/l total Copper

1.1 Methods



Copper, free with Tablet

0.05 – 5 mg/l Cu



1. Fill a clean vial (24 mm Ø) with **10 ml of water sample**, close the vial with the cap tightly.

2. Place the vial in the sample chamber making sure that the \times marks are aligned.

**prepare Zero
press ZERO**

3. Press **ZERO** key.

4. Remove the vial from the sample chamber.

5. Add **one COPPER No. 1 tablet** straight from the foil to the water sample and crush the tablet using a clean stirring rod.

6. Close the vial with the cap tightly and swirl the vial several times until the tablet is dissolved.

7. Place the vial in the sample chamber making sure that the \times marks are aligned.

**Zero accepted
prepare Test
press TEST**

8. Press **TEST** key.


The result is shown in the display in mg/l free Copper.

1.1 Methods

1 5 0

Copper, total with Tablet

0.05 – 5 mg/l Cu



1. Fill a clean vial (24 mm Ø) with **10 ml of water sample**, close the vial with the cap tightly.
2. Place the vial in the sample chamber making sure that the X marks are aligned.

**prepare Zero
press ZERO**

3. Press **ZERO** key.
4. Remove the vial from the sample chamber.
5. Add **one COPPER No. 1 tablet and one COPPER No. 2 tablet** straight from the foil to the water sample and crush the tablets using a clean stirring rod.
6. Close the vial with the cap tightly and swirl the vial several times until the tablets are dissolved.
7. Place the vial in the sample chamber making sure that the X marks are aligned.

**Zero accepted
prepare Test
press TEST**

8. Press **TEST** key.

The result is shown in the display in mg/l total Copper.

1.1 Methods


1 5 3

Copper, free (Note 1) with Vario Powder Pack

0.05 – 5 mg/l Cu



prepare Zero
press ZERO



Zero accepted
prepare Test
press TEST

Count-Down
2:00

1. Fill a clean vial (24 mm Ø) with **10 ml of water sample**, close the vial with the cap tightly.
 2. Place the vial in the sample chamber making sure that the Σ marks are aligned.
 3. Press **ZERO** key.
 4. Remove the vial from the sample chamber.
 5. Add **one VARIO Cu 1 F10 powder pack** straight from the foil to the water sample.
 6. Close the vial with the cap tightly and swirl the vial several times to mix the contents (Note 3).
 7. Place the vial in the sample chamber making sure that the Σ marks are aligned.
 8. Press **TEST** key.
- Wait for a **reaction period of 2 minutes**.

After the reaction period is finished the reading starts automatically.

The result is shown in the display in mg/l Copper

1.1 Methods

Notes:

1. For determination of total Copper digestion is required.
2. Extremely acid water samples (pH 2 or less) must be adjusted between pH 4 and pH 6 before the reagent is added (with 8 mol/l Potassium hydroxide solution KOH).
3. Accuracy is not affected by undissolved powder.
4. Interferences:

Cyanide, CN ⁻	Cyanide prevents full colour development. Add 0.2 ml Formaldehyde to 10 ml water sample and wait for a reaction time of 4 minutes (Cyanide is masked). After this perform test as described. Multiply the result by 1.02 to correct the sample dilution by Formaldehyde.
Silver, Ag ⁺	If a turbidity remains and turns black, silver interferences is likely. Add 10 drops of saturated Potassium chloride solution to 75 ml of water sample. Filtrate through a fine filter. Use 10 ml of the filtered water sample to perform test.

1.1 Methods

1

5

7

Cyanide with Reagent Test

0.01 – 0.5 mg/l CN



prepare Zero
press ZERO

1. Fill a clean vial (24 mm Ø) with **2 ml of water sample and 8 ml of deionized water**, close the vial with the cap tightly.
2. Place the vial in the sample chamber making sure that the Σ marks are aligned.
3. Press **ZERO** key.
4. Remove the vial from the sample chamber.
5. Add **two level grey (No. 4) spoon Cyanide-11** to the prepared water sample, replace the cap tightly and invert the vial several times to mix the contents.
6. **Add two level grey (No. 4) spoon Cyanide-12**, replace the cap tightly and invert the vial several times to mix the contents.
7. Fill the vial with drops of the same size by holding the bottle vertically and squeeze slowly:

3 drops of Cyanide-13

8. Close the vial with the cap tightly and invert the vial several times to mix the contents.
9. Place the vial in the sample chamber making sure that the Σ marks are aligned.
10. Press **TEST** key.

Wait for a reaction **period of 10 minutes**.

After the reaction period is finished the reading starts automatically.

The result is shown in the display in mg/l Cyanide.

Zero accepted
prepare Test
press TEST

Count-Down
10:00

1.1 Methods

Notes:

1. Only free Cyanide and Cyanides that can be destroyed by Chlorine are determined by this test.
2. In the present of Thiocyanate, heavy metal complexes, colorants or aromatic amines, the cyanide must be separated out by distillation before analysis is performed.
3. **Store the reagents in closed containers at a temperature of + 15°C to + 25°C.**

1.1 Methods



Cyanuric acid with Tablet

2 – 160 mg/l Cys



1. Fill a clean vial (24 mm Ø) with **5 ml water sample** and **5 ml deionised water (Note 1)**, close the vial with the cap tightly.

2. Place the vial in the sample chamber making sure that the Σ marks are aligned.

prepare Zero
press ZERO

3. Press **ZERO** key.

4. Remove the vial from the sample chamber.

5. Add **one CYANURIC ACID tablet** straight from the foil to the prepared water sample and crush the tablet using a clean stirring rod.

6. Close the vial with the cap tightly and swirl the vial several times until the tablet is dissolved (Note 2, 3).

7. Place the vial in the sample chamber making sure that the Σ marks are aligned.

Zero accepted
prepare Test
press TEST

8. Press **TEST** key.

The result is shown in the display in mg/l Cyanuric acid.

1.1 Methods

Notes:


1. Use deionised water or tap water free of Cyanuric acid.
2. If Cyanuric acid is present a cloudy solution will be given.
Single particles are uncaused necessarily by Cyanuric acid.
3. Dissolve the tablet completely (therefore swirl the vial approx. 1 minute).
Not dissolved particles of the tablet can cause too high results.

1.1 Methods

1 6 5

DEHA (N,N-Diethylhydroxylamine) with Tablet and Liquid Reagent

20 – 500 µg/l DEHA / 0.02 – 0.5 mg/l DEHA



prepare Zero
press ZERO

1. Fill a clean vial (24 mm Ø) with **10 ml of water sample**, close the vial with the cap tightly (Note 2).

2. Place the vial in the sample chamber making sure that the X marks are aligned.

3. Press **ZERO** key.

4. Remove the vial from the sample chamber.

5. Fill the vial with drops of the same size by holding the bottle vertically and squeeze slowly:

6 drops (0.25ml) of DEHA solution

6. Close the vial with the cap tightly and swirl the vial several times to mix the contents.

7. Add **one DEHA tablet** straight from the foil to the same water sample and crush the tablet using a clean stirring rod.

8. Close the vial with the cap tightly and swirl the vial several times until the tablet is dissolved.

9. Place the vial in the sample chamber making sure that the X marks are aligned.

Zero accepted
prepare Test
press TEST

10. Press TEST key.

Wait for a **reaction period of 10 minutes**.

After the reaction period is finished the reading starts automatically.

The result is shown in the display as DEHA.

Count-Down
10:00

1.1 Methods

Notes:

1. Application: Testing of residual corrosion inhibitors (Oxygen scavengers) in boiler feed water or condensate.
2. Before using clean the vials with Hydrochloric acid (approx. 20%). Rinse thoroughly with deionized water.
3. Keep the sample dark during colour development time. UV-light (sunlight) causes too high measurement results.
4. Ideal temperature for full colour development is $20^{\circ}\text{C} \pm 2^{\circ}\text{C}$.
5. Interferences:
 - Iron (II) interferes at all concentrations:
Repeat the test procedure but without adding the DEHA solution. If the displayed result is above $20\ \mu\text{g/l}$ subtract this value from the DEHA test result.
 - Substances which reduce Iron (III) interfere. Substances which complex iron strongly may interfere also.
 - Substances which may interfere when present in concentrations at:

Borate (as $\text{Na}_2\text{B}_4\text{O}_7$)	500 mg/l
Cobalt	0.025 mg/l
Copper	8.0 mg/l
Hardness (as CaCO_3)	1000 mg/l
Lignosulfonates	0.05 mg/l
Manganese	0.8 mg/l
Molybdenum	80 mg/l
Nickel	0.8 mg/l
Phosphate	10 mg/l
Phosphonates	10 mg/l
Sulfate	1000 mg/l
Zinc	50 mg/l

6. There is an option to change the unit from mg/l to $\mu\text{g/l}$.
The unit mg/l is rounded, e.g.: $25\ \mu\text{g/l} = 0.025\ \text{mg/l} \rightarrow$ display 0.03 mg/l.

▲ mg/l

▼ $\mu\text{g/l}$

1.1 Methods


1 6 7

DEHA (N,N-Diethylhydroxylamin) with Vario Powder Pack and Liquid Reagent

20 – 500 µg/l DEHA / 0.02 – 0.5 mg/l DEHA



Use two clean vials (24 mm Ø) and mark one as blank for zeroing (Note 2).



1. Fill a clean vial with **10 ml deionized water** (this is the blank).

2. Fill the second clean vial with **10 ml water sample** (this is the sample).

3. Add **one VARIO OXYSCAV 1 Rgt powder pack** straight from the foil into each vial.

4. Close the vials with the caps tightly and swirl the vials several times to mix the contents.

5. Add **0.20 ml VARIO DEHA 2 Rgt Solution** to each vial (Note 4).

6. Close the vials with the caps tightly and swirl the vials several times to mix the contents.

Count-Down 1
10:00
start: ↵

7. Press **[↵]** key.

Wait for a reaction **period of 10 minutes** (Note 5).

After reaction period is finished proceed as follows:

8. Place the vial (the blank) in the sample chamber making sure that the **X** marks are aligned.

prepare Zero
press ZERO

9. Press **ZERO** key.

10. Remove the vial from the sample chamber.

11. Place the vial (the sample) in the sample chamber making sure that the **X** marks are aligned.

Zero accepted
prepare Test
press TEST

12. Press **TEST** key.

The result is shown in the display as DEHA.

1.1 Methods

Notes:

1. Application: Testing of residual corrosion inhibitors (Oxygen scavengers) in boiler feed water or condensate.
2. Before using clean the vials with Hydrochloric acid (approx. 20%). Rinse thoroughly with deionized water.
3. Ideally temperature for full colour development is $25^{\circ}\text{C} \pm 3^{\circ}\text{C}$.
4. Volume should always be metered by using suitable pipette (class A).
5. Keep blank and sample dark during colour development time. UV-light (sunlight) causes too high measurement results.
6. Interferences:
 - Iron (II) interferes at all concentrations:
Repeat the test procedure but without adding the VARIO DEHA Rgt 2 solution. If the displayed result is above $20\ \mu\text{g/l}$ subtract this value from the DEHA test result.
 - Substances which reduce Iron (III) interfere. Substances which complex iron strongly may interfere also.
 - Substances who may interfere when present in concentrations at:

Borate (as $\text{Na}_2\text{B}_4\text{O}_7$)	500 mg/l
Cobalt	0.025 mg/l
Copper	8.0 mg/l
Hardness (as CaCO_3)	1000 mg/l
Lignosulfonates	0.05 mg/l
Manganese	0.8 mg/l
Molybdenum	80 mg/l
Nickel	0.8 mg/l
Phosphate	10 mg/l
Phosphonates	10 mg/l
Sulfate	1000 mg/l
Zinc	50 mg/l

7. There is an option to change the unit from mg/l to $\mu\text{g/l}$.
The unit mg/l is rounded, e.g.: $25\ \mu\text{g/l} = 0.025\ \text{mg/l} \rightarrow$ display 0.03 mg/l.

▲ mg/l

▼ $\mu\text{g/l}$

1.1 Methods

1 7 0

Fluoride with Liquid Reagent

0.05 – 2 mg/l F



Regard notes!

1. Fill a clean vial (24 mm Ø) with **exact 10 ml of water sample** (Note 4), close the vial with the cap tightly.
2. Place the vial in the sample chamber making sure that the Σ marks are aligned.
3. Press **ZERO** key.
4. Remove the vial from the sample chamber.
5. Add **exact 2 ml SPADNS reagent solution** (Note 4) to the water sample.
Caution: Vial is filled up to the top! (Note 8)
6. Close the vial with the cap tightly and swirl the vial several times to mix the contents.
7. Place the vial in the sample chamber making sure that the Σ marks are aligned.

prepare Zero
press ZERO

Zero accepted
prepare Test
press TEST

Press **TEST** key.

The result is shown in the display in mg/l Fluoride.

1.1 Methods

Notes:

1. The same batch of SPADNS reagent solution must be used for adjustment and test.
The adjustment process needs to be performed for each new batch of SPANDS reagent solution (see Standard Methods 20th, 1998, APHA, AWWA, WEF 4500 F D., S. 4-82).
The procedure is described in chapter 2.4.5 "Calibration" Mode 40" on page 215.
2. During adjustment and test the same vial should be used for zeroing and test, as different vials may exhibit minor tolerances.
3. The calibration solution and the water samples to be tested should have the same temperature ($\pm 1^\circ\text{C}$).
4. As the test result is highly dependent on exact sample and reagent volumes, the sample and reagent volumes should always be metered by using a 10 ml resp. 2 ml volumetric pipette (class A).
5. The accuracy of the test methods decreases above a level of 1.2 mg/l Fluoride. Although the results are sufficiently accurate for most applications, even more exact results can be achieved by 1:1 dilution of the sample prior to use and subsequent multiplication of the result by 2.
6. SPADNS reagent solution contains Arsenite.
Chlorine concentrations up to 5 mg/l do not interfere.
7. Seawater and wastewater samples must be distilled.
8. It is convenient to use special vials with larger volume.

1.1 Methods

1 9 0

Hardness, Calcium with tablet reagent

50 – 900 mg/l CaCO₃



1. Fill a clean vial (24 mm Ø) with **10 ml deionized water**.
2. Add **one CALCHECK tablet** straight from the foil to the deionised water and crush the tablet using a clean stirring rod.
3. Close the vial with the cap tightly and swirl the vial several times until the tablet is dissolved.
4. Place the vial in the sample chamber making sure that the X marks are aligned.

prepare Zero
press ZERO

Countdown
2:00

5. Press **ZERO** key.
Wait for a **reaction period of 2 minutes**.
After the reaction period is finished the reading starts automatically.
6. Remove the vial from the sample chamber.
7. Add **2 ml water sample** to the prepared vial.
Caution: Vial is filled up to the top! (Note 4)
8. Close the vial with the cap tightly and swirl the vial several times (5x) to mix the contents.
9. Place the vial in the sample chamber making sure that the X marks are aligned.

Zero accepted
prepare Test
press TEST

10. Press **TEST** key.
The result is shown in the display as Calcium Hardness.

1.1 Methods

Notes:

1. Strong alkaline or acidic water samples must be adjusted between pH 4 and pH 10 before the tablet is added (use 1 mol/l Hydrochloric acid resp. 1 mol/l Sodium hydroxide).
2. The tolerance of the method is increasing with higher concentrations. When diluting samples, this should be take in account, always measuring in the first third of the range.
3. This method was developed from a volumetric procedure for the determination of calcium. Due to undefined conditions, the deviations from the standardised method may be greater.
4. It is convenient to use special vials with larger volume.
5. ▲ CaCO₃
°dH
°eH
°fH
▼ °aH

1.1 Methods

2 0 0

Hardness, total with tablet reagent

2 – 50 mg/l CaCO₃



prepare Zero
press ZERO

1. Fill a clean vial (24 mm Ø) with **10 ml of water sample**, close the vial with the cap tightly.

2. Place the vial in the sample chamber making sure that the X marks are aligned.

3. Press **ZERO** key.

4. Remove the vial from the sample chamber.

5. Add **one HARDCHECK P tablet** straight from the foil to the water sample and crush the tablet using a clean stirring rod.

6. Close the vial with the cap tightly and swirl the vial several times until the tablet is dissolved.

7. Place the vial in the sample chamber making sure that the X marks are aligned.

Zero accepted
prepare Test
press TEST

8. Press **TEST** key.

Wait for a **reaction period of 5 minutes**.

Countdown
5:00

After the reaction period is finished the reading starts automatically.

The result is shown in the display as total Hardness.

1.1 Methods

Notes:

1. Strong alkaline or acidic water samples must be adjusted between pH 4 and pH 10 before the tablet is added (use 1 mol/l Hydrochloric acid resp. 1 mol/l Sodium hydroxide).
2. Conversion table:

	mg/l CaCO ₃	°dH	°fH	°eH
1 mg/l CaCO ₃	----	0,056	0,10	0,07
1 °dH	17,8	----	1,78	1,25
1 °fH	10,0	0,56	----	0,70
1 °eH	14,3	0,80	1,43	----

3. ▲ CaCO₃
°dH
°eH
°fH
▼ °aH

1.1 Methods

2 0 1

Hardness, total HR with Tablet

20 – 500 mg/l CaCO₃



Ø 24 mm

1. Fill a clean vial (24 mm Ø) with **1 ml of water sample** and **9 ml of deionised water**, close the vial with the cap tightly.
2. Place the vial in the sample chamber making sure that the Σ marks are aligned.

prepare Zero
press ZERO

3. Press **ZERO** key.
4. Remove the vial from the sample chamber.
5. Add **one HARDCHECK P tablet** straight from the foil to the water sample and crush the tablet using a clean stirring rod.
6. Close the vial with the cap tightly and swirl the vial several times until the tablet is dissolved.
7. Place the vial in the sample chamber making sure that the Σ marks are aligned.

Zero accepted
prepare Test
press TEST

Countdown
5:00

8. Press **TEST** key.
Wait for a **reaction period of 5 minutes**.

After the reaction period is finished the reading starts automatically.

The result is shown in the display as total Hardness.

1.1 Methods

Notes:

1. Strong alkaline or acidic water samples must be adjusted between pH 4 and pH 10 before the tablet is added (use 1 mol/l Hydrochloric acid resp. 1 mol/l Sodium hydroxide).
2. Conversion table:

	mg/l CaCO ₃	°dH	°fH	°eH
1 mg/l CaCO ₃	----	0,056	0,10	0,07
1 °dH	17,8	----	1,78	1,25
1 °fH	10,0	0,56	----	0,70
1 °eH	14,3	0,80	1,43	----

3. ▲ CaCO₃
°dH
°eH
°fH
▼ °aH

1.1 Methods

2 0 5

Hydrazine with powder reagent

0.05 – 0.5 mg/l N_2H_4 / 50 – 500 μ g/l N_2H_4



1. Fill a clean vial (24 mm Ø) with **10 ml of water sample** (Note 1, 2), close the vial with the cap tightly.
2. Place the vial in the sample chamber making sure that the Σ marks are aligned.

prepare Zero
press ZERO

3. Press **ZERO** key.
4. Remove the vial from the sample chamber.
5. Add **1 g HYDRAZINE test powder** (Note 3) to the water sample.
6. Close the vial with the cap tightly and swirl the vial several times to mix the contents.

Count-Down
10:00
start: ↵

7. Press [\downarrow] key.
Wait for a **reaction period of 10 minutes**.
After reaction period is finished proceed as follows:
8. The slight turbidity occurring when the reagent is added must be removed by filtration (Note 4).
9. Place the vial in the sample chamber making sure that the Σ marks are aligned.

Zero accepted
prepare Test
press TEST

10. Press **TEST** key.
The result is shown in the display as Hydrazine.

1.1 Methods

Notes:

1. If the water sample is cloudy, you must filter it before performing the zero calibration.
2. The temperature of the water sample should not exceed 21°C.
3. Using the Hydrazine spoon: 1 g is equivalent to one level spoon.
4. Qualitative folded filter papers for medium precipitates are recommend.
5. In order to check whether the reagent has aged (if it has been stored for a lengthy period), perform the test as described above using tap water. If the result is above the detection limit of 0.05 mg/l, you should only use the reagent with reservations (major result deviation).
6. There is an option to change the unit from mg/l to µg/l.
The unit mg/l is rounded, e.g.: 25 µg/l = 0.025 mg/l → display 0.03 mg/l.

▲ mg/l

▼ µg/l

1.1 Methods

2 0 6

Hydrazine with Vario Liquid Reagent

0.01 – 0.6 mg/l N₂H₄ / 5 – 600 µg/l N₂H₄



Ø 24 mm

Use two clean vials (24 mm Ø) and mark one as blank for zeroing.

1. Fill a clean vial with **10 ml deionized water** (this is the blank).
2. Add **1 ml VARIO Hydra 2 Rgt Solution** into the vial (Note 3).
3. Close the vial with the cap tightly and swirl the vial several times to mix the contents.
4. Place the vial (the blank) in the sample chamber making sure that the X marks are aligned.
5. Press **ZERO** key.
6. Remove the vial from the sample chamber.
7. Fill the second clean vial with **10 ml water sample** (this is the sample).
8. Add **1 ml VARIO Hydra 2 Rgt Solution** into the vial.
9. Close the vial with the cap tightly and swirl the vial several times to mix the contents.
10. Place the vial (the blank) in the sample chamber making sure that the X marks are aligned.

prepare Zero
press ZERO

Zero accepted
prepare Test
press TEST

Countdown
12:00

11. Press **TEST** key.
Wait for a **reaction period of 12 minutes**.

After the reaction period is finished the reading starts automatically.

The result is shown in the display as Hydrazine.

1.1 Methods

Notes:

1. Samples cannot be preserved and must be analysed immediately.
2. Sample temperature should be $21^{\circ}\text{C} \pm 4^{\circ}\text{C}$.
3. Caused by the reagent itself the blank may develop a faint yellow colour.
4. Interferences:
 - Ammonium causes no interferences up to 10 mg/l.
At a concentration of 20 mg/l it is possible that the test result increase up to 20%.
 - Morpholine does not interfere up to 10 mg/l.
 - Highly coloured or turbid samples:
Mix 1 part deionised water with 1 part household bleach. Add 1 drop of these mixture into 25 ml water sample and mix. Use 10 ml prepared sample in place of deionised water in point 1.
Note: at point 7 use the unprepared water sample.
Principle: Hydrazine is oxidised by the household bleach. The interference by colour will be eliminated by zeroing.
5. There is an option to change the unit from mg/L to $\mu\text{g/L}$.
The unit mg/L is rounded, e.g.: $25 \mu\text{g/L} = 0.025 \text{ mg/L} \rightarrow \text{display } 0.03 \text{ mg/L}$.

▲ mg/l
▼ $\mu\text{g/l}$

1.1 Methods

2 0 7



Hydrazine with Vacu-vials® K-5003 (see Notes)

0.01 – 0.7 mg/l N₂H₄ / 10 – 700 µg/l N₂H₄

Insert the adaptor for 13 mm Ø vials.

1. Place the blank in the sample chamber. The blank is part of the test kit.

prepare Zero
press ZERO



2. Press **ZERO** key.
3. Remove the blank from the sample chamber.
4. Fill the sampler to the 25 ml mark with the water sample.
5. Place one Vacu-vial® in the sampler. Snap the tip by pressing the vial against the side of the sampler. The Vacu-vial® breaks at the neck and the vial fills automatically. A small volume of inert gas remains in the Vacu-vial®.
6. Mix the content of the Vacu-vial® by inverting it several times, allowing the bubble to move from one end to the other. Dry the outside of the vial.
7. Place the Vacu-vial® in the sample chamber.

Zero accepted
prepare Test
press TEST

Count-Down
10:00

8. Press **TEST** key.
Wait for a **reaction period of 10 minutes**.
After the reaction period is finished the reading starts automatically.
The result is shown in the display as Hydrazine.

1.1 Methods

Notes:

1. This method is adapted from CHEMetrics.
2. Read the original test instruction and the MSDS (delivered with the test) before performing the test. MSDS also at www.chemetrics.com available.
3. Vacu-vials® is a registered trade mark of the company CHEMetrics, Inc. / Calverton, U.S.A.
4. There is an option to change the unit from mg/l to µg/l.
The unit mg/l is rounded, e.g.: 25 µg/l = 0.025 mg/l → display 0.03 mg/l.

▲ mg/l

▼ µg/l

1.1 Methods



Hydrogen peroxide with tablet reagent

0.03 – 3 mg/l H₂O₂



prepare Zero
press ZERO

1. Fill a clean vial (24 mm Ø) with **10 ml of water sample**, close the vial with the cap tightly.
2. Place the vial in the sample chamber making sure that the Σ marks are aligned.
3. Press **ZERO** key.
4. Remove the vial from the sample chamber, **empty the vial leaving a few drops in**.
5. Add **one HYDROGENPEROXIDE LR tablet** straight from the foil and crush the tablet using a clean stirring rod.
6. Add water sample to the 10 ml mark.
7. Close the vial with the cap tightly and swirl the vial several times until the tablet is dissolved.
8. Place the vial in the sample chamber making sure that the Σ marks are aligned.

Zero accepted
prepare Test
press TEST

Countdown
2:00

9. Press **TEST** key.
Wait for a **reaction period of 2 minutes**.

After the reaction period is finished the reading starts automatically.

The result is shown in the display in mg/l Hydrogen peroxide.

1.1 Methods

Notes:

1. Vial cleaning:

As many household cleaners (e.g. dishwasher detergent) contain reducing substances, the subsequent determination of Hydrogen peroxide may show lower results. To avoid any measurement errors, only use glassware free of Chlorine consumption.

Preparation: Put all applicable glassware into Sodium hypochlorite solution (0.1 g/l) for one hour, then rinse all glassware thoroughly with deionized water.

2. Preparing the sample:

When preparing the sample, the escape of Hydrogen peroxide gases, e.g. by pipetting or shaking, must be avoided. The analysis must take place immediately after taking the sample.

3. The DPD colour development is carried out at a pH value of 6.3 to 6.5. The reagent tablet therefore contains a buffer for the pH adjustment.


Strong alkaline or acidic water samples must be adjusted between pH 6 and pH 7 before the tablet is added (use 0.5 mol/l Sulfuric acid resp. 1 mol/l Sodium hydroxide).

4. Exceeding of the measuring range:

Concentrations above 5 mg/l Hydrogen peroxide can produce results within the measuring range up to 0 mg/l. In this event, the water sample must be diluted with water free of Hydrogen peroxide. 10 ml of the diluted sample will be mixed with the reagent and the measurement repeated.


Oxidizing agents such as Chlorine, Ozone etc. interfere as they react like Hydrogen peroxide.

1.1 Methods



Iodine with Tablet

0.05 – 3.6 mg/l I



**prepare Zero
press ZERO**

1. Fill a clean vial (24 mm Ø) with **10 ml of water sample**, close the vial with the cap tightly.
2. Place the vial in the sample chamber making sure that the Σ marks are aligned.
3. Press **ZERO** key.
4. Remove the vial from the sample chamber, **empty the vial leaving a view drops in.**
5. Add **one DPD No. 1 tablet** straight from the foil to the water sample and crush the tablet using a clean stirring rod.
6. Add water sample to the 10 ml mark.
7. Close the vial with the cap tightly and swirl the vial several times until the tablet is dissolved.
8. Place the vial in the sample chamber making sure that the Σ marks are aligned.

**Zero accepted
prepare Test
press TEST**

9. Press **TEST** key.

The result is shown in the display in mg/l Iodine.

1.1 Methods

Notes:

1. Oxidising reagents, such as Chlorine, Bromine, etc. interfere as they react like Iodine.

1.1 Methods

2 2 0

Iron with Tablet

0.02 – 1 mg/l Fe

*Determination of total dissolved Iron Fe^{2+} and Fe^{3+}

2 2 2

Iron with Vario Powder Pack

0.02 – 3 mg/l Fe

*Determination of all soluble iron and most insoluble forms of iron.

2 2 3

Iron, total with Vario Powder Pack

0.02 – 1.8 mg/l Fe


*Determination of all soluble iron and most insoluble forms of iron; most insoluble iron oxides are recovered by the reagent.

*This information refer to analysis of the water sample without digestion.

Further information you will find at the method notes.

1.1 Methods

Notes:



Digestion procedure for the determination of total soluble and insoluble iron.

1. Add 1 ml of concentrated sulfuric acid to 100 ml water sample. Heat and boil for 10 minutes or until all particles are dissolved. After cooling down the sample is set to a pH-value of 3 to 6 by using ammonia solution. Refill with deionised water to the previous volume of 100 ml and mix well. 10 ml of this pre-treated solution is used for the following analysis. Perform as described at the selected test method.
2. Water which has been treated with organic compounds like corrosion inhibitors must be oxidised where necessary to break down the iron. Therefore add 1 ml concentrated sulfuric acid and 1 ml concentrated nitric acid to 100 ml water sample and boil to approx. half volume. After cooling down proceed as described above.

1.1 Methods

2 2 0

Iron (Note 1) with Tablet

0.02 – 1 mg/l Fe



prepare Zero
press ZERO

1. Fill a clean vial (24 mm Ø) with **10 ml of water sample**, close the vial with the cap tightly.
2. Place the vial in the sample chamber making sure that the Σ marks are aligned.
3. Press **ZERO** key.
4. Remove the vial from the sample chamber.
5. Add **one IRON LR tablet** straight from the foil to the water sample and crush the tablet using a clean stirring rod.
6. Close the vial with the cap tightly and swirl the vial several times until the tablet is dissolved.
7. Place the vial in the sample chamber making sure that the Σ marks are aligned.

Zero accepted
prepare Test
press TEST

8. Press **TEST** key.
Wait for a **reaction period of 5 minutes**.

Countdown
5:00

After the reaction period is finished the reading starts automatically.

The result is shown in the display in mg/l Iron.

1.1 Methods

Notes:

1. This method determines the total dissolved Iron as Fe^{2+} and Fe^{3+} .
2. The IRON (II) LR tablet is used for differentiation – as described above – instead of the IRON LR tablet.
3. For the determination of total dissolved and undissolved iron digestion is required. An example is described on page 97.

1.1 Methods



Iron (Note 1) with Vario Powder Pack

0.02 – 3 mg/l Fe



Ø 24 mm


1. Fill a clean vial (24 mm Ø) with **10 ml of water sample**, close the vial with the cap tightly.

2. Place the vial in the sample chamber making sure that the marks \times are aligned.

prepare Zero
press ZERO

3. Press **ZERO** key.

4. Remove the vial from the sample chamber.



5. Add **one Vario Ferro F10 powder pack straight** from the foil to the water sample.

6. Close the vial with the cap tightly and swirl the vial several times to mix the contents (Note 4).

7. Place the vial in the sample chamber making sure that the \times marks are aligned.

Zero accepted
prepare Test
press TEST

8. Press **TEST** key.

Wait for a **reaction period of 3 minutes (Note 5)**.

Countdown
3:00

After the reaction period is finished the reading starts automatically.

The result is shown in the display in mg/l Iron.

1.1 Methods

Notes:

1. The reagent reacts with all soluble iron and most insoluble forms of iron in the water sample.
2. Iron oxide requires a prior digestion, use mild, vigorous or Digesdahl digestion (e.g. for digestion with acid see page 97).
3. Very strong alkaline or acidic water samples must be adjusted to a ph-Value between 3 and 5 before analysis.
4. Accuracy is not affected by undissolved powder.
5. Water samples containing visible rust should be allowed to react at least five minutes.

1.1 Methods




Iron, total (TPTZ, Note 1) with Vario Powder Pack

0.02 – 1.8 mg/l Fe



Use two clean vials (24 mm Ø) and mark one as blank for zeroing.

1. Fill a clean vial with **10 ml deionized water** (this is the blank).



2. Fill the second clean vial with **10 ml water sample** (this is the sample).

3. Add **one Vario IRON TPTZ F10 powder pack** straight from the foil into each vial.

4. Close the vials with the caps tightly and swirl the vials several times to mix the contents.

Count-Down
3:00
start: ↵

5. Press **[↵]** key.
Wait for a reaction **period of 3 minutes**.

After reaction period is finished proceed as follows:

6. Place the vial (the blank) in the sample chamber making sure that the **Σ** marks are aligned.

prepare Zero
press ZERO

7. Press **ZERO** key.

8. Remove the vial from the sample chamber.

9. Place the vial (the sample) in the sample chamber making sure that the **Σ** marks are aligned.

Zero accepted
prepare Test
press TEST

10. Press **TEST** key.
The result is shown in the display in mg/l Iron.

1.1 Methods

Notes:

1. For determination of total Iron digestion is required.
TPTZ reagent recovers most insoluble iron oxides without digestion.
2. Rinse all glassware with 1:1 Hydrochloric acid solution first and then rinse with deionised water to remove iron deposits that can cause slightly high results.
3. Strong alkaline or acidic water samples must be adjusted between pH 3 and pH 8 before the reagent is added (use 0.5 mol/l Sulfuric acid resp. 1 mol/l Sodium hydroxide).
4. Interferences:
When interferences occurred, the colour development was inhibited or a precipitate was formed.

The values below refer to a standard with an iron concentration of 0.5 mg/l.

The following substances do not interfere when present up to the levels given:

Substance	no interference to
Cadmium	4.0 mg/l
Chromium ⁽³⁺⁾	0.25 mg/l
Chromium ⁽⁶⁺⁾	1.2 mg/l
Cobalt	0.05 mg/l
Copper	0.6 mg/l
Cyanide	2.8 mg/l
Manganese	50 mg/l
Mercury	0.4 mg/l
Molybdenum	4.0 mg/l
Nickel	1.0 mg/l
Nitrite Ion	0.8 mg/l

1.1 Methods

2 4 0

Manganese with Tablet

0.2 – 4 mg/l Mn



1. Fill a clean vial (24 mm Ø) with **10 ml of water sample**, close the vial with the cap tightly.
2. Place the vial in the sample chamber making sure that the \times marks are aligned.

prepare Zero
press ZERO

3. Press **ZERO** key.
4. Remove the vial from the sample chamber.
5. Add **one MANGANESE LR 1 tablet** straight from the foil to the water sample and crush the tablet using a clean stirring rod and dissolve the tablet.
6. Add **one MANGANESE LR 2 tablet** straight from the foil to the same water sample and crush the tablet using a clean stirring rod.
7. Close the vial with the cap tightly and swirl the vial several times until the tablets are dissolved.
8. Place the vial in the sample chamber making sure that the marks \times are aligned.

Zero accepted
prepare Test
press TEST

Count-Down
5:00

9. Press **TEST** key.
Wait for a **reaction period of 5 minutes**.
After the reaction period is finished the reading starts automatically.
The result is shown in the display in mg/l Manganese.

1.1 Methods

Note:

1. ▲ Mn
MnO₄
▼ KMnO₄

1.1 Methods


2 4 2

Manganese LR with Vario Powder Pack

0.01 – 0.7 mg/l Mn



Use two clean vials (24 mm Ø) and mark one as blank for zeroing (Note 1).



1. Fill a clean vial with **10 ml of deionized water** (this is the blank).
2. Fill the second clean vial with **10 ml of water sample** (this is the sample).

3. Add **one Vario Ascorbic Acid powder pack** straight from the foil into each vial. (Note 2)

4. Close the vials with the caps tightly and swirl the vials several times to mix the contents.

5. Fill each vial with drops of the same size by holding the bottle vertically and squeeze slowly (Note 3):
15 drops of Alkaline Cyanide reagent solution

6. Close the vials with the caps tightly and swirl the vials several times to mix the contents.

7. Fill each vial with drops of the same size by holding the bottle vertically and squeeze slowly:
21 drops of PAN Indicator solution

8. Close the vials with the caps tightly and swirl the vials several times to mix the contents.

Countdown 1
2:00
start: ↵

9. Press **[↵]** key.
Wait for a **reaction period of 2 minutes** (Note 4).

After reaction period is finished proceed as follows:

9. Place the vial (the blank) in the sample chamber making sure that the marks are **∞** aligned.

prepare Zero
press ZERO

10. Press **ZERO** key.

11. Remove the vial from the sample chamber.

Zero accepted
prepare Test
press TEST

12. Place the vial (the sample) in the sample chamber making sure that the marks are **∞** aligned.

13. Press **TEST** key.

The result is shown in the display in mg/l Manganese.

1.1 Methods

Notes:

1. Rinse all glassware with 1:1 Nitric acid solution first and then rinse with deionised water.
2. Water samples that contain more than 300 mg/l CaCO_3 hardness: After adding the Vario Ascorbic Acid powder pack add additionally 10 drops of Rochelle Salt Solution.
3. After addition of the reagent solution "Alkaline-Cyanide" a cloudy or turbid solution may form in some water samples. The turbidity should disappear after point 7.
4. Water samples containing more than 5 mg/l iron should be allowed to react at least 10 minutes.
5. Conversion:
 $\text{mg/l MnO}_4 = \text{mg/l Mn} \times 2.17$
6. ▲ Mn
 MnO_4
 ▼ KMnO_4

1.1 Methods

2 4 3


Manganese HR with Vario Powder Pack

0.1 – 18 mg/l Mn



prepare Zero
press ZERO

1. Fill a clean vial (24 mm Ø) with **10 ml of water sample**, close the vial with the cap tightly.
2. Place the vial in the sample chamber making sure that the Σ marks are aligned.
3. Press **ZERO** key.
4. Remove the vial from the sample chamber.



5. Add **one VARIO Citrat powder pack** straight from the foil to the water sample.
6. Close the vial with the cap tightly and swirl the vial several times to mix the contents.
7. **Add one VARIO Sodium periodate powder pack** straight from the foil to the same water sample.
8. Close the vial with the cap tightly and swirl the vial several times to mix the contents.
9. Place the vial in the sample chamber making sure that the Σ marks are aligned
10. Press **TEST** key.

Zero accepted
prepare Test
press TEST

Countdown
2:00

Wait for a **reaction period of 2 minutes**.

After the reaction period is finished the reading starts automatically.

The result is shown in the display in mg/l Manganese.

1.1 Methods

Notes:

1. This test is applicable for the determination of soluble Manganese in water and wastewater.
2. Highly buffered or water sample with extreme pH-values may exceed the buffering capacity of the reagents and requires sample pre-treatment.
If samples were acidified for storing, adjust the pH between 4 and 5 with 5 mol/l (5 N) Sodium hydroxide before test. Do not exceed pH 5, as manganese may precipitate.
3. Interferences:

Interfering substance	Interference level
Calcium	greater than 700 mg/l
Chloride	greater than 70 000 mg/l
Iron	greater than 5 mg/l
Magnesium	greater than 100 000 mg/l

4. ▲ Mn
MnO₄
▼ KMnO₄

1.1 Methods


2 5 0

Molybdate with Tablet

1 – 50 mg/l MoO_4 / 0.6 – 30 mg/l Mo



prepare Zero
press ZERO



Zero accepted
prepare Test
press TEST

1. Fill a clean vial (24 mm Ø) with **10 ml of water sample**, close the vial with the cap tightly.
2. Place the vial in the sample chamber making sure that the X marks are aligned.
3. Press **ZERO** key.
4. Remove the vial from the sample chamber and **empty the vial**.
5. Fill **20 ml water sample** in a 100 ml beaker.
6. Add **one MOLYBDATE HR No. 1 tablet** straight from the foil to the water sample and crush the tablet using a clean stirring rod.
7. Add **one MOLYBDATE HR No. 2 tablet** straight from the foil to the same water sample and crush the tablet using a clean stirring rod.
8. Dissolve the tablets using a clean stirring rod.
9. Rinse out the vial with the prepared water sample and then fill to the 10 ml mark.
10. Close the vial with the cap tightly.
11. Place the vial in the sample chamber making sure that the X marks are aligned.
12. Press **TEST** key.

The result is shown in the display in mg/l Molybdate / Molybdenum.

1.1 Methods

Notes:


1. The tablets must be added in the correct sequence.
2. Under test conditions (pH 3.8 – 3.9) iron does not interfere nor do other metals at levels likely to be found in industrial water systems.
3. Conversions:
 $\text{mg/l Mo} = \text{mg/l MoO}_4 \times 0.6$
 $\text{mg/l Na}_2\text{MoO}_6 = \text{mg/l MoO}_4 \times 1.3$
4. ▲ MoO₄
Mo
▼ Na₂MoO₄

1.1 Methods


2 5 2

Molybdate / Molybdenum HR with Vario Powder Pack

0.5 – 66 mg/l MoO₄ / 0.3 – 40 mg/l Mo



prepare Zero
press ZERO



1. Fill a clean vial (24 mm Ø) with **10 ml of water sample**, close the vial with the cap tightly.
2. Place the vial in the sample chamber making sure that the X marks are aligned.
3. Press **ZERO** key.
4. Remove the vial from the sample chamber.
5. Add **one Vario Molybdenum HR 1 F10 powder pack** straight from the foil to the water sample.
6. Close the vial with the cap tightly and swirl the vial several times to mix the contents.
7. Add **one Vario Molybdenum HR 2 F10 powder pack** straight from the foil to the same water sample.
8. Close the vial with the cap tightly and swirl the vial several times to mix the contents.
9. Add **one Vario Molybdenum HR 3 F10 powder pack** straight from the foil to the same water sample.
10. Close the vial with the cap tightly and swirl the vial several times to mix the contents.
11. Place the vial in the sample chamber making sure that the X marks are aligned.
12. Press **TEST** key.

Zero accepted
prepare Test
press TEST

Countdown
5:00

Wait for a **reaction period of 5 minutes**.

After reaction period is finished the reading starts automatically.

The result is shown in the display in mg/l Molybdate / Molybdenum.

1.1 Methods

Notes:

1. Filter turbid water samples using filter paper and funnel before analysis.
2. Highly buffered water samples or extreme pH values should be adjusted to a pH of nearly 7 with 1 mol/l Nitric acid or 1 mol/l Sodium hydroxide.
3. Concentration from 10 mg/l Cu causes too high test values if the described reaction time of 5 minutes is increased. So it is very important to perform the test procedure continuously.
4. Substances who may interfere when present in concentrations at:

Aluminium	50 mg/l
Chromium	1000 mg/l
Iron	50 mg/l
Nickel	50 mg/l
Nitrite	all levels



5. ▲ MoO_4
Mo
▼ Na_2MoO_4

1.1 Methods

2 6 5

Nitrate with Tube Test

1 – 30 mg/l N



Countdown
5:00
start: ↵

1. Open the white cap of one vial (Reagent A) and add **1 ml deionised water** (this is the blank).
2. Open the white cap of the other vial (Reagent A) and add **1 ml water sample** (this is the sample).
3. Add **one Vario Nitrate Chromotropic powder pack straight** from the foil into each vial.
4. Close the vials with the caps tightly and invert the vials gently several times (10 x) to mix the contents (Note 1).

5. Press **[↵]** key.
Wait for a **reaction period of 5 minutes**.
6. After reaction period is finished proceed as follows:
7. Place the vial (the blank) in the sample chamber making sure that the marks are **Δ** aligned.
Place the cover on the adapter.

prepare Zero
press ZERO

8. Press **ZERO** key.
9. Remove the vial from the sample chamber.
10. Place the vial (the sample) in the sample chamber making sure that the marks are **Δ** aligned.
Place the cover on the adapter.

Zero accepted
prepare Test
press TEST

11. Press **TEST** key.
The result is shown in the display in mg/l Nitrate.

1.1 Methods

Notes:

1. Some solids may not dissolve.

2. Conversion:

$$\text{mg/l NO}_3 = \text{mg/l N} \times 4.43$$

3. ▲ N

▼ NO₃

1.1 Methods

2 7 0

Nitrite with Tablet

0.01 – 0.5 mg/l N



prepare Zero
press ZERO

1. Fill a clean vial (24 mm Ø) with **10 ml of water sample**, close the vial with the cap tightly.
2. Place the vial in the sample chamber making sure that the marks X are aligned.
3. Press **ZERO** key.
4. Remove the vial from the sample chamber.
5. Add **one NITRITE LR tablet** straight from the foil to the water sample and crush the tablet using a clean stirring rod.
6. Close the vial with the cap tightly and swirl the vial several times until the tablet is dissolved.
7. Place the vial in the sample chamber making sure that the marks X are aligned.

Zero accepted
prepare Test
press TEST

8. Press **TEST** key.
Wait for a **reaction period of 10 minutes**.

Count-Down
10:00

After the reaction period is finished the reading starts automatically.

The result is shown in the display in mg/l Nitrite.

1.1 Methods

Notes:

1. The following ions can produce interferences since under the reaction conditions they cause precipitation:
Antimony (III), Iron (III), Lead, Mercury (I), Silver, Chloroplatinate, Metavanadate and Bismuth.
Copper (II)-ions may cause lower test results as they accelerate the decomposition of the Diazonium salt.
It is improbable in practice that these interfering ions will occur in such high concentrations that they cause significant reading errors.
2. Conversion:
 $\text{mg/l NO}_2 = \text{mg/l N} \times 3.29$
3. ▲ N
▼ NO₂

1.1 Methods

2 7 2

Nitrite LR with Vario Powder Pack

0.01 – 0.3 mg/l N




1. Fill a clean vial (24 mm Ø) with **10 ml of water sample**, close the vial with the cap tightly.

2. Place the vial in the sample chamber making sure that the \times marks are aligned.

prepare Zero
press ZERO

3. Press **ZERO** key.

4. Remove the vial from the sample chamber.



5. Add **one VARIO Nitri 3 powder pack** straight from the foil to the water sample.

6. Close the vial with the cap tightly and swirl the vial several times to mix the contents.

7. Place the vial in the sample chamber making sure that the \times marks are aligned.

Zero accepted
prepare Test
press TEST

8. Press **TEST** key.

Wait for a **reaction period of 20 minutes**.

Countdown
20:00

After the reaction period is finished the reading starts automatically.

The result is shown in the display in mg/l Nitrite.

1.1 Methods

Notes:

1. Interferences:

- Strong oxidizing and reducing substances interfere.
- Cupric and ferrous ions cause low results.
- Antimonous, Auric, Bismuth, Chloroplatinate, Ferric, Lead, Mercurous, Metavanadate, Silver ions interfere by causing precipitation.
- In samples with very high concentrations of Nitrate (> 100 mg/L N) a small amount of Nitrite will be found. Such high levels of Nitrate appear to undergo a slight amount of reduction to Nitrite, either spontaneously or during the reaction time of the test.

2. ▲ N

▼ NO₂

1.1 Methods


2 8 0

Nitrogen, total LR with Vario Tube Test

0.5 – 25 mg/l N



Ø 16 mm



Countdown

3:00

start: ↵

Countdown

2:00

start: ↵

1. **Open two TN Hydroxide LR digestion vials** and add **one Vario TN Persulfate Rgt. powder pack** (Note 2, 3).
2. Add **2 ml deionised water** to the prepared vial (this is the blank, Note 4, 5).
3. Add **2 ml water sample** to the other prepared vial (this is the sample).
4. Close the vials with the caps and shake to mix the contents (at least 30 seconds, Note 6).
5. Heat the vials for **30 minutes** in the preheated reactor at a temperature of **100°C** (Note 7).
6. After 30 minutes remove the vials from the reactor. **(CAUTION: The vials are hot!)** Allow the vials to cool to room temperature.
7. Open the cooled down digestion vials and add **one Vario TN Reagent A Powder Pack** to each vial (Note 2).
8. Close the vials with the caps and shake to mix the contents (at least 15 seconds).
9. Press [**↵**] key.
Wait for a **reaction period of 3 minutes**.
After reaction period is finished proceed as follows:
10. Open the digestion vials and add **one Vario TN Reagent B powder pack** to each vial (Note 2).
11. Close the vials with the caps and shake to mix the contents (at least 15 seconds, Note 8).
12. Press [**↵**] key.
Wait for a **reaction period of 2 minutes**.
After reaction period is finished proceed as follows:
13. Open **two TN Acid LR/HR (Reagent C) vials** and add **2 ml of the digested, treated blank** to one vial (this is the blank).
14. Add **2 ml of the digested, treated water sample** to the other TN Acid LR/HR vial (this is the sample).
15. Close the vials with the caps and swirl the vials gently several times to mix the contents (10 x, Note 9). **(CAUTION: Vials warm up).**

1.1 Methods

**prepare Zero
press ZERO**

**Count-Down
5:00**

**Zero accepted
prepare Test
press TEST**

16. Place the vial (the blank) in the sample chamber making sure that the marks Δ are aligned.
17. Press **ZERO** key.
Wait for a **reaction period of 5 minutes**.
After the reaction period is finished the reading starts automatically.
18. Remove the vial from the sample chamber.
19. Place the vial (the sample, Note 10) in the sample chamber making sure that the marks Δ are aligned.
20. Press **TEST** key.
The result is shown in the display in mg/l Nitrogen.

Notes:



1. Appropriate safety precautions and a good lab technique should be used during the whole procedure.
2. Use a funnel to add the reagent.
3. Wipe off any Persulfate reagent that may get on the lid or the tube threads.
4. Volumes for samples and blank should always be metered by using 2 ml volumetric pipettes (class A).
5. One blank is sufficient for each set of samples.
6. The reagent may not dissolve completely.
7. It is very important to remove the vials from the reactor after exactly 30 minutes.
8. The reagent will not completely dissolve.
9. Hold the vial in a vertical position with the cap pointing up. Turn the vial upside-down. Wait for all of the solution to flow down to the cap. Return the vial to the upright position. Wait for all the solution to flow to the bottom of the vial. This process is one inversion; 10 inversions = approx. 30 seconds.
10. After zero calibration with the blank it is possible to measure several samples.
11. Great quantities of nitrogen free, organic compounds which are included in some water samples may reduce the effectiveness of the digestion by reacting with the Persulfate reagent. Samples which are well known to contain great quantities of organic compounds must be diluted and digestion and measurement must be repeated for checking the effectiveness of the digestion.
12. Application: for water, wastewater and seawater
13. Interferences:
Interfering substances that resulted in a concentration change of 10%:
Bromide more than 60 mg/l and Chloride more than 1000 mg/l produce positive interferences.
TN = Total Nitrogen
14. \blacktriangle N
 NH₄
 \blacktriangledown NH₃

1.1 Methods

2 **8** **1**

Nitrogen, total HR with Vario Tube Test

5 – 150 mg/l N



Countdown
3:00
start: ↵

Countdown
2:00
start: ↵

1. **Open two TN Hydroxide HR digestion vials** and add **one Vario TN Persulfate Rgt. powder pack** (Note 2, 3).
2. Add **0.5 ml deionised water** to the prepared vial (this is the blank, Note 4, 5).
3. Add **0.5 ml water sample** to the other prepared vial (this is the sample).
4. Close the vials with the caps and shake to mix the contents (at least 30 seconds, Note 6).
5. Heat the vials for **30 minutes** in the preheated reactor at a temperature of **100°C** (Note 7).
6. After 30 Minutes remove the vials from the reactor. **(CAUTION: The vials are hot!)** Allow the vials to cool to room temperature.
7. Open the cooled down digestion vials and add **one Vario TN Reagent A Powder Pack** to each vial (Note 2).
8. Close the vials with the caps and shake to mix the contents (at least 15 seconds).
9. Press [↵] key.
Wait for a **reaction period of 3 minutes**. After reaction period is finished proceed as follows:
10. Open the digestion vials and add **one Vario TN Reagent B powder pack** to each vial (Note 2).
11. Close the vials with the caps and shake to mix the contents (at least 15 seconds, Note 8).
12. Press [↵] key.
Wait for a **reaction period of 2 minutes**. After reaction period is finished proceed as follows:
13. Open **two TN Acid LR/HR (Reagent C) vials** and add **2 ml of the digested, treated blank** to one vial (this is the blank).
14. Add **2 ml of the digested, treated water sample** to the other TN Acid LR/HR vial (this is the sample).
15. Close the vials with the caps and swirl the vials gently several times to mix the contents (10 x, Note 9). **(CAUTION: Vials warm up).**

1.1 Methods

**prepare Zero
press ZERO**

**Count-Down
5:00**

**Zero accepted
prepare Test
press TEST**

16. Place the vial (the blank) in the sample chamber making sure that the ▲ marks are aligned.

17. Press **ZERO** key.
Wait for a **reaction period of 5 minutes**.

After the reaction period is finished the reading starts automatically.

18. Remove the vial from the sample chamber.

19. Place the vial (the sample, Note 10) in the sample chamber making sure that the ▲ marks are aligned.

20. Press **TEST** key.
The result is shown in the display in mg/l Nitrogen.

Notes:

1. Appropriate safety precautions and a good lab technique should be used during the whole procedure.
2. Use a funnel to add the reagent.
3. Wipe off any Persulfate reagent that may get on the lid or the tube threads.
4. Volumes for samples and blank should always be metered by using suitable pipettes (class A).
5. One blank is sufficient for each set of samples.
6. The reagent may not dissolve completely.
7. It is very important to remove the vials from the reactor after exactly 30 minutes.
8. The reagent will not completely dissolve.
9. Hold the vial in a vertical position with the cap pointing up. Turn the vial upside-down. Wait for all of the solution to flow down to the cap. Return the vial to the upright position. Wait for all the solution to flow to the bottom of the vial. This process is one inversion; 10 inversions = approx. 30 seconds.
10. After zero calibration with the blank it is possible to measure several samples.
11. Great quantities of nitrogen free, organic compounds which are included in some water samples may reduce the effectiveness of the digestion by reacting with the Persulfate reagent. Samples which are well known to contain great quantities of organic compounds must be diluted and digestion and measurement must be repeated for checking the effectiveness of the digestion.
12. Application: for water, wastewater and seawater
13. Interferences:
Interfering substances that resulted in a concentration change of 10%:
Bromide more than 60 mg/l and Chloride more than 1000 mg/l produce positive interferences.
TN = Total Nitrogen
14. ▲ N
 NH₄
 ▼ NH₃

1.1 Methods

2 9 0

Oxygen, active* with Tablet

0.1 – 10 mg/l O₂



prepare Zero
press ZERO

1. Fill a clean vial (24 mm Ø) with **10 ml of water sample**, close the vial with the cap tightly.
2. Place the vial in the sample chamber making sure that the \times marks are aligned.
3. Press **ZERO** key.
4. Remove the vial from the sample chamber.
5. Add **one DPD No. 4 tablet** straight from the foil to the water sample and crush the tablet using a clean stirring rod.
6. Close the vial with the cap tightly and swirl the vial several times until the tablet is dissolved.
7. Place the vial in the sample chamber making sure that the \times marks are aligned.

Zero accepted
prepare Test
press TEST

8. Press **TEST** key.
Wait for a **reaction period of 2 minutes**.

Countdown
2:00

After the reaction period is finished the reading starts automatically.

The result is shown in the display in mg/l active Oxygen.

1.1 Methods

Notes:

* **Active Oxygen is a synonym for a common disinfectant (based on "Oxygen") in Swimming Pool Treatment.**

1. When preparing the sample, the escape of Oxygen gases, e.g. by pipetting or shaking, must be avoided.
2. The analysis must take place immediately after taking the sample.

1.1 Methods

2 9 2

Oxygen, dissolved with Vacu-vials® K-7553 (see Notes)


10 – 800 µg/l O₂

Insert the adaptor for 13 mm Ø round vials.

1. Place the blank in the sample chamber. The blank is part of the test kit.

**Zero vorbereiten
ZERO drücken**

2. Press **ZERO** key.
3. Remove the blank from the sample chamber.
4. Water should flow through the special sampler for several minutes to remove any air bubbles sticking at the surface.



The water must flow from the bottom to the top.

5. When the sampler is bubble-free press one Vacu-vial® into the lower edge of the sampler. The Vacu-vial® breaks at the neck and the vial fills automatically.

A small volume of inert gas remains in the Vacu-vial®.

6. Remove the Vacu-vial® point downwards from the sampler immediately.

As the content of the vial has a higher density than water, it is important to remove the vial from the sampler within 5 seconds to prevent any loss of reagent.

7. The Vacu-vial® is closed with one finger (covered with a glove) to prevent entry of air. Invert the vial several times. Dry the outside of the vial.

8. Place the Vacu-vial® in the sample chamber.

**Zero akzeptiert
Test vorbereiten
TEST drücken**

9. Press **TEST** key.

The result is shown in the display in µg/l Oxygen.

1.1 Methods

Notes:

1. This method is adapted from CHEMetrics.
2. Read the original test instruction and the MSDS (delivered with the test) before performing the test. MSDS also at www.chemetrics.com available.
3. Vacu-vials® is a registered trade mark of the company CHEMetrics, Inc. / Calverton, U.S.A.

1.1 Methods

3 **0** **0**

Ozone with Tablet

0.02 – 1 mg/l O₃

Ozon
>> **with Cl**
without Cl

The following selection is shown in the display:

>> **with Cl**

for the determination of Ozone in the presence of Chlorine.

>> **without Cl**

for the determination of Ozone in the absence of Chlorine.

Select the desired method with the arrow keys
[▲] and [▼]. Confirm with [↵] key.

1.1 Methods

Notes:

1. Vial cleaning:

As many household cleaners (e.g. dishwasher detergent) contain reducing substances, the subsequent determination of Ozone may show lower results. To avoid any measurement errors, only use glassware free of Chlorine consumption.

Preparation: Put all applicable glassware into Sodium hypochlorite solution (0.1 g/l) for one hour, then rinse all glassware thoroughly with deionized water.

2. Preparing the sample:

When preparing the sample, the escape of Ozone gases, e.g. by pipetting or shaking, must be avoided. The analysis must take place immediately after taking the sample.

3. The DPD colour development is carried out at a pH value of 6.3 to 6.5. The reagent tablet therefore contains a buffer for the pH adjustment.

Strong alkaline or acidic water samples must be adjusted between pH 6 and pH 7 before the tablet is added (use 0.5 mol/l Sulfuric acid resp. 1 mol/l Sodium hydroxide).

4. Turbidity (lead to errors):

The use of the DPD No. 1 tablet in samples with high Calcium ion content* and/or high conductivity* can lead to turbidity of the sample and therefore incorrect measurements.

** it is not possible to give exactly values, because the development of turbidity depends on nature and ingredients of the sample.*

5. Exceeding of the measuring range:

Concentrations above 5 mg/l Ozone can produce results within the measuring range up to 0 mg/l. In this event, the water sample must be diluted with water free of Ozone. 10 ml of the diluted sample will be mixed with the reagent and the measurement repeated.

6. If ??? is displayed at the differentiated test result see page 242.


Oxidizing agents such as Bromine, Chlorine etc. interfere as they react like Ozone.

1.1 Methods

3 0 0

Ozone, in the presence of Chlorine with Tablet

0.02 – 1 mg/l O₃



prepare Zero
press ZERO

1. Fill a clean vial (24 mm Ø) with **10 ml of water sample**, close the vial with the cap tightly.
 2. Place the vial in the sample chamber making sure that the X marks are aligned.
 3. Press **ZERO** key.
 4. Remove the vial from the sample chamber, **empty the vial leaving a few drops in**.
 5. Add **one DPD No. 1 tablet** and **one DPD No. 3 tablet** straight from the foil and crush the tablets using a clean stirring rod.
 6. Add water sample to the 10 ml mark.
 7. Close the vial with the cap tightly and swirl the vial several times until the tablets are dissolved.
 8. Place the vial in the sample chamber making sure that the X marks are aligned.
 9. Press **TEST** key.
Wait for a **reaction period of 2 minutes**.
- After the reaction period is finished the reading starts automatically.
10. Remove the vial from the sample chamber, empty the vial, rinse vial and cap several times. Fill the vial with **a few drops of water sample**.
 11. Add **one DPD No. 1 tablet** and **one DPD No. 3 tablet** straight from the foil and crush the tablets using a clean stirring rod.

Zero accepted
prepare T1
press TEST

Countdown
2:00

1.1 Methods

12. **Fill a second clean vial with 10 ml of water sample.**
13. Add **one GLYCINE tablet** straight from the foil and crush the tablet using a clean stirring rod.
14. Close the vial with the cap tightly and swirl the vial several times until the tablet is dissolved.
15. **Transfer the content of the second vial into the prepared vial.**
16. Close the vial with the cap tightly and swirl the vial several times until the tablets are dissolved.
17. Place the vial in the sample chamber making sure that the \times marks are aligned.

T1 accepted
prepare T2
press TEST

Count-Down
2:00

..* mg/l O₃
..* mg/l total Cl

18. Press **TEST** key.
Wait for a **reaction period of 2 minutes.**

After the reaction period is finished the reading starts automatically.

The result is shown in the display in:

mg/l Ozone
mg/l total Chlorine

Notes:

See page 129.

1.1 Methods

3 0 0

Ozone, in absence of Chlorine with Tablet

0.02 – 1 mg/l O₃



Ø 24 mm

prepare Zero
press ZERO

1. Fill a clean vial (24 mm Ø) **with 10 ml of water sample**, close the vial with the cap tightly.
2. Place the vial in the sample chamber making sure that the Σ marks are aligned.
3. Press **ZERO** key.
4. Remove the vial from the sample chamber, **empty the vial leaving a few drops in**.
5. Add **one DPD No. 1 tablet** and **one DPD No. 3 tablet** straight from the foil and crush the tablets using a clean stirring rod.
6. Add water sample to the 10 ml mark.
7. Close the vial with the cap tightly and swirl the vial several times until the tablets are dissolved.
8. Place the vial in the sample chamber making sure that the Σ marks are aligned.
9. Press **TEST** key.
Wait for a **reaction period of 2 minutes**.

Zero accepted
prepare Test
press TEST

Countdown
2:00

After the reaction period is finished the reading starts automatically.

The result is shown in the display in
mg/l Ozone.

Notes:

See page 129.

1.1 Methods

7

0

PHMB (Biguanide) with Tablet

2 – 60 mg/l PHMB



prepare Zero
press ZERO

1. Fill a clean vial (24 mm Ø) with **10 ml of water sample**, close the vial with the cap tightly.

2. Place the vial in the sample chamber making sure that the Σ marks are aligned.

3. Press **ZERO** key.

4. Remove the vial from the sample chamber.

5. Add **one PHMB PHOTOMETER tablet** straight from the foil to the water sample and crush the tablet using a clean stirring rod.

6. Close the vial with the cap tightly and swirl the vial several times until the tablet is dissolved.

7. Place the vial in the sample chamber making sure that the Σ marks are aligned.

Zero accepted
prepare Test
press TEST

8. Press **TEST** key.

The result is shown in the display in mg/l PHMB.

1.1 Methods

Notes:

1. Clean vials with the brush after analysis directly.
2. Using vials and stirring rods for a longer time it is possible that they turn blue. In this case clean vials and stirring rods with a laboratory detergent (see chapter 1.2.2 Cleaning of vials and accessories for analysis). Rinse vials and caps thoroughly with tap water and than with deionized water.
3. The test result is influenced by Hardness and Total Alkalinity.
The calibration of this method was done using water of the following concentration:

Ca-Hardness: 200 mg/l CaCO_3

Total Alkalinity: 120 mg/l CaCO_3

1.1 Methods

- 3 2 0 Phosphate, ortho LR with Tablet**
0.05 – 4 mg/l PO₄
Determination of ortho-Phosphate ions
- 3 2 1 Phosphate, ortho HR with Tablet**
1 – 80 mg/l PO₄
Determination of ortho-Phosphate ions
- 3 2 3 Phosphate, ortho with Vario Powder Pack**
0.06 – 2.5 mg/l PO₄
Determination of ortho-Phosphate ions
- 3 2 4 Phosphate, ortho with Vario Tube Test**
0.06 – 5 mg/l PO₄
Determination of ortho-Phosphate ions
- 3 2 7 Phosphat 1, ortho with Vacu-vials®**
5 – 40 mg/l PO₄
Determination of ortho-Phosphate ions
- 3 2 8 Phosphat 2, ortho with Vacu-vials®**
0.05 – 5 mg/l PO₄
Determination of ortho-Phosphate ions
- 3 2 5 Phosphate, acid hydrolizable with Vario Tube Test**
0.02 – 1.6 mg/l P
Determination of ortho-Phosphate ions + condensed, inorganic Phosphates
- 3 2 6 Phosphate, total with Vario Tube Test**
0.02 – 1.1 mg/l P
Determination of ortho-Phosphate ions + condensed, inorganic Phosphates + organically combined Phosphates

1.1 Methods

General:

Ortho-Phosphate ions react with the reagent to a intense blue colour (methods **320**, **323**, **324**, **325** and **326**).

Phosphate in organic and condensed inorganic forms (meta-, pyro- and polyphosphates) must be converted to ortho-Phosphate ions before analysis.

Pretreatment of the sample with acid and heat provides the conditions for hydrolysis of the condensed inorganic forms. Organically combined phosphates are converted to ortho-Phosphate ions by heating with acid and persulfate.

The amount of organically combined phosphates can be calculated:

mg/l Phosphate, organic = mg/l Phosphate, total – mg/l Phosphate, acid hydrolysable

In methods **321** and **327** the ortho-Phosphat ions react with the Vanadate-molybdate-reagent under acid conditions to a yellow coloured product.

Notes - only for tube tests and tests with powder packs:

323, 324, 325, 326

1. Application: for water, wastewater and seawater.
2. Highly buffered samples or samples with extreme pH Values should be adjusted between pH 2 and pH 10 before analysis (with 1 mol/l Hydrochloric acid or 1 mol/l Sodium hydroxide).
3. Interferences:
Large amounts of turbidity may cause inconsistent results.

Interfering substance

Aluminium
Arsenate
Chromium
Copper
Iron
Nickel
Silica (Silicium dioxide)
Silicate
Sulfide
Zinc

Interference level:

greater than 200 mg/l
at any level
greater than 100 mg/l
greater than 10 mg/l
greater than 100 mg/l
greater than 300 mg/l
greater than 50 mg/l
greater than 10 mg/l
at any level
greater than 80 mg/l

Phosphate, ortho \triangleq Phosphorus, reactive

1.1 Methods

3 2 0

Phosphate, ortho LR with Tablet

0.05 – 4 mg/l PO₄



prepare Zero
press ZERO

1. Fill a clean vial (24 mm Ø) with **10 ml of water sample**, close the cap tightly.
2. Place the vial in the sample chamber making sure that the marks \times are aligned.
3. Press **ZERO** key.
4. Remove the vial from the sample chamber.
5. Add **one PHOSPHATE No. 1 LR tablet** straight from the foil to the water sample and crush the tablet using a clean stirring rod.
6. Add **one PHOSPHATE No. 2 LR tablet** straight from the foil to the same water sample and crush the tablet using a clean stirring rod.
7. Close the vial with the cap tightly and swirl the vial several times until the tablet is dissolved.
8. Place the vial in the sample chamber making sure that the marks \times are aligned.

Zero accepted
prepare Test
press TEST

9. Press **TEST** key.
Wait for a **reaction period of 10 minutes**.

Countdown
10:00

After the reaction period is finished the reading starts automatically.

The result is shown in the display in mg/l ortho-Phosphate.

1.1 Methods

Notes:

1. Only ortho-Phosphate ions react.
2. The tablets must be added in the correct sequence.
3. The test sample should have a pH-Value between 6 and 7.
4. Interferences:
Higher concentrations of Cu, Ni, Cr (III), V (V) and W (VI) interfere due to their colour.
Silicates do not interfere (masked by Citric acid in the tablets).
5. see also page 137
6. Conversion:
 $\text{mg/l P} = \text{mg/l PO}_4 \times 0.33$
 $\text{mg/l P}_2\text{O}_5 = \text{mg/l PO}_4 \times 0.75$
7. ▲ PO₄
P
▼ P₂O₅

1.1 Methods



Phosphate HR, ortho with Tablet

1 – 80 mg/l PO₄



prepare Zero
press ZERO

1. Fill a clean vial (24 mm Ø) with **10 ml of water sample**, close the vial with the cap tightly.
2. Place the vial in the sample chamber making sure that the Σ marks are aligned.
3. Press **ZERO** key.
4. Remove the vial from the sample chamber.
5. Add **one PHOSPHATE HR P1 tablet** straight from the foil to the water sample and crush the tablet using a clean stirring rod.
6. Add **one PHOSPHATE HR P2 tablet** straight from the foil to the same water sample and crush the tablet using a clean stirring rod.
7. Close the vial with the cap tightly and swirl the vial several times until the tablets are dissolved.
8. Place the vial in the sample chamber making sure that the Σ marks are aligned.

Zero accepted
prepare Test
press TEST

Count-Down
10:00

9. Press **TEST** key.
Wait for a **reaction period of 10 minutes**.
After the reaction period is finished the reading starts automatically.
The result is shown in the display in mg/l ortho-Phosphate.

1.1 Methods

Notes:

1. For samples under 5 mg/l PO_4 it is recommendable to analyse the water sample with e.g. method 320 "Posphate LR, ortho with Tablet".
2. Only ortho-Phosphate ions react.
3. see also page 137
4. Conversions:
 $\text{mg/l P} = \text{mg/l PO}_4 \times 0.33$
 $\text{mg/l P}_2\text{O}_5 = \text{mg/l PO}_4 \times 0.75$
5. ▲ PO_4
P
▼ P_2O_5

1.1 Methods




Phosphate, ortho with Vario Powder Pack

0.06 – 2.5 mg/l PO₄



prepare Zero
press ZERO



Zero accepted
prepare Test
press TEST

Count-Down
2:00

1. Fill a clean vial (24 mm Ø) with **10 ml of water sample**, close the vial with the cap tightly.

2. Place the vial in the sample chamber making sure that the Σ marks are aligned.

3. Press **ZERO** key.

4. Remove the vial from the sample chamber.

5. Add **VARIO Phosphate Rgt. F10 powder pack** straight from the foil to the water sample.

6. Close the vial with the cap tightly and swirl the vial several times to mix the contents (approx. 10-15 sec., Note 1).

7. Place the vial in the sample chamber making sure that the Σ marks are aligned.

8. Press **TEST** key.

Wait for a **reaction period of 2 minutes**.

After the reaction period is finished the reading starts automatically.

The result is shown in the display in mg/l ortho-Phosphate.

1.1 Methods

Notes:

1. The reagent dissolves not completely.

2. see also page 137

3. Conversions:

$$\text{mg/l P} = \text{mg/l PO}_4 \times 0.33$$

$$\text{mg/l P}_2\text{O}_5 = \text{mg/l PO}_4 \times 0.75$$

4. ▲ PO₄

P

▼ P₂O₅

1.1 Methods




Phosphate, ortho with Vario Tube Test

0.06 – 5 mg/l PO₄



Ø 16 mm

prepare Zero
press ZERO



1. Open the white cap of one **tube PO₄-P Dilution** and add **5 ml water sample**.
2. Place the vial in the sample chamber making sure that the **X** marks are aligned.
3. Press **ZERO** key.
4. Remove the vial from the sample chamber.
5. Add **one VARIO Phosphate Rgt. F10 powder pack** straight from the foil to the water sample (Note 1).
6. Close the vial with the cap tightly and swirl the vial several times to mix the contents (approx. 10-15 sec., Note 2).
7. Place the vial in the sample chamber making sure that the **X** marks are aligned.

Zero accepted
prepare Test
press TEST

Count-Down
2:00

8. Press **TEST** key.
Wait for a **reaction period of 2 minutes**.
After the reaction period is finished the reading starts automatically.
The result is shown in the display in mg/l ortho-Phosphate.

1.1 Methods

Notes:

1. Use a funnel to add the reagent.
2. The reagent dissolves not completely.
3. see also page 137
4. Conversions:
 $\text{mg/l P} = \text{mg/l PO}_4 \times 0.33$
 $\text{mg/l P}_2\text{O}_5 = \text{mg/l PO}_4 \times 0.75$
5. ▲ PO₄
P
▼ P₂O₅

1.1 Methods

3 **2** **7**


Phosphate 1, ortho with Vacu-vials® K-8503 (see Notes)

5 – 40 mg/l PO₄

Insert the adaptor for 13 mm Ø vials.

1. Place the blank in the sample chamber. The blank is part of the test kit.

**prepare Zero
press ZERO**




2. Press **ZERO** key.

3. Remove the blank from the sample chamber.

4. Fill the sampler to the 25 ml mark with the water sample.

5. Place one Vacu-vial® in the sampler. Snap the tip by pressing the vial against the side of the sampler.



The Vacu-vial® breaks at the neck and the vial fills automatically.

A small volume of inert gas remains in the Vacu-vial®.

6. Mix the content of the Vacu-vial® by inverting it several times, allowing the bubble to move from one end to the other. Dry the outside of the vial.

7. Place the Vacu-vial® in the sample chamber.

**Zero accepted
prepare Test
press TEST**

8. Press **TEST** key.
Wait for a **reaction period of 5 minutes**.

**Count-Down
5:00**

After the reaction period is finished the reading starts automatically.

The result is shown in the display in mg/l ortho-Phosphate.

1.1 Methods

Notes:

1. This method is adapted from CHEMetrics.
2. Read the original test instruction and the MSDS (delivered with the test) before performing the test. MSDS also at www.chemetrics.com available.
3. Vacu-vials® is a registered trade mark of the company CHEMetrics, Inc. / Calverton, U.S.A.
4. Only ortho-Phosphate ions react.
5. Sulfide, Thiosulfate und Thiocyanate cause low test results.
6. ▲ PO_4
P
▼ P_2O_5

1.1 Methods

3 **2** **8**

Phosphate 2, ortho with Vacu-vials® K-8513 (see Notes)


0.05 – 5 mg/l PO₄

Insert the adaptor for 13 mm Ø vials.

1. Place the blank in the sample chamber. The blank is part of the test kit.

prepare Zero
press ZERO

2. Press **ZERO** key.



3. Remove the blank from the sample chamber.

4. Fill the sampler to the 25 ml mark with the water sample.

5. Fill the sampler with drops of the same size by holding the bottle vertically and squeeze slowly:




2 drops A-8500 Activator Solution

6. Close the sampler with the cap tightly and swirl several times to mix the contents.



7. Place one Vacu-vial® in the sampler. Snap the tip by pressing the vial against the side of the sampler. The Vacu-vial® breaks at the neck and the vial fills automatically. A small volume of inert gas remains in the Vacu-vial®.



8. Mix the content of the Vacu-vial® by inverting it several times, allowing the bubble to move from one end to the other. Dry the outside of the vial.

9. Place the Vacu-vial® in the sample chamber.

Zero accepted
prepare Test
press TEST

10. Press **TEST** key.

Wait for a **reaction period of 3 minutes**.

Count-Down
3:00

After the reaction period is finished the reading starts automatically.

The result is shown in the display in mg/l ortho-Phosphate.

1.1 Methods

Notes:


1. This method is adapted from CHEMetrics.
2. Read the original test instruction and the MSDS (delivered with the test) before performing the test. MSDS also at www.chemetrics.com available.
3. Vacu-vials® is a registered trade mark of the company CHEMetrics, Inc. / Calverton, U.S.A.
4. Only ortho-Phosphate ions react.
5. Sulfide, Thiosulfate und Thiocyanate cause low test results.
6. ▲ PO_4
P
▼ P_2O_5

1.1 Methods

3 **2** **5**


Phosphate, acid hydrolyzable with Vario Tube Test

0.02 – 1.6 mg/l P (≙ 0.06 – 5 mg/l PO₄)



1. Open the white cap of one **digestion tube PO4-P Acid reagent** and add **5 ml water sample**.
2. Close the vial with the cap tightly. Invert the vial gently several times to mix the contents.
3. Heat the vials for **30 minutes** in the preheated reactor at a temperature of **100°C**.
4. After 30 minutes remove the vial from the reactor. **(CAUTION: The vials are hot!)** Allow the vials to cool to room temperature.
5. Open the cooled down digestion vial and add **2 ml 1.00 N Sodium hydroxide solution** to the vial.
6. Close the vial with the cap and invert the vial gently several times to mix the contents.
7. Place the vial in the sample chamber making sure that the Δ marks are aligned.
8. Press **ZERO** key.
9. Remove the vial from the sample chamber.
10. Add **one VARIO Phosphate Rgt. F10 powder pack** straight from the foil to the vial (Note 2).
11. Close the vial with the cap tightly and swirl the vial several times to mix the contents (approx. 10-15 sec., Note 3).
12. Place the vial in the sample chamber making sure that the Δ marks are aligned.
13. Press **TEST** key.

prepare Zero
press ZERO



Zero accepted
prepare Test
press TEST

Count-Down
2:00

Wait for a **reaction period of 2 minutes**.

After the reaction period is finished the reading starts automatically.

The result is shown in the display in mg/l acid hydrolyzable Phosphate.

1.1 Methods

Notes:


1. Use appropriate safety precautions and a good lab technique should be used during the whole procedure.
2. Use a funnel to add the reagent.
3. The reagent dissolves not completely.
4. see also page 137
5. Conversions:
 $\text{mg/l PO}_4 = \text{mg/l P} \times 3.07$
 $\text{mg/l P}_2\text{O}_5 = \text{mg/l P} \times 2.29$
6. ▲ PO_4
P
▼ P_2O_5

1.1 Methods

3 **2** **6**


Phosphate, total with Vario Tube Test

0.02 – 1.1 mg/l P (≙ 0.06 – 3.5 mg/l PO₄)



1. Open the white cap of one **digestion tube PO4-P Acid reagent** and add **5 ml water sample**.
2. Add **one Vario Potassium Persulfate F10 powder pack** straight from the foil to the vial (Note 2).
3. Close the vial with the cap tightly. Invert the vial several times to mix the contents.
4. Heat the vials for **30 minutes** in the preheated reactor at a temperature **of 100°C**.
5. After 30 minutes remove the vial from the reactor. **(CAUTION: The vials are hot!)** Allow the vials to cool to room temperature.
6. Open the cooled down digestion vial and add **2 ml 1.54 N Sodium hydroxide solution** to the vial.
7. Close the vial with the cap and invert the vial gently several times to mix the contents.
8. Place the vial in the sample chamber making sure that the Δ marks are aligned.
9. Press **ZERO** key.
10. Remove the vial from the sample chamber.
11. Add **one VARIO Phosphate Rgt. F10 powder pack** straight from the foil to the vial (Note 2).
12. Close the vial with the cap tightly and swirl the vial several times to mix the contents (approx. 10-15 sec., Note 3).
13. Place the vial in the sample chamber making sure that the Δ marks are aligned.
14. Press **TEST** key.
Wait for a **reaction period of 2 minutes**.
After the reaction period is finished the reading starts automatically.
The result is shown in the display in mg/l total Phosphate.

**prepare Zero
press ZERO**



**Zero accepted
prepare Test
press TEST**

**Count-Down
2:00**

1.1 Methods

Notes:

1. Use appropriate safety precautions and a good lab technique should be used during the whole procedure.
2. Use a funnel to add the reagent.
3. The reagent dissolves not completely.
4. see also page 137
5. Conversions:
 $\text{mg/l PO}_4 = \text{mg/l P} \times 3.07$
 $\text{mg/l P}_2\text{O}_5 = \text{mg/l P} \times 2.29$
6. ▲ P
 PO₄
 ▼ P₂O₅

1.1 Methods



pH-value LR 5.2 – 6.8 with Tablet



prepare Zero
press ZERO

1. Fill a clean vial (24 mm Ø) with **10 ml of water sample**, close the vial with the cap tightly.
2. Place the vial in the sample chamber making sure that the Σ marks are aligned.
3. Press **ZERO** key.
4. Remove the vial from the sample chamber.
5. Add **one BROMOCRESOLPURPLE PHOTOMETER tablet** straight from the foil to the water sample and crush the tablet using a clean stirring rod.
6. Close the vial with the cap tightly and swirl the vial several times until the tablet is dissolved.
7. Place the vial in the sample chamber making sure that the Σ marks are aligned.

Zero accepted
prepare Test
press TEST

8. Press **TEST** key.

The result is shown in the display as pH-value.

1.1 Methods

Notes:

1. For photometric determination of pH-Values only use BROMOCRESOLPURPLE tablets in black printed foil pack and marked with PHOTOMETER.
2. pH-Values below 5.2 and above 6.8 can produce results inside the measuring range. A plausibility test (pH-meter) is recommend.
3. The accuracy of the colorimetric determination of pH-values depends on various boundary conditions (buffer capacity of the sample, salt content etc.).
4. Salt error
Correction of test results (average values) for samples with salt content of:

Indicator	Salt content		
Bromcresolpurple	1 molar - 0.26	2 molar - 0.33	3 molar - 0.31

The values of Parson and Douglas (1926) are based on the use of Clark and Lubs buffers.

1 Mol NaCl = 58.4 g/l = 5.8 %

1.1 Methods

3 3 0

pH-value 6.5 – 8.4 with Tablet



1. Fill a clean vial (24 mm Ø) with **10 ml of water sample**, close the vial with the cap tightly.
2. Place the vial in the sample chamber making sure that the Σ marks are aligned.

prepare Zero
press ZERO

3. Press **ZERO** key.
4. Remove the vial from the sample chamber.
5. Add **one PHENOL RED PHOTOMETER tablet** straight from the foil to the water sample and crush the tablet using a clean stirring rod.
6. Close the vial with the cap tightly and swirl the vial several times until the tablet is dissolved.
7. Place the vial in the sample chamber making sure that the Σ marks are aligned.

Zero accepted
prepare Test
press TEST

8. Press **TEST** key.

The result is shown in the display as pH-value.

1.1 Methods

Notes:

1. For photometric determination of pH-values only use PHENOL RED tablets in black printed foil pack and marked with PHOTOMETER.
2. Water samples with low values of Alkalinity-m (below 35 mg/l CaCO_3) may give wrong pH readings.
3. pH-values below 6.5 and above 8.4 can produce results inside the measuring range. A plausibility test (pH-meter) is recommended.
4. The accuracy of the colorimetric determination of pH-values depends on various boundary conditions (buffer capacity of the sample, salt content etc.).
5. Salt error

Correction of test results (average values) for samples with salt content of:

Indicator	Salt content		
Phenol red	1 molar - 0.26	2 molar - 0.33	3 molar - 0.31

The values of Parson and Douglas (1926) are based on the use of Clark and Lubs buffers. 1 Mol NaCl = 58.4 g/l = 5.8 %

1.1 Methods



pH-value 6.5 – 8.4 with Liquid Reagent



prepare Zero
press ZERO

1. Fill a clean vial (24 mm Ø) with **10 ml of water sample**, close the vial with the cap tightly.
2. Place the vial in the sample chamber making sure that the Σ marks are aligned.

3. Press **ZERO** key.

4. Remove the vial from the sample chamber.

5. Fill the vial with drops of the same size by holding the bottle vertically and squeeze slowly:

6 drops of PHENOL RED solution

6. Close the vial with the cap tightly and swirl the vial several times to mix the contents.

7. Place the vial in the sample chamber making sure that the Σ marks are aligned.

Zero accepted
prepare TEST
press Test

8. Press **TEST** key.

The result is shown in the display as pH-value.

1.1 Methods

Notes:

1. When testing chlorinated water the residual chlorine content can influence the colour reaction of the liquid reagent. This can be avoided (without interfering the pH measurement) by adding a small crystal of Sodiumthiosulfate ($\text{Na}_2\text{S}_2\text{O}_3 \times 5 \text{H}_2\text{O}$) to the sample before adding the PHENOL RED solution. PHENOL RED tablets already contain Thiosulfate.
2. Due to differing drop size results can show a discrepancy in accuracy by comparison with tablets. This can be minimised by using a pipette (0.18 ml PHENOLRED solution is equivalent to 6 drops).
3. After use replace the bottle cap securely.
4. **Store the reagent in a cool, dry place ideally at between 6°C and 10°C.**

1.1 Methods



pH-value HR 8.0 – 9.6 with Tablet



prepare Zero
press ZERO

1. Fill a clean vial (24 mm Ø) with **10 ml of water sample**, close the vial with the cap tightly.
2. Place the vial in the sample chamber making sure that the \times marks are aligned.
3. Press **ZERO** key.
4. Remove the vial from the sample chamber.
5. Add **one THYMOLBLUE PHOTOMETER tablet** straight from the foil to the water sample and crush the tablet using a clean stirring rod.
6. Close the vial with the cap tightly and swirl the vial several times until the tablet is dissolved.
7. Place the vial in the sample chamber making sure that the \times marks are aligned.

Zero accepted
prepare TEST
press Test

8. Press **TEST** key.

The result is shown in the display as pH-value.

1.1 Methods

Notes:

1. For photometric determination of pH-values only use THYMOLBLUE tablets in black printed foil pack and marked with PHOTOMETER.
2. pH-Values below 8.0 and above 9.6 can produce results inside the measuring range. A plausibility test (pH-meter) is recommend.
3. The accuracy of the colorimetric determination of pH-values depends on various boundary conditions (buffer capacity of the sample, salt content etc.).
4. Salt error

Correction of test results (average values) for samples with salt content of:

Indicator	Salt content		
Thymolblue	1 molar - 0.22	2 molar - 0.29	3 molar - 0.34

The values of Parson and Douglas (1926) are based on the use of Clark and Lubs buffers. 1 Mol NaCl = 58.4 g/l = 5.8 %

1.1 Methods

3 4 0

Potassium with Tablet

0.7 – 12 mg/l K



prepare Zero
press ZERO

1. Fill a clean vial (24 mm Ø) with **10 ml of water sample**, close the vial with the cap tightly.
2. Place the vial in the sample chamber making sure that the Σ marks are aligned.
3. Press **ZERO** key.
4. Remove the vial from the sample chamber.
5. Add **one Potassium T tablet** straight from the foil to the water sample and crush the tablet using a clean stirring rod.
6. Close the vial with the cap tightly and swirl the vial several times until the tablet is dissolved.
7. Place the vial in the sample chamber making sure that the Σ marks are aligned.

Zero accepted
prepare Test
press TEST

8. Press **TEST** key.

The result is shown in the display in mg/l Potassium.

1.1 Methods

Notes:

1. If Potassium is present a cloudy solution will be given.
Single particles are uncaused necessarily by Potassium.

1.1 Methods

3 5 0

Silica/Silicon dioxide with Tablet

0.05 – 4 mg/l SiO₂



Ø 24 mm

prepare Zero
press ZERO

Countdown
5:00
start: ↓

Zero accepted
prepare Test
press TEST

Countdown
1:00

1. Fill a clean vial (24 mm Ø) with **10 ml of water sample**, close the vial with the cap tightly.
2. Place the vial in the sample chamber making sure that the X marks are aligned.
3. Press **ZERO** key.
4. Remove the vial from the sample chamber.
5. Add **one SILICA No. 1 tablet** straight from the foil to the water sample and crush the tablet using a clean stirring rod.
6. Close the vial with the cap tightly and swirl the vial several times until the tablet is dissolved.
7. Press [↓] key.
Wait for a **reaction period of 5 minutes**.

After reaction period is finished proceed as follows:

8. Add **one SILICA PR tablet** straight from the foil to the same water sample and crush the tablet using a clean stirring rod.
9. Add **one SILICA No. 2 tablet** straight from the foil to the same water sample and crush the tablet using a clean stirring rod.
10. Close the cap tightly and swirl the vial several times until the tablets are dissolved.
11. Place the vial in the sample chamber making sure that the X marks are aligned.
12. Press **TEST** key.
Wait for a **reaction period of 1 minute**.

After the reaction period is finished the reading starts automatically.

The result is shown in the display in mg/l Silica.

1.1 Methods

Notes:

1. The tablets must be added in the correct sequence.
2. Phosphate ions does not interfere under the given reaction conditions.
3. If Phosphate is known to be absent, the addition of the SILICA PR tablet may be omitted.
4. Conversion:
 $\text{mg/l Si} = \text{mg/l SiO}_2 \times 0.47$
5. ▲ SiO₂
▼ Si

1.1 Methods

3 5 1

Silica LR / Silicon dioxide LR with Vario Powder Pack and Liquid Reagent

0.1 – 1.6 mg/l SiO₂




Use two clean vials (24 mm Ø) and mark one as blank for zeroing.

1. Fill each vial with **10 ml water sample**.
2. Add **0.5 ml Vario Molybdate 3 reagent solution** into each vial.
3. Close the vials with the caps tightly and swirl the vials several times to mix the contents (Note 1).
4. Press [**↵**] key.

Countdown
4:00
start: ↵

Wait for a **reaction period of 4 minutes** (Note 2).



After reaction period is finished proceed as follows:

5. Add **one Vario Silica Citric Acid F10 Powder Pack** straight from the foil into each vial.
6. Close the vials with the caps tightly and swirl the vials several times to mix the contents.
7. Press [**↵**] key.

Countdown
1:00
start: ↵

Wait for a **reaction period of 1 minute** (Note 3).

After reaction period is finished proceed as follows:

8. Place the vial (the blank) in the sample chamber making sure that the **X** marks are aligned.
9. Add **one Vario LR Silica Amino Acid F F10 Powder Pack** straight from the foil into the vial (the sample).
10. Close the vial with the cap tightly and swirl the vial several times to mix the contents.

1.1 Methods

**prepare Zero
press ZERO**

11. Press **ZERO** key (blank is already placed in the sample chamber - see point 8).

**Count-Down
2:00**

Wait for a **reaction period of 2 minutes**.

After the reaction period is finished the zero-reading starts automatically.

12. Remove the vial from the sample chamber.
13. Place the vial (the sample) in the sample chamber making sure that the Σ marks are aligned.

**Zero accepted
prepare Test
press TEST**

14. Press **TEST** key.

The result is shown in the display in mg/l Silica.

Notes:

1. Close the vials with the cap directly after adding the Vario Molybdate 3 reagent solution, otherwise it can result in minimum findings.
2. The given reaction time of 4 minutes refers to a water sample temperature of 20°C.
At 30°C a reaction time of 2 minutes, at 10°C a reaction time of 8 minutes is required.
3. The given reaction time of 1 minute refers to a water sample temperature of 20°C.
At 30°C a reaction time of 30 seconds, at 10°C a reaction time of 2 minutes is required.
4. Interferences:

Substance	Interference
Iron	large amounts interfere
Phosphate	does not interfere at concentrations less than 50 mg/l PO ₄ at 60 mg/l PO ₄ the interference is approx. – 2% at 75 mg/l PO ₄ the interference is approx. – 11%
Sulfide	interfere at all levels

Occasionally water samples contain silica forms which reacts very slowly with Molybdate. The nature of these forms is not known.

A pre-treatment with Sodium hydrogencarbonate and then with Sulfuric Acid will make these forms reactive to Molybdate (pre-treatment is given in "Standard Methods for the Examination of Water and Wastewater" under "Silica-Digestion with Sodium Bicarbonate").

5. ▲ SiO₂
▼ Si

1.1 Methods


3 5 2

Silica HR / Silicon dioxide HR with Vario Powder Pack

1 – 90 mg/l SiO₂



prepare Zero
press ZERO



Countdown
10:00
start: ↵

Zero accepted
prepare Test
press TEST

Count-Down
2:00

1. Fill a clean vial (24 mm Ø) with **10 ml of water sample** (Note 1), close the vial with the cap tightly.
2. Place the vial in the sample chamber making sure that the X marks are aligned.
3. Press **ZERO** key.
4. Remove the vial from the sample chamber.
5. Add **one Vario Silica Molybdate F10 powder pack** straight from the foil to the water sample.
6. Close the vial with the cap tightly and swirl the vial several times to mix the contents.
7. Add **one Vario Silica HR Acid Rgt. F10 powder pack** straight from the foil to the same water sample (Note 2).
8. Close the vial with the caps tightly and swirl the vials several times to mix the contents.
9. Press [↵] key.

Wait for a **reaction period of 10 minutes**.

After reaction period is finished proceed as follows:

10. Add **one Vario Silica Citric Acid F10 powder pack** straight from the foil to the water sample (Note 3).
11. Close the vial with the cap tightly and swirl the vial several times to mix the contents.
12. Place the vial in the sample chamber making sure that the X marks are aligned.
13. Press **TEST** key.

Wait for a **reaction period of 2 minutes**.

After reaction period is finished the reading starts automatically.

The result is shown in the display in mg/l Silica.

1.1 Methods

Notes:

1. Temperature of the sample should be 15°C – 25°C.
2. If Silica or Phosphate is present a yellow colour is developed.
3. In this step any yellow colour due to Phosphate is removed.
4. Interferences:

Substance	Interference
Iron	large amounts interfere
Phosphate	does not interfere at concentrations less than 50 mg/l PO ₄ at 60 mg/l PO ₄ the interference is approx. – 2% at 75 mg/l PO ₄ the interference is approx. – 11 %
Sulfide	interfere at all levels

Occasionally water samples contain silica forms which reacts very slowly with Molybdate. The nature of these forms is not known.

A pre-treatment with Sodium hydrogencarbonate and then with Sulfuric Acid will make these forms reactive to Molybdate (pre-treatment is given in "Standard Methods for the Examination of Water and Wastewater" under "Silica-Digestion with Sodium Bicarbonate").


5. ▲ SiO₂
▼ Si

1.1 Methods




Sodium hypochlorite (Soda bleaching lye) with Tablet

0.2 – 16 % w/w NaOCl



Preparation:

1. Fill a 5 ml plastic syringe with the test solution, ensuring that all air bubbles are expelled. Fill the 5 ml test solution slowly into a 100 ml beaker and dilute to the 100 ml mark with chlorine-free water. Mix thoroughly.
2. Fill a 5 ml plastic syringe with the diluted test solution (step 1) to the 1 ml mark, ensuring that all air bubbles are expelled. Fill the 1 ml test solution slowly into a 100 ml beaker and dilute to the 100 ml mark with chlorine-free water. Mix thoroughly.



Ø 24 mm

Performing test procedure:

1. Fill a clean vial (24 mm Ø) with **10 ml of the prepared water sample**, close the vial with the cap tightly.
2. Place the vial in the sample chamber making sure that the Σ marks are aligned.
3. Press **ZERO** key.
4. Remove the vial from the sample chamber.
5. Add **one CHLORINE HR (KI) tablet** straight from the foil to the water sample and crush the tablet using a clean stirring rod.
6. Add **one ACIDIFYING GP tablet** straight from the foil to the same water sample and crush the tablet using a clean stirring rod.
7. Close the vial with the cap tightly and swirl the vial several times until the tablets are dissolved.

prepare Zero
press ZERO

1.1 Methods

8. Place the vial in the sample chamber making sure that the \times marks are aligned.

**Zero accepted
prepare Test
press TEST**

9. Press **TEST** key.

The result is shown in the display in % w/w as available chlorine present in the original sample of Sodium hypochlorite.

Notes:

1. Please pay attention by handling with sodium hypochlorite. The material has a very strong alkalinity and can cause corruptions. The contact with eyes, skin and clothes etc. has to be avoided. It is necessary to look at the detailed information the producer has given about the product.
2. The tablets must be added in the correct sequence.
3. This method gives you the opportunity of a fast and simple test. The test can be arranged on the premises but the result will not give you a detailed specification like a laboratory method.
4. By following the strict order of procedure an exactness of +/- 1 weight % can be reached.

1.1 Methods


3

5

5

Sulfate with Tablet

5 – 100 mg/l SO₄



1. Fill a clean vial (24 mm Ø) with **10 ml of water sample**, close the vial with the cap tightly.

2. Place the vial in the sample chamber making sure that the Σ marks are aligned.

**prepare Zero
press ZERO**

3. Press **ZERO** key.

4. Remove the vial from the sample chamber.

5. Add **one SULFATE T tablet** straight from the foil to the water sample and crush the tablet using a clean stirring rod.

6. Close the vial with the cap tightly and swirl the vial several times until the tablet is dissolved.

7. Place the vial in the sample chamber making sure that the Σ marks are aligned.

**Zero accepted
prepare Test
press TEST**

8. Press **TEST** key.

The result is shown in the display in mg/l Sulfate.

1.1 Methods

Notes:

1. If Sulfate is present a cloudy solution will be given.

1.1 Methods

3 6 0

Sulfate with Vario Powder Pack

5 – 100 mg/l SO₄




1. Fill a clean vial (24 mm Ø) with **10 ml of water sample**, close the vial with the cap tightly.

2. Place the vial in the sample chamber making sure that the \times marks are aligned.

prepare Zero
press ZERO

3. Press **ZERO** key.



4. Remove the vial from the sample chamber.

5. Add **one VARIO Sulpha 4 / F10** powder pack straight from the foil to the water sample.

6. Close the vial with the cap tightly and swirl the vial several times to mix the contents.

7. Place the vial in the sample chamber making sure that the \times marks are aligned.

Zero accepted
prepare Test
press TEST

8. Press **TEST** key.
Wait for a **reaction period of 5 minutes**.

Countdown
5:00

After the reaction period is finished the reading starts automatically.

The result is shown in the display in mg/l Sulfate.

1.1 Methods

Note:

1. If Sulfate ions are present a cloudy solution will be given.

1.1 Methods



Sulfide with Tablet

0.04 – 0.5 mg/l S



Ø 24 mm

1. Fill a clean vial (24 mm Ø) with **10 ml of water sample**, close the vial with the cap tightly.

2. Place the vial in the sample chamber making sure that the Σ marks are aligned.

**prepare Zero
press ZERO**

3. Press **ZERO** key.

4. Remove the vial from the sample chamber.

5. Add **one SULFIDE No. 1 tablet** to the water sample and crush the tablet using a clean stirring rod and dissolve the tablet.

6. Add **one SULFIDE No. 2 tablet** to the same water sample and crush the tablet using a clean stirring rod.

7. Close the vial with the cap tightly and swirl the vial several times until the tablets are dissolved.

8. Place the vial in the sample chamber making sure that the Σ marks are aligned.

**Zero accepted
prepare Test
press TEST**

9. Press **TEST** key.

Wait for a **reaction period of 10 minutes**.

**Countdown
10:00**

After the reaction period is finished the reading starts automatically.

The result is shown in the display in mg/l Sulfide.

1.1 Methods

Notes:

1. The tablets must be added in the correct sequence.
2. Chlorine and other oxidizing agents which react with DPD do not interfere in the test.
3. To avoid loss of Sulfide collect the sample carefully with a minimum of aeration. It is essential to test the sample immediately after collection.
4. The temperature of test performance should be 20°C. Difference to this temperature can lead to higher or lower results.
5. Conversion:
$$\text{H}_2\text{S} = \text{mg/l S} \times 1.06$$
6. ▲ S
▼ H₂S

1.1 Methods

3 7 0

Sulfite with Tablet

0.1 – 5 mg/l SO₃



1. Fill a clean vial (24 mm Ø) with **10 ml of water sample**, close the vial with the cap tightly.

2. Place the vial in the sample chamber making sure that the X marks are aligned.

prepare Zero
press ZERO

3. Press **ZERO** key.

4. Remove the vial from the sample chamber.

5. Add **one SULFITE LR tablet** straight from the foil to the water sample and crush the tablet using a clean stirring rod.

6. Close the vial with the cap tightly and swirl the vial several times until the tablet is dissolved.

7. Place the vial in the sample chamber making sure that the X marks are aligned.

Zero accepted
prepare Test
press TEST

8. Press **TEST** key.

Wait for a **reaction period of 5 minutes**.

Count-Down
5:00

After the reaction period is finished the reading starts automatically.

The result is shown in the display in mg/l Sulfite.

1.1 Methods

Notes:

1. ▲ SO_3
▼ Na_2SO_3

1.1 Methods

3 9 0

Urea with Tablet and Liquid Reagent

0.1 – 3 mg/l (NH₂)₂CO / mg/l Urea



prepare Zero
press ZERO

1. Fill a clean vial (24 mm Ø) with **10 ml of water sample**, close the vial with the cap tightly.
2. Place the vial in the sample chamber making sure that the X marks are aligned.
3. Press **ZERO** key.

4. Remove the vial from the sample chamber.

5. Add **2 drops of Urea reagent 1** to the water sample (Note 8).

6. Close the vial with the cap tightly and swirl the vial several times to mix the contents.

Countdown
5:00
start: ↓

7. Add **1 drop of Urea Reagent 2** (Urease) to the same water sample (Note 8).

8. Close the vial with the cap tightly and swirl the vial several times to mix the contents.

9. Press [**↓**] key.

Wait for a **reaction period of 5 minutes**.

After reaction period is finished proceed as follows:

10. Add **one AMMONIA No. 1 tablet** straight from the foil to the prepared water sample and mix to dissolve with a clean stirring rod.

11. Add **one AMMONIA No. 2 tablet** straight from the foil to the same water sample and mix to dissolve with a clean stirring rod.

1.1 Methods

12. Close the vial with the cap tightly and swirl the vial several times until the tablets are dissolved.

13. Place the vial in the sample chamber making sure that the X marks are aligned.

Zero accepted
prepare Test
press TEST

14. Press **TEST** key.
Wait for a **reaction period of 10 minutes**.

Countdown
10:00

After the reaction period is finished the reading starts automatically.

The result is shown in the display in mg/l Urea.

Notes:

1. The sample temperature should be between 20°C and 30°C.
2. Determination at the latest one hour after sample taking.
3. The tablets must be added in the correct sequence.
4. **Store reagent 2 (Urease) in the refrigerator at a temperature of 4°C to 8°C.**
5. The AMMONIA No. 1 tablet will only dissolve completely after the AMMONIA No. 2 tablet has been added.
6. Ammonium and chloramines are also measured during urea measurement.
7. Before analysing seawater samples, a measuring spoon of Ammonia Conditioning Powder must be added to the sample and swirled to dissolve before AMMONIA No. 1 tablet is added.
8. Fill the vial with drops of the same size by holding the bottle vertically and squeeze slowly.

1.1 Methods

4 0 0

Zinc with Tablet

0.02 – 1 mg/l Zn



Ø 24 mm

1. Fill a clean vial (24 mm Ø) with **10 ml of water sample**.
2. Add **one COPPER / ZINC LR tablet** straight from the foil to the water sample, crush the tablet using a clean stirring rod.
3. Close the vial with the cap tightly and swirl the vial several times until the tablet is dissolved.
4. Place the vial in the sample chamber making sure that the X marks are aligned.

prepare Zero
press ZERO

Count-Down
5:00

5. Press **ZERO** key.
Wait for a **reaction period of 5 minutes**.
After the reaction period is finished the reading starts automatically.
6. Remove the vial from the sample chamber.
7. Add **one EDTA tablet** straight from the foil to the prepared vial and crush the tablet using a clean stirring rod.
8. Close the vial with the cap tightly and swirl the vial several times until the tablet is dissolved.
9. Place the vial in the sample chamber making sure that the X marks are aligned.

Zero accepted
press ZERO
press TEST

10. Press **TEST** key.
The result is shown in the display in mg/l Zinc.

1.1 Methods

Notes:

1. The tablets must be added in the correct sequence.
2. In the case of high levels of residual chlorine, perform the analysis with a dechlorinated water sample. To dechlorinate add one DECHLOR tablet to the water sample (point 1). Crush and mix to dissolve the tablet. Then add the COPPER / ZINC LR tablet (point 2) and continue with the test procedure as described above.

1.2 Important notes

1.2.1 Correct use of reagents

The reagents must be added in the correct sequence.

Tablet reagents:

The tablet reagents should be added to the water sample straight from the foil without touching them with the fingers.


Liquid reagents:

Add drops of the same size to the water sample by holding the bottle vertically and squeezing slowly.

After use replace the bottle caps securely noting the colour coding.

Note recommendation for storage (e.g. cool and dry).


Powder Packs:



Vacu-vials® of CHEMetrics:

Vacu-vials® should be stored dark and at room temperature.

For further information see MSDS.



1.2.2 Cleaning of vials and accessories for analysis

Vials, caps and stirring rods should be cleaned thoroughly **after each analysis** to prevent influences.

Procedure:


Clean vials and accessories after each analysis as soon as possible.

- a. Clean vials and accessories with laboratory detergent (e.g. Extran® MA 02 (neutral, phosphatic), Extran® MA 03 (alkaline, phosphate-free) from Merck KGaA).
- b. Rinse with tap water thoroughly.
- c. On demand (see Notes) perform special cleaning at this point, e.g.: rinse with diluted Hydrochloric acid solution.
- d. Rinse with deionized water thoroughly.


1.2.3 Guidelines for photometric measurements

1. Vials, caps and stirring rods should be cleaned thoroughly after each analysis to prevent influences. Even minor reagent residues can cause errors in the test result.
2. The outside of the vial must be clean and dry before starting the analysis. Clean the outside of the vials with a towel. Fingerprints or other marks will be removed.
3. If there is no defined vial for the blank, the zeroing and the test must be carried out with the same vial as there may be slight differences in optical performance between vials.
4. The vials must be positioned in the sample chamber for zeroing and test with the Δ mark on the vial aligned with the ∇ mark on the instrument.

Correct position of the vial (Ø 24 mm):




Correct position of the vial (Ø 16 mm):



5. Always perform zeroing and test with closed vial cap. Only use cap with sealing ring.
6. Bubbles on the inside wall of the vial leads to incorrect measurements. To prevent this, remove the bubbles by swirling the vial before performing the test.
7. Avoid spillage of water in the sample chamber. If water should leak into the instrument housing, it can destroy electronic components and cause corrosion.
8. Contamination of the lens in the sample chamber can result in errors. Check at regular intervals and – if necessary – clean the light entry surfaces of the sample chamber using a moist cloth or cotton buds.
9. Large temperature differences between the instrument and the environment can lead to errors – e.g. due to the formation of condensation in the area of the lens or on the vial.
10. To avoid errors caused by stray-light do not use the instrument in bright sunlight.

Correct filling of the vial:



1.2.4 Sample dilution techniques

Proceed as follows for accurate dilutions:

Pipette the water sample (see table) into a 100 ml volumetric flask and fill up to 100 ml mark with deionized water. Swirl to mix the contents.

Water sample [ml]	Multiplication factor
1	100
2	50
5	20
10	10
25	4
50	2

Pipette the required volume of the diluted sample into the vial and proceed as described in the test methods.

Caution:

1. Dilution decreases accuracy.
2. Do not dilute water samples for measurement of pH-values. This will lead to incorrect test results. If there is displayed "Overrange" use another instrument (e.g. pH-meter).

1.2.5 Correcting for volume additions

If a larger volume of acid or base is used to pre-adjust the pH-value, a volume correction of the displayed result is necessary.

Example:

For adjusting the pH-value of a 100 ml water sample 5 ml of acid had to be added. The corresponding displayed result is 10 mg/l.

Total volume = 100 ml + 5 ml = 105 ml

Correction factor = 105 ml / 100 ml = 1.05

Corrected result = 10 mg/l x 1.05 = 10.5 mg/l

Part 2

Operating manual

2.1 Operation

2.1.1 Commissioning

Before working with the photometer insert the rechargeable batteries and the Lithium battery (content of delivery). The rechargeable batteries are not charged. See chapter 2.1.2 Saving data – Important Notes, 2.1.3 Replacement of rechargeable batteries resp. Lithium battery, and 2.1.4 Charging the rechargeable batteries.

Before using the photometer select language (mode 10), select mode 34 and perform “Delete Data”. Set date and time (see chapter 2.4 Photometer settings).

2.1.2 Saving data – Important Notes

The Lithium battery saves data (stored results and photometer setting) if there is no power from the power supply from the rechargeable batteries or the mains adapter.

Recommendation: Exchange of the lithium battery every 5 years.

Note: When neither mains adapter nor batteries supply energy to the instrument, all stored data and settings will be lost, if the lithium battery is taken out.

Recommendation: Keep the instrument connected to mains adapter supply while changing the lithium battery.

2.1.3 Replacement of rechargeable batteries resp. Lithium-battery

1. Switch the instrument off.
2. If necessary remove vial from the sample chamber.
3. Place the instrument upside down on a clean and even surface.
4. Unscrew the two screws (A) of the battery compartment cover (B).
5. Lift battery compartment cover off.
6. If necessary remove old rechargeable batteries (C) and/or the Lithium-battery (D) (See 2.1.4).
7. Place 7 new rechargeable batteries and/or the Lithium-battery.
Ensuring the correct polarity!
8. Replace the battery compartment cover.
9. Tighten the screws carefully.

CAUTION

Dispose of used rechargeable batteries and Lithium-batteries in accordance with all federal, state and local regulations.

2.1.4 Charging the rechargeable batteries

The rechargeable batteries are uncharged in the instrument. As soon as the photometer is connected with the mains adapter to the mains the rechargeable batteries are charged. Empty rechargeable batteries should be charged in the instrument for at least 5 days. 10 charging and discharging cycles are necessary before the rechargeable batteries obtain their full capacity.


2.1.5 Fuse

The instrument contains a fuse (E) (type: 1 A, inert, 20 mm). If a replacement is necessary proceed as described in "Replacement of rechargeable batteries resp. Lithium-battery". If the instrument can be operated with the mains adapter but not with the rechargeable batteries, the fuse could be defect (try new rechargeable batteries first).

2.1.6 Protective caps:

If not used protect the two connections against damage (e.g. corrosion) caused by environmental influences (e.g. dust or splashing) keep the protective caps in place (G).

- (A) screws
- (B) battery compartment cover
- (C) rechargeable batteries 7 Ni-MH-rechargeable batteries (Type AA, 1100 mAh)
- (D) battery Lithium-battery (Type CR 2032, 3V)
- (E) fuse 1 A, inert, 20 mm
- (F) instrument















2.2 Overview of function keys




Attention:

With the software-update V012.002.3.003.001 an "ESC-function" is implemented. If your keypad doesn't show an [Esc]-key please note that the grey key without a print (lowest key on the left) has the "ESC-function".

2.2.1 Overview

	Switching the photometer on or off
	Returning to selection of methods or previous menu
	Function key: description in the text if key available
	Function key: description in the text if key available
	Function key: description in the text if key available
	Confirming
	Menu of photometer settings and further functions
	Moving the cursor ">>" up resp. down
	Storing of displayed test result
	Performing Zero
	Performing Test
	Displaying date and time / user-countdown

2.2.2 Displaying time and date:

	Press ["clock"] key.
	The display shows: After 15 seconds the photometer reverts to the previous display automatically or press [] key or [Esc].

2.2.3 User-countdown

With this function the operator is able to define his own countdown.



Press [“clock”] key.

19.30.20 2006-06-15

The display shows time and date:



Press [“clock”] key.

Count-Down
mm : ss
99 : 99

The display shows:

Either press [↵] key to accept the last used user-countdown.

or


press any number key to start entering a new value

The entering comprises two digits each.

Enter minutes and seconds,

e.g.: 2 minutes, 0 seconds = [0][2][0][0].

Confirm with [↵] key.

0 2 0 0


Count-Down
02:00
start: ↵

The display shows:

Start count down with [↵] key.

After countdown has finished the photometer reverts to the previous display automatically.

2.3 Operation mode



Switch the photometer on by pressing the ON/OFF key.

Autotest ...

The photometer performs an electronic self-test.

2.3.1 Automatic switch off

The instrument switches off automatically after 20 minutes. This is indicated 30 seconds before by a beeper. Press any key to avoid the instrument switching off. As long as the instrument is working (for example countdown or printing) the automatic switch off is inactive.

2.3.2 Selecting a method

>> **30 Alkalinity-m**
35 Alkalinity-p
40 Aluminium

The display shows a selection:

There are two possibilities to select the required method:



a) enter method-number directly

e.g.: [8] [0] to select Bromine



b) press arrow key [▼] or [▲] to select the required method from the displayed list.



Confirm with [↵] key.

2.3.2.1 Method-Information (F1)

Use F1 key to switch between the compact and the detailed list for method selection.

100 Chlorine
0.02-6 mg/l Cl₂
Tablet
24 mm
DPD No 1
DPD No 3

Example:

Line 1: Method number, Method name

Line 2: Range

Line 3: Kind of reagent


Line 4: Vial

Line 5-7: Used reagent

tube = reagent vial contained in tube test

2.3.2.2 Chemical Species Information


Pressing the F2 key the display shows a list with available chemical species and corresponding ranges. Changing chemical species see chapter 2.3.7 page 198.



2.3.3 Differentiation

Chlorine >> diff free total
--

Differentiation is possible in some methods (e.g. Chlorine). The photometer then requires the type of determination.



Press arrow key [▼] or [▲] to select the required determination.

Confirm with [↵] key.

2.3.4 Performing Zero

prepare Zero press ZERO
--

The display shows:

Zero

Prepare a clean vial as described in "Method" and place the vial in the sample chamber making sure that the X marks are aligned.

Press **ZERO** key.

Zero accepted prepare Test press TEST
--

The display shows:

2.3.5 Performing Tests

When zero calibration is complete, remove the vial from the sample chamber and perform the tests as described under "Method".


When the results have been displayed:

- at some methods you can change between different chemical species
- you can store and/or print out the results
- perform further analysis with the same zero
- select a new method

2.3.6 Ensuring reaction periods (countdown)

For the compliance with reaction periods there is incorporated a time delay, the countdown.

There are two kinds of countdowns:



- Press **↵** key.
Prepare water sample, start countdown with **↵** key and proceed as described in the mode description.
The vial must not be placed in the sample chamber.



- Press **TEST** key.
Prepare the water sample as described in the method description and place the vial in the sample chamber. The display shows the countdown by pressing the **TEST** key and the countdown is started automatically. After the reaction period is finished the reading starts automatically.



Notes:

1. It is possible to finish the working countdown by pressing the **↵** key. Reading starts immediately. In this case the operator is responsible for ensuring the necessary reaction period by himself.

Non-compliance with reaction periods lead to incorrect test results.

2. The time remaining is displayed continuously.
The beeper indicates the last 10 seconds.

2.3.7 Changing chemical species

For some methods there is a possibility to change the chemical species of the test result. If the test result is displayed press arrow key [▲] or [▼].

Example:

320 Phosphate LR T	-----[▼]----->	320 Phosphate LR T	<----- [▼] -----	320 Phosphate LR T
0.05-4 mg/l PO ₄		0.02-1.3 mg/l P		0.04-3 mg/l P ₂ O ₅
	<----- [▲] -----		----- [▲] ----->	
1.00 mg/l PO ₄		0.33 mg/l P		0.75 mg/l P ₂ O ₅


If the special species of a test result is changed the displayed range is adjusted automatically. For an already stored result it is not possible to change the chemical species. The last displayed chemical species is kept by the instrument and will be displayed if this method is used the next time. If there is the possibility to change the chemical species for a method it is described in the manual. The arrows with the possible chemical species are printed below the notes of the method:

- ▲ PO₄
- P
- ▼ P₂O₅


2.3.8 Storing results



Press **STORE** during the test result is displayed.



The display shows:



- We advise you to enter a numeric code (up to 6 places). (A Code-No. can contain references to the operator or the sample-taking place.)



After entering confirm with [←] key.

- If a code number is not necessary confirm by pressing [←] directly. (The assignment for the Code-No. is then 0 automatically.)

The entire data set is stored with date, time, Code-No., method and test result.



The display shows:

The test result is then shown again.

Storage: 900
free records left

Note:

The display shows the number of free data sets.

Storage: only 29
free records left

If there are less than 30 data sets free the display shows:

Clear the memory as soon as possible (see "Deleting stored results"). If memory capacity is used up it would be impossible to save additional test results.

2.3.9 Printing results

If a printer is installed and switched on, it is possible to print out the test results (without saving it before).

F3

Press **F3** key.

The entire data set is printed with date, time, Code-No., method and test result. Printing example:

```
100 Chlorine T  
0.02-6 mg/l Cl2  
Profi-Mode: no  
2006-07-01 14:53:09  
Test No.: 1  
Code-Nr.: 007  
4.80 mg/l Cl2
```

The test No. is an internal number that is set automatically if a test result is stored. It appears only at the print out.

2.3.10 Perform additional measurements

Test

To perform additional tests using the same method:

Zero accepted
prepare Test
press TEST

- Press **TEST** key

The display shows:

Test

Confirm with **TEST** key

or

Zero

- Press **ZERO** key to perform a new zero calibration.

prepare Zero
press ZERO

The display shows:

2.3.11 Selecting a new method



Press [ESC] key to return to method selection.



Or enter the required method number directly, e.g. [1] [6] [0] for Cyanuric acid.



Confirm with [↵] key.

2.3.12 Measure absorbance

Range: -2600 mAbs to +2600 mAbs

Method-No.	Title
900	mAbs 430 nm
910	mAbs 530 nm
920	mAbs 560 nm
930	mAbs 580 nm
940	mAbs 610 nm
950	mAbs 660 nm

Select the desired wavelength from the method list or by entering the corresponding method-number directly.

900 mAbs 430 nm
-2600 mAbs - + 2600 mAbs
prepare Zero
press ZERO

The display shows e.g.:

Perform zeroing always with a filled (e.g. deionised water) vial.

Zero accepted
prepare Test
press TEST

The display shows:

Perform measurement of the sample.

500 mAbs

The display shows e.g.:

TIP: To ensure reaction times the User-Countdown may be helpful (chapter 2.2.3, page 194).

2.4 Photometer settings <MODE-Menu>

Table of Mode-Functions

MODE-Function	No.	Description	Page
Calibration	40	Performance of fluoride	215
Clear calibration	46	Deleting user calibration	220
Clock	12	Setting date and time	203
Countdown	13	Switching the countdown on/off to ensure reaction times	204
Delete data	34	Deleting all stored results	215
Key beep	11	Switching the acoustic signal on/off to indicate key-pressing	203
Langelier	70	Calculation of Langelier saturation Index (Water Balance)	232
Language	10	Selecting language	202
LCD contrast	80	Setting the display contrast	234
Method list	60	User method list, adaption	222
M list all on	61	User method list, switching on all methods	223
M list all off	62	User method list, switching off all methods	223
Print	20	Printing all stored results	206
Print, code-Nr.	22	Print only results of a selected Code-No. range	208
Print, date	21	Print only results of a selected time period	207
Print, method	23	Print only results of one selected method	209
Printing parameters	29	Setting of printing options	210
Profi-Mode	50	Switching the detailed operator instructions on/off	221
Signal beep	14	Switching the acoustic signal on/off to indicate end of reading	205
Storage	30	Displaying all stored results	211
Stor., code	32	Displaying only results of a selected Code-No. range	213
Stor., date	31	Displaying only results of a selected time period	212
Stor., method	33	Displaying only results of one selected method	214
System-info	91	Information about the instrument e.g. current software-version	234
Temperature	71	Selection of °C or °F for Langelier Mode 70	233
User calibration	45	Storage user calibration	219

User concentration	64	Entering of the data that are necessary to run a user concentration method	224
User polynoms	65	Entering of the data that are necessary to run a user polynomial	226
User methods clear	66	Delete all data of a user polynomial or of a concentration method	229
User methods print	67	Print out all data that are stored with mode 64 (concentration) or mode 65 (polynomial)	230
User methods init	69	Initialise the user-method system (polynomial and concentration)	231

The selected settings are kept by the photometer also after it was switched off. To change photometer settings a new setting is required.

2.4.1 blank because of technical requirements

2.4.2 Instrument basic settings 1


Selecting a language



Press [MODE] [1] [0] keys.



Confirm with [↵] key.




The display shows:

Press arrow key [▼] or [▲] to select the required language from the displayed list.



Confirm with [↵] key.

Key-beep



Press [MODE] [1] [1] keys.



Confirm with [↵] key.

<Key-Beep>
ON: 1 OFF: 0

The display shows:



- Press [0] key to switch the key beep off.



- Press [1] key to switch the key beep on.



Confirm with [↵] key.

Note:

In the case of methods with reaction periods, an acoustic signal still sounds during the last 10 seconds of the countdown even if the key-beep is switched off.

Setting date and time



Press [MODE] [1] [2] keys.



Confirm with [↵] key.

<clock>
yy-mm-dd hh:mm
--:-- --:--

The display shows:

The entering comprises two digits each.

yy-mm-dd hh:mm
06-05-14 --:--

Enter year, month and day,
e.g.: 14. Mai 2006 = [0][6][0][5][1][4]

yy-mm-dd hh:mm
06-05-14 15:07

Enter hours and minutes
e.g.: 3.07 p.m. = [1][5][0][7]



Confirm with [↵] key.

Note:

While conforming date and time with [↵] key the seconds are adjusted to zero automatically.

Countdown (Ensuring reaction periods)


Some methods require a reaction period. This reaction period is incorporated in the method as standard by the countdown function.

It is possible to switch the countdown off for all methods:

   Press [MODE] [1] [3] keys.



Confirm with [↵] key.



The display shows:



- Press [0] key to switch the countdown off.



- Press [1] key to switch the countdown on.



Confirm with [↵] key.

Notes:

1. It is possible to finish the working countdown by pressing the [↵] key (application e.g. serial analysis).
The "user-countdown" is also available if the countdown is switched off.
2. If the countdown function is switched off, the operator is responsible for ensuring the necessary reaction period by himself.

Non-compliance with reaction periods lead to incorrect test results.

Signal-beep

Performing a zero or a measurement takes 8 seconds. The photometer indicates the end of zeroing or measuring by a short beep.



Press [MODE] [1] [4] keys.



Confirm with [↵] key.



The display shows:



- Press [0] key to switch the signal-beep off.



- Press [1] key to switch the signal-beep on.



Confirm with [↵] key.

Note:

In the case of methods with reaction periods, an acoustic signal still sounds during the last 10 seconds of the countdown even if the key-beep is switched off.

2.4.3 Printing of stored results


Printing all results



Press [MODE] [2] [0] keys.



Confirm with [↵] key.



The display shows:

Press [↵] key for printing out all stored test results.



The display shows e.g.:

After printing the photometer goes back to <Mode-Menu> automatically.

Note:

It is possible to cancel the entering by [ESC].
All stored data are printed out.

Printing results of a selected time period



Press [MODE] [2] [1] keys.



Confirm with [↵] key.

<Print>
sorted: date
from yy-mm-dd
_ _ - _ - _

The display shows:

Enter year, month and day for the first day of the required period, e.g.: 14 Mai 2006 = [0][6][0][5][1][4]



Confirm with [↵] key.

to yy-mm-dd
_ _ - _ - _

The display shows:

Enter year, month and day for the last day of the required period, e.g.: 19 Mai 2006 = [0][6][0][5][1][9]



Confirm with [↵] key.

from 2006-05-14
to 2006-05-19
Start: ↵
cancel: ESC

The display shows:

Press [↵] key and all stored results in the selected date range are printed.

After printing the photometer goes back to mode menu automatically.

Note:

It is possible to cancel the entering by [ESC].

If you want to print only results of one day enter the same date twice to characterise the period.

Printing results of a selected Code-No. range



Press [MODE] [2] [2] keys.



Confirm with [↵] key.

<Print>
sorted: Code-No.
from _____

The display shows:

Enter numeric code number (up to 6 places) for the first required Code-No., e.g.: [1].



Confirm with [↵] key.

to _____

The display shows:

Enter numeric code number (up to 6 places) for the last required Code-No., e.g.: [1] [0].



Confirm with [↵] key.

from 000001
to 000010
Start: ↵
cancel: ESC

The display shows:

Press [↵] key and all stored results in the selected Code-Number range are printed.

After printing the photometer goes back to mode menu automatically.

Note:

It is possible to cancel the entering by [ESC].

If you want to print only results of one Code-Number enter the same Code-Number twice.

If you want to print all results without Code-No. (Code-Nr. is 0) enter Zero [0] twice.

Printing results of one selected method



Press [MODE] [2] [3] keys.



Confirm with [↵] key.

```
<Print>
>>20 Acid demand
  30 Alkalinity-tot
  40 Aluminium T
```

The display shows:

Select the required method from the displayed list or enter the method-number directly.



Confirm with [↵] key.

In case of differentiated methods select the required kind of determination and confirm with [↵] key.

```
<Print>
method
30 Alkalinity-tot
Start:  ↵
cancel: ESC
```

The display shows:

Press [↵] key and all stored results of the selected method are printed.

After printing the photometer goes back to mode menu automatically.

Printing Parameter



Press [MODE] [2] [9] keys.



Confirm with [↵] key.

```
<printing parameter>
1: Flow control
2: Baud rate

cancel:          ESC
```

The display shows:



Press [1] key to select "Flow control".

```
<Flow Control>
is: Hardware
select:  ↑ ↓
save:    ↵
cancel:  ESC
```

The display shows:



Press arrow key [▼] or [▲] to select the required Protocol.
(Xon/Xoff, Hardware, no control)



Confirm with [↵] key.



Finish with [ESC] key.
Flow Control will be set to the selection displayed at "is".



Press [2] key to select "Baud rate".

```
<Baud rate>
is: 19200
select:  ↑ ↓
save:    ↵
cancel:  ESC
```

The display shows:



Press arrow key [▼] or [▲] to select the required baud rate.
(600, 1200, 2400, 4800, 9600, 14400, 19200)



Confirm with [↵] key.



End with [ESC] key.

Back to Mode-Menu with [ESC] key.

Back to method selection with [ESC] key.

Note:

Select "Hardware" as Protocol and "19200" as baud rate if you use the printer **DP 1012**.
Select "Hardware" as Protocol and "9600" as baud rate if you use the printer **DPN 2335**.
For setting of the printer see chapter 2.5.1 Connection to a printer.

2.4.4 Recall / delete stored results

Recall all stored results



Press [MODE] [3] [0] keys.



Confirm with [↵] key.

```
<Storage>
display all data
Start:  ↵ cancel:  ESC
print:  F3
print all: F2
```

The display shows:

The stored data sets are displayed in chronological order, started with the latest stored test result. Press [↵] key and all stored results are displayed.

- Press [F3] key to print the displayed result.
- Press [F2] key to print all results.
- End with [ESC].
- Press arrow key [▼] to display the following test result.
- Press arrow key [▲] to display the previous test result.



```
no data
```

If there are no test results in memory the display shows:

Recall results of a selected time period



Press [MODE] [3] [1] keys.



Confirm with [↵] key.

```
<Storage>
sorted: date
from yy-mm-dd
_ _ - _ - _
```

The display shows:

Enter year, month and day for the for the first day of the required period, e.g.: 14 Mai 2006 = [0][6][0][5][1][4]



Confirm with [↵] key.

```
to yy-mm-dd
_ _ - _ - _
```

The display shows:

Enter year, month and day for the last day of the required period, e.g.: 19 Mai 2006 = [0][6][0][5][1][9]



Confirm with [↵] key.

```
from 2006-05-14
to 2006-05-19
Start: ↵ cancel: ESC
print: F3
print all: F2
```

The display shows:

- Press [↵] key and all stored results in the selected date range are displayed.
- Press [F3] key to print the displayed result.
- Press [F2] key to print all selected results.
- End with [ESC].

Note:

It is possible to cancel the entering by [ESC].

If you want to recall only results of one day enter the same date twice to characterise the time period.

Recall results of a selected Code-No. range



Press [MODE] [3] [2] keys.



Confirm with [↵] key.

```
<Storage>
sorted: Code-No.
from _ _ _ _ _
```

The display shows:

Enter numeric code number (up to 6 places) for the first required Code-No., e.g.: [1].



Confirm with [↵] key.

```
to _ _ _ _ _
```

The display shows:

Enter numeric code number (up to 6 places) for the last required Code-No., e.g.: [1] [0].



Confirm with [↵] key.

```
from 000001
to 000010
Start: ↵ cancel: ESC
print: F3
print all: F2
```

The display shows:

- Press [↵] key and all stored results in the selected Code-No. range are displayed.
- Press [F3] key to print the displayed result.
- Press [F2] key to print all selected results.
- End with [ESC].

Note:

It is possible to cancel the entering by [ESC].

If you want to recall only results of one Code-Number enter the same Code-Number twice.

If you want to recall all results without Code-No. (Code-Nr. is 0) enter Zero [0] twice.

Recall results of one selected method



Press [MODE] [3] [3] keys.



Confirm with [↵] key.

```
<Storage>
>>20 Acid demand
  30 Alkalinity-tot
  40 Aluminium T
```

The display shows:

Select the required method from the displayed list or enter the method-number directly.



Confirm with [↵] key.

In case of differentiated methods select the required kind of determination and confirm with [↵] key.

```
<Storage>
method
30 Alkalinity-tot
Start: ↵ cancel: ESC
print: F3
print all: F2
```

The display shows:

- Press [↵] key and all stored results of the selected method are displayed.
- Press [F3] key to print the displayed result.
- Press [F2] key to print all selected results.
- End with [ESC].

Delete stored results



Press [MODE] [3] [4] keys.



Confirm with [↵] key.

```
<Delete data>
Delete all data?
YES : 1 NO : 0
```

The display shows:



- Press [0] key to retain the data sets in memory.



- After pressing key [1] the following acknowledgment is displayed:

```
<Delete data>
Delete data ↓
Do not delete: ESC
```

Press [↵] key to delete.



ATTENTION:
All stored test results are deleted

or cancel without deleting data by pressing [ESC] key.

Note:

All stored test results are deleted.

2.4.5 Calibration

Calibration Fluoride



Regard notes!

Press [MODE] [4] [0] keys.



Confirm with [↵] key.

```
<Calibration>
170 Fluoride
Zero: deionised water
press ZERO
```

The display shows:

1. Fill a clean vial (24 mm Ø) with exact **10 ml of deionised water**, close the vial with the cap tightly.
2. Place the vial in the sample chamber making sure that the marks X are aligned.

Zero accepted
T1: 0 mg/l F
press TEST

T1 accepted
T2: 1 mg/l F
press TEST

Calibration
accepted



Esc

1

7

0



Error, absorbance
T2>T1

3. Press **ZERO** key.
4. Remove the vial from the sample chamber.
5. Add **exact 2 ml SPADNS reagent solution** to the water sample. **Caution: Vial is filled up to the top!**
6. Close the vial with the cap tightly and swirl the vial gently several times to mix the contents.
7. Place the vial in the sample chamber making sure that the Σ marks are aligned.
8. Press **TEST** key.
9. Remove the vial from the sample chamber, empty the vial, rinse vial and cap several times and fill the vial with exact **10 ml Fluoride standard** (Concentration 1 mg/l F).
10. Add **exact 2 ml SPADNS reagent solution** to the Fluoride standard.
Caution: Vial is filled up to the top!

11. Place the vial in the sample chamber making sure that the Σ marks are aligned.
12. Press **TEST** key.

The display shows:

Confirm with [\downarrow] key.

Back to method selection with **ESC** key.

Select method Fluoride with keys [1][7][0] und [\downarrow].

Note:

The same batch of SPADNS reagent solution must be used for adjustment and test. The adjustment process needs to be performed for each new batch of SPANDS reagent solution (see Standard methods 20th, 1998, APHA, AWWA, WEF 4500 F D., S. 4-82).

As the test result is highly dependent on exact sample and reagent volumes, the sample and reagent volumes should always be metered by using a 10 ml resp. 2 ml volumetric pipette (class A).

If there appears an error message please repeat adjustment.

User-Calibration

If a test method is user calibrated the method name is displayed inverse.

Procedure:

- Prepare a standard of known concentration and use this standard instead of the sample according to the test procedure.
- It is recommend to use well known standards which are formulated according to DIN EN, ASTM or other international norms or to use certified standards which are commercially available.
- After measuring this standard solution it is possible to change the displayed results to the required value.
- If a method use a mathematic equation for the calculation of the result, it is only possible to calibrate the basic tests since all the other tests use the same polynom.
- The same applies for some test procedures which use a polynom of another test procedure.

Return to factory calibration:

If the user calibration is deleted the factory calibration is automatically activated.

Remarks:

The method "Fluoride" cannot be calibrated with mode 45 since the test requires a calibration related to the batch of the liquid reagent (SPADNS) (mode 40, chapter "calibration (fluoride)").

Table

No.	Method	Recommended range for user user-calibration
20	Acid demand	1-3 mmol/l
35	Alkalinity-p	100-300 mg/l CaCO ₃
30	Alkalinity-total	50-150 mg/l CaCO ₃
40	Aluminium T	0.1-0.2 mg/l Al
50	Aluminium PP	0.1-0.2 mg/l Al
60	Ammonium T	0.3-0.5 mg/l N
62	Ammonium PP	0.3-0.5 mg/l N
65	Ammonium LR TT	1 mg/l N
66	Ammonium HR TT	20 mg/l N
85	Boron	1 mg/l B
80	Bromine	Calibration with basic test 100 Chlorine free
90	Chloride	10-20 mg/l Cl
100	Chlorine T	0.5-1.5 mg/l Cl
101	Chlorine L	Calibration with basic test 100 Chlorine free
110	Chlorine PP	0.5-1 mg/l Cl ₂
105	Chlorine (KI) HR	70-150 mg/l Cl
120	Chlorine dioxide	Calibration with basic test 100 Chlorine free
130	COD LR	100 mg/l O ₂
131	COD MR	500 mg/l O ₂
132	COD HR	5 g/l O ₂ = 5000mg/lO ₂
150	Copper T	0.5-1.5 mg/l Cu

No.	Method	Recommended range for user user-calibration
153	Copper PP	0.5-1.5 mg/l Cu
157	Cyanide	0.1-0.3 mg/l CN
160	Cyanuric acid	30-60 mg/l Cys
165	DEHA T	200-400 µg/l DEHA
167	DEHA PP	200 µg/l DEHA
170	Fluoride	Calibration with 0 und 1 mg/l F through Mode 40
190	Hardness, Calcium	100-200 mg/l CaCO ₃
200	Hardness, total T	15-25 mg/l CaCO ₃
201	Hardness, total HR T	Calibration with basic test 200 Chlorine free
205	Hydrazine P	0.2-0.4 mg/l N ₂ H ₄
206	Hydrazine L	0.2-0.4 mg/l N ₂ H ₄
207	Hydrazine C	0.2-0.4 mg/l N ₂ H ₄
210	Hydrogen peroxide	Calibration with basic test 100 Chlorine free
215	Iodine	Calibration with basic test 100 Chlorine free
220	Iron T	0.3-0.7 mg/l Fe
222	Iron PP	0.1-4 mg/l Fe
223	Iron (TPTZ) PP	0.3-0.7 mg/l Fe
240	Manganese T	1-2 mg/l Mn
242	Manganese PP	0.1-0.4 mg/l Mn
243	Manganese HR PP	4-6 mg/l Mn
250	Molybdate T	5-15 mg/l Mo
252	Molybdate HR PP	10-30 mg/l Mo
265	Nitrate TT	10 mg/l N
270	Nitrite T	0.2-0.3 mg/l N
272	Nitrite LR PP	0.1-0.2 mg/l N
280	Nitrogen, total LR	10 mg/l N
281	Nitrogen, total HR	50-100 mg/l N
300	Ozone (DPD)	Calibration with basic test 100 Chlorine free
290	Oxygen, active	Calibration with basic test 100 Chlorine free
292	Oxygen, dissolved	possible against meter for dissolved oxygen
280	Nitrogen, total LR	10 mg/l N
281	Nitrogen, total HR	50-100 mg/l N
329	pH-Value LR	6.0-6.6
330	pH-Value T	7.6-8.0
331	pH-Value L	7.6-8.0
332	pH-Value HR	8.6-9.0
70	PHMB	15-30 mg/l
320	Phosphate LR T	1-3 mg/l PO ₄
321	Phosphate HR T	30-50 mg/l PO ₄
323	Phosphate, ortho PP	0.1-2 mg/l PO ₄
324	Phosphate, ortho TT	3 mg/l PO ₄
327	Phosphate 1, ortho C	20-30 mg/l PO ₄
328	Phosphate 2, ortho C	1-3 mg/l PO ₄
325	Phosphate, total TT	0.3-6 mg/l P
326	Phosphate, hydr. TT	0.3-0.6 mg/L P
340	Potassium	3 mg/l K
350	Silica	0.5-1.5 mg/l SiO ₂

No.	Method	Recommended range for user user-calibration
351	Silica LR PP	1 mg/l SiO ₂
352	Silica HR PP	50 mg/l SiO ₂
212	Sodium hypochlorite	8 %
360	Sulfate PP	50 mg/l SO ₄
355	Sulfate T	50 mg/l SO ₄
365	Sulfide	0.2-0.4 mg/l S
370	Sulfite	3-4 mg/l SO ₃
390	Urea	1-2 mg/l CH ₄ N ₂ O
400	Zinc	0.2-0.4 mg/L Zn

Store user-calibration

100 Chlorine T
0.02-6 mg/l Cl₂
0.90 mg/l free Cl₂

Perform the required method as described in the manual using a standard of known concentration instead of the water sample.



If the test result is displayed press [MODE] [4] [5] keys and confirm with [↵] key.



The display shows:

<user calibration>
100 Chlorine T
0.02-6 mg/l Cl₂
0.90 mg/l free Cl₂
up: ↑, down: ↓
save: ↵

Pressing the arrow key [▲] once increases the displayed result.

Pressing the arrow key [▼] once decreases the displayed result.

Press keys till the displayed result corresponds to the value of the standard.



Confirm with [↵] key to store the new calibration factor.

Cancel user calibration by pressing [ESC] key.

Jus Factor
saved

The display shows:

100 Chlorine T
0.02-6 mg/l Cl₂
1.00 mg/l free Cl₂

Now the method name is displayed inverse and the test result is calculated with the new calibration factor.

Delete user-calibration

This chapter only applies for methods which can be user-calibrated.

100 Chlorine T
0.02-6 mg/l Cl2

Select the required method.

prepare ZERO
press ZERO

Instead of zeroing the instrument press [MODE] [4] [6] keys and confirm with [↵] key.

Mode **4** **6**



<user calibration>
100 Chlorine T
0.02-6 mg/l Cl2
clear user
calibration?
YES: 1, NO: 0

The display shows:

1

Press [1] key to delete user-calibration.

0

Press [0] key to keep the valid user-calibration.

The instrument goes back to Zero-query automatically.

2.4.6 Lab function

Reduced operator guidance => "Profi-Mode"

This function may be used for routine analyses with many samples of one method.
The following information is always stored in the methods:

- Method
- Range
- Date and time
- Differentiation of results
- Detailed operator instruction
- Compliance with reaction periods

If the Profi-Mode is active, the photometer provides only a minimum of operator instructions.
The criteria specified above d, e, f are not longer included.

   Press [MODE] [5] [0] keys in succession.



Confirm with [↵] key.

<Profi-Mode>
ON : 1 OFF : 0

The display shows:



- Press [0] key to switch the Profi-Mode off.



- Press [1] key to switch the Profi-Mode on.

switched off

The display shows:

or

switched on



Confirm with [↵] key.

Note:

Storage of test results is possible. In case of stored test results the display shows "Profi-Mode" additionally.

The selected settings are kept by the photometer also after it was switched off. To change photometer setting a new setting is required.

2.4.7 User operations

User-method list

After switching on the instrument a scroll list of all available methods is automatically shown in the display. To shorten this list according to the requirements of the user it is possible to create a user defined scroll list.

The program structure requires that this list must have at least one active (switched on) method. For this reason it is necessary to activate first all required methods and than to switch off the automatic activated one if this one is not required.

User-method list, adaptation



Press [MODE] [6] [0] keys.



Confirm with [↵] key.

```
<Method list>
selected: •
toggle: F2
save: ↵
cancel: ESC
```

The display shows:



Start with [↵] key.

```
<Method list>
>> 30•Alkalinity-tot
    40•Aluminium
    50•Ammonium
....
```

The complete method list is displayed.

Methods with a point [•] behind the method number will be displayed in the method selection list. Methods without a point will not be displayed in the method selection list.

```
>> 30•Alkalinity-tot
```

Press key [▲] or [▼] to select the required method from the displayed list.

```
[F2]
```

Switch with [F2] key between "active" [•] and "inactive" [].

```
>> 30 Alkalinity-tot
```

Select next method, activate or inactivate it and so on.

```
[F2]
```

```
>> 30•Alkalinity-tot
```

Confirm with [↵] key.



Cancel without storing by pressing [ESC] key.

Recommendation:

If only a few methods are required it is recommendable to perform Mode 62 first, followed by Mode 60.

All user-Polynomials (1-25) and -Concentrations (1-10) are displayed in the method list, although they are not programmed by the user. Non-programmed user-methods can't be activated!


User-method list, switch all methods on

This mode function activates all methods. After switching on the instrument a scroll list of all available methods is automatically shown in the display.



Press [MODE] [6] [1] keys.



Confirm with [] key.

**<Mlist all on>
switch on all
methods
YES: 1, NO: 0**

The display shows:

1

- Press [1] key to display all methods in the method selection list.

0

- Press [0] key to keep the valid method selection list.

The instrument goes back to mode-menu automatically.


User-method list, switch all methods off

The program structure requires that the method list must have at least one active (switched on) method. For this reason the instrument activates one method automatically.



Press [MODE] [6] [2] keys.



Confirm with [] key.

**<Mlist all off>
switch off all
methods
YES: 1, NO: 0**

The display shows:

1

- Press [1] key to display only one method in the method selection list.

0

- Press [0] key to keep the valid method selection list.

The instrument goes back to mode-menu automatically.

User-Concentration-Methods

It is possible to enter and store up to 10 User-Concentration-Methods.

Therefore you need 2 to 14 standards of known concentration and one blank (deionised water or reagent blank value). The Standards should be measured with increasing concentrations and from the brightest to the darkest colouration.

The measuring range for „Underrange“ and „Ovrange“ is defined with -2600 mAbs* and +2600 mAbs*. After selection of a method the concentration of the lowest and highest used standard is displayed as measuring range. The operation range should be within these range to achieve best results.

*1000 mAbs = 1 Abs = 1 E

Entering a User Concentration:



Press [MODE] [6] [4] keys.



Confirm with [↵] key.

```
< User concentr.>
choose no.: ____
(850-859)
```

The display shows:

[8] [5] [0]

Enter a method-number in the range from 850 to 859, e.g.: [8] [5] [0]



Confirm with [↵] key.

```
Overwrite conc. meth.?
YES: 1, NO: 0
```

Note:

if the entered number has already been used to save a concentration the display shows the query:

- Press [0] or [ESC] key to go back to method-No. query.
- Press [1] key to start entry-mode.

```
wavelength:
1: 530 nm    4: 430 nm
2: 560 nm    5: 580 nm
3: 610 nm    6: 660 nm
```

Enter the required wavelength, e.g.: [2] for 560 nm.

[2]

```
choose unit:
>>
    mg/l
    g/l
    mmol/l
    mAbs
    µg/l
    E
    A
    %
```

Press [▲] or [▼] keys to select the required unit.



choose resolution
 1: 1
 2: 0.1
 3: 0.01
 4: 0.001

Confirm with [↵] key.

Press the appropriate numerical key to select the required resolution.

Note:
 Please enter the required resolution according to the instrument presetting:

range	max. resolutions
0.000 ...9.999	0.001
10.00 ...99.99	0.01
100.0... 999.9	0.1
1000 ...9999	1

< User concentr.>
 prepare Zero
 press ZERO



Measurement procedure with standards of known concentration:

The display shows:

Prepare Zero and press [Zero] key.

Note:
 Use deionised water or reagent blank value.

The display shows:

< User concentr.>
 Zero accepted
 S1: +_____
 ↵ | ESC | F1

0 . 0 5

Enter the concentration of the first standard; e.g.: 0.05

- One step back with [ESC].
- Press [F1] key to reset numerical input.



Confirm with [↵] key.

< User concentr.>
 S1: 0.05 mg/l
 prepare
 press TEST



The display shows:

Prepare the first standard and press [Test] key.

The display shows the input value and the measured extinction value. Confirm with [↵] key.

S1: 0.05 mg/l
 mAbs: 12 ↵

Enter the concentration of the second standard; e.g.: 0.1

- One step back with [ESC].
- Press [F1] key to reset numerical input.

S1 accepted
 S2: +_____
 ↵ | ESC | F1

0 . 1

Confirm with [↵] key.



S2: 0.10 mg/l
prepare
press TEST

Prepare the second standard and press [Test] key.

S2: 0.10 mg/l
mAbs: 150 ↵

The display shows the input value and the measured extinction value. Confirm with [↵] key.

S2 accepted
S3: + _____
↵ | ESC | F1 | Store

Note:

- Perform as described above to measure further standards.
- The minimum of measured standards is 2.
- The maximum of measured standards is 14 (S1 to S14).

Store

If all required standards or the maximum value of 14 standards are measured press [Store] key.

stored!

The display shows:

The instrument goes back to the mode menu automatically.

Now the concentration is stored in the instrument and can be recalled by entering its method number or selecting it from the displayed method list.

TIP:

Save all your concentration data in a written form because in case of power outage (e.g. changing the battery) all concentration data will be lost and must be entered again. You might want to use Mode 67 to transfer all concentration data to a PC.

User-Polynomials

It is possible to enter and store up to 25 User-Polynomials.

The program allows the user to apply a Polynomial up to the 5th degree:

$$y = A + Bx + Cx^2 + Dx^3 + Ex^4 + Fx^5$$

If only a Polynomial of a lower degree is necessary the other coefficients are specified as zero (0), e.g.: for the 2nd degree is D, E, F = 0.

The values of the coefficients A, B, C, D, E, F must be entered in an academic notation with maximal 6 decimal places, e.g.: 121,35673 = 1,213567E+02

Entering a User-Polynomial:

Mode 6 5

Press [MODE] [6] [5] keys.

↵

Confirm with [↵] key.

<User polynoms>
choose no.: ____
(800-824)

The display shows:

8 0 0

Enter a method-number in the range from 800 to 824, e.g.: [8] [0] [0]



Confirm with [↵] key.

Overwrite polynom?
YES: 1, NO: 0

Note:
if the entered number has already been used to save a polynomial the display shows the query:

- Press [0] or [ESC] key to go back to method-No. query.
- Press [1] key to start entry-mode.

wavelength:
 1: 530 nm 4: 430 nm
 2: 560 nm 5: 580 nm
 3: 610 nm 6: 660 nm

Enter the required wavelength, e.g.: [2] for 560 nm.



< User polynoms >
 $y = A+Bx+Cx^2+Dx^3+Ex^4+Fx^5$
 A: + _____

- Press [▲] or [▼] key to change between plus and minus sign
- Enter data of the coefficient A including decimal point, e.g.: 1.32



Confirm with [↵] key.

A: 1.32 ___ E+ ___

- Press [▲] or [▼] key to change between plus and minus sign
- Enter the exponent of the coefficient A, e.g.: 3



Confirm with [↵] key.

B: + _____

Successively the instrument queries the data for the other coefficients (B, C, D, E and F).

Note:
If zero [0] is entered for the value of the coefficient, the input of the exponent is omitted automatically.



Confirm every input with [↵] key.

measurement range
 Min mAbs: + _____
 Max mAbs: + _____

Enter measurement ranges from –2600 to +2600 mAbs.

- Press [▲] or [▼] key to change between plus and minus sign.
- Enter the values in Absorbance (mAbs) for the upper limit (Max) and the lower limit (Min).



Confirm every input with [↵] key.

choose unit:
 >>
 mg/l
 g/l
 mmol/l
 mAbs
 µg/l
 E
 A
 %

Press [▲] or [▼] keys to select the required unit.



Confirm with [↵] key.

choose resolution
 1: 1
 2: 0.1
 3: 0.01
 4: 0.001

Press the appropriate numerical key to select the required resolution.

Note:

Please enter the required resolution according to the instrument presetting:

range	max. resolutions
0.000 ...9.999	0.001
10.00 ...99.99	0.01
100.0... 999.9	0.1
1000 ...9999	1

stored!

The display shows:

The instrument goes back to the mode menu automatically.

Now the polynomial is stored in the instrument and can be recalled by entering its method number or selecting it from the displayed method list.

TIP:

Save all your polynomial data in a written form because in case of power outage (e.g. changing the battery) all polynomial data will be lost and must be entered again. You might want to use Mode 67 to transfer all polynomial data to a PC.

Delete User-Methods (Polynomial or Concentration)

In principle a valid user-method can be overwritten.

An existing user-method (Polynomial or Concentration) can be totally deleted as well and is removed out of the method selection list:



Press [MODE] [6] [6] keys.



Confirm with [↵] key.

```
<User m. clear>
choose no.: _____
(800-824), (850-859)
```

The display shows:

8 0 0

Enter the number of the User-Method you want to delete (in the range from 800 to 824 or 850 to 859), e.g.: 800



Confirm with [↵] key.

```
M800
delete?
YES: 1, NO: 0
```

There is displayed the query:

1

- Press [1] key to delete the selected User-Method.

0

- Press [0] key to keep the valid User-Method.

The instrument goes back to mode menu automatically.

Print Data of User-Methods (Polynomials & Concentration)

With these Mode function all data (e.g. wavelength, unit ...) of stored user-polynomial and concentration methods can be printed out or transferred with HyperTerminal to a PC.



Press [MODE] [6] [7] keys.



Confirm with [↵] key.

```
<User m. print>
Start: ↵
```

The display shows:



Press [↵] key to print out the data (e.g. wavelength, unit, ...) of all stored User-Methods.

```
M800
M803
...
```

The display shows e.g.:

After data transfer the photometer goes back to mode menu automatically.

Initialise User-Method-System (Polynomials & Concentration)

Power loss at the storage device will cause incoherent data. The user-method system must be initialised with this mode function to set it to a predefined state.

ATTENTION:

All stored user-methods (polynomial & concentration) are deleted with initialisation.



Press [MODE] [6] [9] keys.



Confirm with [↵] key.

<User m. init>
Start: ↵

The display shows:



Confirm with [↵] key.

Initialising?
YES: 1, NO: 0

There is displayed the query:



- Press [1] key to start initialisation.



- Press [0] key to to cancel without initialisation.

The instrument goes back to mode menu automatically.

2.4.8 Special functions

Langelier Saturation Index (Water Balance)


For calculation the following tests are required:

- pH-value
- Temperature
- Calcium hardness
- Total Alkalinity
- TDS (Total Dissolved Solids)

Run the test separately and note the results.

Calculate the Langelier Saturation Index as described:

Calculation of Langelier Saturation Index



With Mode 71 (see below) it is possible to select between degree Celsius or degree Fahrenheit.

Press [MODE] [7] [0] keys.

Confirm with [↵] key.

<Langelier>
temperature °C:
3°C <=T<=53°C
+ _ _ _ _

The display shows:



Enter the temperature value (T) in the range between 3 and 53°C and confirm with [↵] key. If °F was selected, enter the temperature value in the range between 37 und 128°F.

calcium hardness
50<=CH<=1000
+ _ _ _ _

The display shows:



Enter the value for Calcium hardness (CA) in the range between 50 and 1000 mg/l CaCO₃ and confirm with [↵] key.

tot. alkalinity
5<=TA<=800
+ _ _ _ _

The display shows:



Enter the value for Total Alkalinity (TA) in the range between 5 and 800 mg/l CaCO₃ and confirm with [↵] key.

total dissol. solids
0<=TDS<=6000
+ _ _ _ _

The display shows:



Enter the value for TDS (Total Dissolved Solids) in the range between 0 und 6000 mg/l and confirm with [↵] key.

pH value
0<=pH<=12
+ _ _ _ _

The display shows:



Enter the pH-value in the range between 0 and 12 and confirm with [↵] key.

<Langelier>
Langelier
saturation index
0.00
Esc ↵

The display shows the Langelier Saturation Index.

Press [↵] key to start new calculation.

Return to mode menu by pressing [ESC] key.

Examples:

Operating error:

Values out of defined range:

CH<=1000 mg/l CaCO3!

The entered value is to high.

CH>=50 mg/l CaCO3!

The entered value is to low.



Confirm display message with [↵] key and enter a value in the defined range.

Selection of temperature unit

Entering the temperature value is possible in degree Celsius or degree Fahrenheit. Therefore the following preselection is (once) required.

Mode [7] [1]

Press [MODE] [7] [1] keys.



Confirm with [↵] key.

<temperature>
1: °C 2: °F

The display shows:

1

Press [1] key to select degree Celsius.

2

Press [2] key to select degree Fahrenheit.


The instrument goes back to mode menu automatically.


2.4.9 Instrument basic settings 2


Adjusting display contrast

   Press [MODE] [8] [0] keys.

 Confirm with [↵] key.

 The display shows:

 • Press arrow key [▲] to increase contrast of the LCD display.

 • Press arrow key [▼] to decrease contrast of the LCD display.


 Confirm with [↵] key.


2.4.10 Instrument special functions /service


Photometer-Information

   Press [MODE] [9] [1] keys.

 Confirm with [↵] key.

 This method informs you about the current software version, about the current detected mains power supply, about the number of performed tests and free memory capacity.

 Press arrow key [▼] to display the number of performed tests and free memory capacity.

 Finish with [ESC] key.

2.5 Data transfer

Switch the photometer and the personal computer or printer off. Connect the photometer (RS232 interface) and the serial interface of the personal computer or printer using a cable in line with the specified assignment (see technical data). The cable for connection to a personal computer is included in delivery contents.

2.5.1 Connection to a printer

Printer with a serial connection are suitable for connection with the photometer (see chapter 3.4 Technical data interface).

A suitable paper label printer is the printer DPN 2335.

Before using the printer **DPN 2335** with the Photometer you should change the following standard adjustments:
(Detailed information of changing the adjustment you will find in the printer manual).

Baud rate:	9600
Parity:	None
Data bits:	8

Note: The printer must be connected and switched on before printing.

Caution: Adjust printing parameter in Mode 29. See chapter 2.4.3 Printing Parameter.

2.5.2 Data transfer to a personal computer

Transferring test results from the photometer to a personal computer requires a transfer program, e.g. HyperTerminal.
Please find detailed information at our homepage on the download-area.

2.5.3 Internet-Updates

It is possible to update new software applications and additional languages via internet. Please find detailed information at our homepage on the download-area.

Remark:

To prevent loss of stored test results store or print out them before performing an Update.

2.6 blank because of technical requirements

Part 3

Enclosure

3.1 Unpacking

Carefully inspect all items to ensure that every part of the list below is present and no visible damage has occurred during shipment. If there is any damage or something is missing, please contact your local distributor immediately.

3.2 Delivery content

Standard content of PC MultiDirect:

-
- 1 Photometer in plastic case
- 1 Adapter for 16 mm Ø vials
- 1 Cap for adapter
- 2 Protective caps for connections
- 1 Rechargeable battery set (7 Ni-MH-cells; Type AA)
- 1 Lithium battery (CR 2032; 3V)
- 1 Mains adapter, 100 – 240 V, 50 – 60 Hz
- 1 Cable for connection to PC
- 3 Round vials with cap, height 48 mm, Ø 24 mm
- 3 Round vials with cap, height 90 mm, Ø 16 mm
- 1 Beaker cup, plastic, 100 ml
- 1 Cleaning brush
- 1 Stirring rod, plastic
- 1 Syringe, plastic, 2 ml
- 1 Syringe, plastic, 5 ml
- 1 Syringe, plastic, 10 ml
- 1 Instruction manual
- 1 Guarantee declaration

Reagent sets are not part of the standard scope of delivery. Please see the General Catalogue for details of available reagent sets.

3.3 Blank because of technical requirements

3.4 Technical data

Display	Graphic-Display (7-line, 21-characters)
Serial Interface	serial RS232 for printer- and PC-connection; 9-pin D-sub-mail connector, data format ASCII, 8-bit Data, no parity, 1 start-bit, 1 stop-bit, baud rate and protocol: adjustable Pin assignation: Pin 1 = free Pin 2 = Rx Data Pin 3 = Tx Data Pin 4 = free Pin 5 = GND Pin 6 = free Pin 7 = RTS Pin 8 = CTS Pin 9 = free
Light source	LEDs and photo sensor amplifier in protected cell compartment. Wavelength ranges: $\lambda_1 = 530 \text{ nm IF } \Delta \lambda = 5 \text{ nm}$ $\lambda_2 = 560 \text{ nm IF } \Delta \lambda = 5 \text{ nm}$ $\lambda_3 = 610 \text{ nm IF } \Delta \lambda = 6 \text{ nm}$ $\lambda_4 = 430 \text{ nm IF } \Delta \lambda = 5 \text{ nm}$ $\lambda_5 = 580 \text{ nm IF } \Delta \lambda = 5 \text{ nm}$ $\lambda_6 = 660 \text{ nm IF } \Delta \lambda = 5 \text{ nm}$ IF = Interference filter
Photometric accuracy*	0.100 Abs \pm 0.008 Abs 1.000 Abs \pm 0.020 Abs
Operation	Acid and solvent resistant touch-sensitive keyboard with integral beeper as acoustic indicator.
Power supply	7 Ni-MH cells (Type AA with 1100 mAh); external main adapter (Input: 100-240 V, 50-60 Hz; Output: 15V=530 mA) Lithium battery (CR 2032, 3V); for keeping data if there is no power supply from the rechargeable batteries or the main adapter
Auto off	20 minutes after last function, 30 seconds acoustical signal before switch off
Charging time	approx. 10 hours
Dimensions	approx. 265 x 195 x 70 mm (unit) approx. 440 x 370 x 140 mm (case)
Weight (unit batteries)	approx. 1000 g (with main adapter and rechargeable batteries)
Working condition (without condensation)	5 – 40°C at max. 30-90% relative humidity
Language options	English, German, French; Spanish, Italian further languages via Internet-Update
Storage capacity	ca. 1000 data sets

Subject to technical modification!


* measured with standard solutions

3.5 Abbreviations

Abbreviation	Definition
°C	degree Celsius (Centigrade)
°F	degree Fahrenheit °F = (°C x 1.8) + 32
°dH	degree German Hardness
°fH	degree French hardness
°eH	degree English Hardness
°aH	degree American Hardness
Abs	Absorption unit (\triangleq Extinction E) 1000 mAbs = 1 Abs \triangleq 1 A \triangleq 1 E
μ g/l	(= ppb) Microgram per liter
mg/l	(= ppm) Milligram per liter
g/l	(= ppt) gram per litre
KI	Potassium iodide
K _{S4.3}	Acid demand to pH 4.3 – this method is similar to the Total Alkalinity but converted into the unit “mmol/l”, as the German DIN 38409 demand.
TDS	Total Dissolved Solids
LR	Low Range
MR	Medium Range
HR	High Range
C	Reagents of Chemetrics®
L	Liquid reagent
P	Powder (-reagent)
PP	Powder Pack
T	Tablet
TT	Tube Test
DEHA	N,N-Diethylhydroxylamine
DPD	Diethyl-p-phenylendiamine
DTNB	Ellmans reagent
PAN	1-(2-Pyridylazo)-2-naphthol
PDMAB	Paradimethylaminobenzaldehyde
PPST	3-(2-Pyridyl)-5,6-bis(4-phenylsulfonic acid)1,2,4-triazine
TPTZ	2,4,6-Tri-(2-Pyridyl)-1,3,5-triazine

3.6 Troubleshooting

3.6.1 Operating messages in the display / error display

Display	Possible Causes	Elimination
Overrange	reading is exceeding the range water sample is too cloudy too much light on the photo cell	if possible dilute sample or use other measuring range filtrate water sample seal on the cap? Repeat measurement with seal on the cap of the vial.
Underrange	result is under the detection limit	indicate result with lower x mg/l x = low end of measuring range; if necessary use other analytical method
Storagesystem error use Mode 34	mains power fails or is not existing	insert or change Lithium battery. Delete data with Mode 34
capacity of rechargeable battery 	full capacity warning signal every 3 minutes warning signal every 12 seconds warning signal, the instrument switches itself off	capacity of the rechargeable battery is too low; charge the rechargeable battery; operate instrument with mains adapter
Jus Overrange E4	The user calibration is out of the accepted range	Please check the standard, reaction time and other possible faults. Repeat the user calibration.
Jus Underrange E4		
Overrange E1	The concentration of the standard is too high/too low, so that during user-calibration the limit of the range was exceeded	Perform the test with a standard of higher/lower concentration
Underrange E1		
E40 user calibration not possible	If the display shows Overrange/Underrange for a test result a user calibration is not possible	Perform the test with a standard of higher/lower concentration
Zero not accepted	Light absorption is too great or too low	Refer to chapter 2.3.4 Performing Zero (page 196) Clean sample chamber. Repeat zeroing.

Display	Possible Causes	Elimination
<p>???</p> <p>Example 1</p> <p>0,60 mg/l free Cl ??? comb Cl 0,59 mg/l total Cl</p> <p>Example 2</p> <p>Underrange ??? comb Cl 1,59 mg/l total Cl</p> <p>Example 3</p> <p>0,60 mg/l free Cl ??? comb Cl Overage</p>	<p>The calculation of a value (e.g. combined Chlorine) is not possible</p>	<p>Test procedure correct? If not – repeat test</p> <p>Example 1: The readings for free and total Chlorine are different, but considering the tolerances of each reading they are the same. For this reason the combined Chlorine is most likely zero.</p> <p>Example 2: The reading for free Chlorine is under the detection limit. The instrument is not able to calculate the combined Chlorine. In this case the combined Chlorine is most likely the same as the total Chlorine.</p> <p>Example 3: The reading for total Chlorine is exceeding the range. The instrument is not able to calculate the combined Chlorine. The test should be repeated with a diluted sample.</p>
<p>Error absorbance e.g.: T2>T1</p>	<p>calibration of Fluoride was not correct</p>	<p>Repeat calibration</p>
<p>Printer „timeout“</p>	<p>printer switched off; no connection</p>	<p>Connect printer Check connections Switch printer on</p>

3.6.2 General problems

Problem	Possible Causes	Elimination
Test result deviates from the expected.	Chemical species not as required.	Press arrow keys to select the required chemical species.
No differentiation: e.g. for the test Chlorine there is no selection between differentiated, free or total.	Profi-Mode is switched on.	Switch Profi-Mode off with Mode 50.
The pre-programmed countdown is not displayed.	Countdown is not activated and/or the Profi-Mode is activated.	Switch the countdown on with Mode 13 and/or switch the Profi-Mode off with Mode 50.
It seems that a method is not available.	Method is not activated in the user method list.	Activate the required method in the user method list with Mode 60.
Instrument can be operated with the mains adapter but not with the rechargeable batteries.	Rechargeable batteries are not charged or defect. Fuse (Type A, inert, 20 mm) may be defect.	Charge rechargeable batteries or change them. If the problem still exists change fuse.

3.7 Declaration of CE-Conformity


The manufacturer: **Tintometer GmbH**
Schleefstraße 8 a
44287 Dortmund
Germany

declares, that this product

Product name: **PC MultiDirect**

Conforms with EN 61 326 for specific defined electromagnetic environment.
Conforms with EN 61 326 (domestic).

Dortmund, 06. August 2003



Cay-Peter Voss, Managing Director



Tintometer GmbH

Lovibond® Water Testing
Schleefstraße 8-12
D-44287 Dortmund
Tel.: (+49) (0)2 31 / 9 45 10-0
Fax: (+49) (0)2 31 / 9 45 10-30
verkauf@tintometer.de
www.tintometer.de
Deutschland

Tintometer AG

Hauptstraße 2
CH-5212 Hausen AG
Tel.: (+41) (0)56 / 4 42 28 29
Fax: (+41) (0)56 / 4 42 41 21
tintometer@bluewin.ch
www.tintometer.ch

Schweiz

The Tintometer Limited

Lovibond House
Solar Way / Solstice Park
Amesbury, SP4 7SZ
Tel.: (+44) 19 80 664 800
Fax: (+44) 19 80 625 412
sales@tintometer.com
www.tintometer.com
UK



QM-System
Certificate No. 5394
ISO 9001 : 2000

Technical changes without notice
Printed in Germany 08/07

Lovibond® and Tintometer®
are registered trademarks of the
Tintometer group of companies