

# HPLC PTFE Manifolds



## Laboratory HPLC Disposal Containers - Laboratory HPLC Disposal Containers

### HPLC PTFE Manifolds

Enables up to six waste lines (per manifold) from HPLC processes to be fed into a single intake fitting. Safety container with manifold installed is OSHA compliant and meets NFPA Code 30. Use on round containers; not for use on oval style can.

#### Colors:

- Grey/White
- Silver/White

#### Sizes:

- 12.7cm x 17.8cm



# Laboratory HPLC Disposal Containers - Laboratory HPLC Disposal Containers

## HPLC PTFE Manifolds

### General

Brand	Justrite
Category	Laboratory HPLC Disposal Containers
Product Number	2817
Product Family	Safety Cans & Containers
Product Group	Laboratory HPLC Disposal Containers
Product Name	HPLC PTFE Manifolds
Shipping Location	Barendrecht 2992SB, NL
Update trigger	20210310-1301

## HPLC Poly Manifold with Stainless Steel Stacker Connector



JCN28179SI  
28179

Art. nr.	Art. Nr. Justrite US Ref (ERP)	Dimensions	Weight (kg)	Height (cm)
JCN28179SI	28179	11.4x13x18	0.1	11.4

## HPLC Poly Manifold With Polypropylene Quick Disconnect



JCN28178GY  
28178

Art. nr.	Art. Nr. Justrite US Ref (ERP)	Dimensions	Weight (kg)	Height (cm)
JCN28178GY	28178	11.4x13x18	0.1	11.4

## HPLC Poly Manifold With Stainless Steel Quick Disconnect



JCN28177SI  
28177

Art. nr.	Art. Nr. Justrite US Ref (ERP)	Dimensions	Weight (kg)	Height (cm)
JCN28177SI	28177	11.4x13x18	0.1	11.4