



TROUSSES D'ANALYSE [KIT TEST



Recharge Oto/Phénol

Ricambio Oto/Phenol



Référence	40061
Description	Recharge Oto/Phénol
Mesure	Ø 2,5 x 7,2 cm
Poids Oto/Phénol	0,018 kg / 0,021 Kg

Codice articolo	40061
Descrizione	Ricambio Oto/Phenol
Misura	Ø 2,5 x 7,2 cm
Peso Oto/Phenol	0,018 kg / 0,021 Kg



EMBALLAGE IMBALLAGGIO

Type Tipo	Blister / Blister
Unité Unità	1
Mesure Misura	9 x 2,3 x 20 cm
Volume Volume	0,00041 m ³
Poids Peso	0,06 kg



EMBALLAGE IMBALLAGGIO

Type Tipo	Bande carton / Cartone
Unité Unità	24 – display
Mesure Misura	14,5 x 38,5 x 38,5 cm
Volume Volume	0,021 m ³
Poids Peso	1,8 kg



MIS SUR PALETTE PALLET

Type Tipo	Européen/Europeo
Unité Unità	1536 und/pallet
Mesure Misura	120 X 80 X 169 cm
Distribution pallet / Distribuzione di pallet	4 couches/strati X 16 emballage /imballi: 64
Volume Volume	1,62 m ³
Poids Peso	138 kg

